

Product datasheet for **RG228360**

AADACL1 (NCEH1) (NM_001146276) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AADACL1 (NCEH1) (NM_001146276) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	NCEH1
Synonyms:	AADACL1; NCEH
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG228360 representing NM_001146276
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGAGCAGCTGCCGCGGGCAGAAAAGTTGCCGGAGGTCTCCGGGTGGTATCGCCCTTCTCTTTGCCAGC
 CCGCTGGCGAGCCGAGCCAGGGCAAGATGAGGTGTCCTGTGTCCTGCTACCGCCCTGGTGGCGCTGGC
 CGCCTATTACGTCTACATCCCGCTGCCTGGCTCCGTGTCCGACCCTGGAAGCTGATGCTGCTGGACGCC
 ACTTTCCGGGGTGACAGCAAGTGAAGTAACTGATCCACTACCTGGGACTGAGCCATCACCTGCTGGCAC
 TGAATTTTATCATTGTTTCTTTGGCAAAAAAGCGCGTGGTCTTCTGCCCAAGTGAAGGTGACCGACAC
 AGACTTTGATGGTGTGGAAGTCAAGTGTGTTGAAGGCCCTCCGAAGCCGAAGAGCCACTGAAACGCAGC
 GTCGTTTATATCCACGGAGGAGCTGGCCTTGGCAAGTGAAGTGCCTGCTGTCACCTCAGATGAAA
 TCAGGTATTATGATGAGCTGTGTACAGCAATGGCTGAGGAATTGAATGCTGTCTTGTGTTCCATTGAATA
 CAGGCTAGTTCCAAAGGTTTATTTTCTGAGCAAATTCATGATGTTGTACGGGCCACAAAGTATTTCTG
 AAGCCAGAAGTCTTACAGAAGTATATGGTTGATCCAGGCAGAATTTGCATTTCTGGTGACAGTGGTGGT
 GAAATCTGGCTGCTGCCCTTGGACAACAGTTTACTCAAGATGCCAGCCTAAAAATAAGCTCAAACCTACA
 AGCTTTAATTTATCCAGTCTTCAAGCTTTAGATTTTAAACACACCATCTTATCAGCAAAATGTGAACACC
 CCAATCTGCCCGCTATGTCATGGTGAAGTATTGGGTGGACTACTTCAAAGGCAACTATGACTTTGTGC
 AGGCAATGATCGTTAACAATCACACTTCACTTGTGTTGGAAGAGGCTGCTGCTGTGAGGGCCCGTCTAAA
 CTGGACATCCCTCTTGCCTGCATCCTTCACAAAGAACTACAAGCCTGTTGTACAGACCACAGGCAATGCC
 AGGATTTCCAGGAGCTTCTCAGTGTGTTGAGTGGCCTCCGCCCACTATTGCAGACCAGGCAGTGC
 TGCAGCTCCTCCAAAGACCTACATTCTGACGTGTGAGCATGATGTCCTCAGAGACGATGGCATCATGTA
 TGCCAAGCGTTTGGAGAGTGCCTGTGGAGGTGACCCTGGATCACTTTGAGGATGGCTTTCACGGATGT
 ATGATTTTCACTAGCTGGCCACCAACTTCTCAGTGGGAATCCGGACTAGGAATAGTTACATCAAGTGGC
 TAGATCAAACCTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG228360 representing NM_001146276
 Red=Cloning site Green=Tags(s)


MSSCRGQKAVAGGLRVVSPFPLCQPAGEPSQGMRSSCVLLTALVALAAYVYIPLPGSVSDPWKLMLLDA
 TFRGAQQVSNLIHYLGLSHLLALNFIIVSFGKSAWSSAQVKVTDTFDGVVRFEGPPKPEEPLKRS
 VVYIHGGWALASASWSPSDEIRYYDELCTAMAEELNAVIVSIEYRLVPKVYFPEQIHDVVRATKYFL
 KPEVLQKYMVDPGRICISGDSAGNLAALGQQFTQDASLKNKLQAL IYPVLQALDFNTPSYQQNVNT
 PILPRYVMVKYVDFYFKGNYDFVQAMIVNNHTSLDVEEAAVRARLNWTSLLPASFTKNYKPVVQTTGNA
 RIVQELPQLLDARSAPLIADQAVLQLLPKTYILTCEHDVLRDDGIMYAKRLESAGVEVTLDFHEDGFHGC
 MIFTSWPTNFSVGIRTRNSYIKWLDQNL

TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:

Plasmid Map:


ACCN: NM_001146276

ORF Size: 1344 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001146276.1</u> , <u>NP_001139748.1</u>
RefSeq Size:	4318 bp
RefSeq ORF:	1251 bp
Locus ID:	57552
UniProt ID:	<u>Q6PIU2</u>
Cytogenetics:	3q26.31
Protein Families:	Transmembrane
Gene Summary:	Hydrolyzes 2-acetyl monoalkylglycerol ether, the penultimate precursor of the pathway for de novo synthesis of platelet-activating factor. May be responsible for cholesterol ester hydrolysis in macrophages, thereby contributing to the development of atherosclerosis. Also involved in organ detoxification by hydrolyzing exogenous organophosphorus compounds. May contribute to cancer pathogenesis by promoting tumor cell migration.[UniProtKB/Swiss-Prot Function]