

Product datasheet for **RG228389**

KRT72 (NM_001146226) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KRT72 (NM_001146226) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	KRT72
Synonyms:	CK-72; K6irs; K6IRS2; K72; KRT6; KRT6IRS2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG228389 representing NM_001146226
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAGCCGCCAACTGACCCATTTCCCCCGGGGAGCGCCTGGGTTTTCAGCGGTTGCTCCGCGGTCTCTCT
 CTGGCGGGATCGGCAGCAGCTCCGCCTCATTCCGGGCCCGGGTCAAGGGCTCGGCCCTCTTTGGCAGCAA
 GAGCCTCTCTGCCTTGGGGGCAGCCGAAGCCTGGCGCTCAGCGCTGCTGCACGGCGGGCGCGGCCGC
 CTGGGCGGCTTCTGTGGGACCGCCTTCGGCAGCGCCGGGCTGGGGCCAAAGTGTCCCTCCGTGTGCCAC
 CCGGGGGCATCCCTCAGGTCACCGTCAACAAGAGCCTCCTGGCCCGCTCAACGTGGAGATGGACCCGA
 GATCCAGAGGGTGCAGCGCCAGGAGCGGGAGCAGATCAAGGCGCTAAACAACAAGTTGCGCTCCTTCATC
 GACAAGGTGCGGTTCTGGAGCAGCAGAATCAGGTGCTAGAGACCAAGTGAACCTCTACAGCAGCTGG
 ACTTGAACAACCTGCAGGAAGAACCTGGAGCCATTTATGAGGGCTACATCAGCAACCTGCAGAAGCAGCT
 GGAGATGCTGTCTGGGGACGGGTTGAGGCTGGATTCGGAGCTGAGGAACATGCAGGATTTGGTGGAGGAC
 TACAAGAAGAGGTATGAGGTGGAGATTAACAGACGCACAGCTGCTGAGAAATGAGTTTGTGGTGTCAAGA
 AGGACGTGGATGCTGCTTACATGAATAAGGTTGAGCTCCAGGCCAAGGTGGACTCCTTGACAGATGAGAT
 TAAATCTTCAAGTGCCTTTATGAAGGGGAGATCACTCAGATCCAGTCCCACATCAGCGACACGTCCATC
 GTCCTGTCAATGGACAACAACCGGGATCTGGACCTGGACAGCATCATTGCCGAGGTCCGTGCCCAGTACG
 AGGAGATTGCCCTAAAGAGCAAGGCCGAGGCTGAGACCTGTACCAGACCAAGTGTGCCGATCTGGAGAC
 GGCCATCGCCGACGCTGAACAGCGGGGGGACTGCGCCCTGAAAGATGCCCGGCCAAGCTGGATGAGCTG
 GAGGGCGCCCTGCACCAGGCCAAGGAGGAGCTGGCAGGATGCTGCGTGAAGTACCAGGAGCTCGTGAGCC
 TGAAGCTGGCCCTGGATATGGAGATCGCCACCTACCGCAAGCTGCTGGAGAGCGAGGAGTGCAGGATGTC
 TGGCGAATATCCAAATCTGTGAGCATCTCCGTCATCAGCAGCACCAATGCTGGGGCAGGAGGGGCTGGC
 TTCAGCATGGGCTTTGGCGCCTCAAGCAGTTATAGCTACAAAACCTGCAGCTGCAGACGTCAAGACAAAAG
 GCAGCTGTGGCAGTGAAGTCAAGGATCCCTTGCCAAAACCTCGGGGAGCAGCTGTGCCACAAAAGGC
 CTCCAGA

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG228389 representing NM_001146226
 Red=Cloning site Green=Tags(s)

MSRQLTHFPRGERLGFSGCSAVLSGGIGSSSASFRARVKGSASFSGSKLSCLGGSRLSALSAARRGGGR
 LGGFVGTAFGSAGLGPKCPSVCPGGIPQVTVNKSLAPLNVEMDPEIQRVRAQEREQIKALNNKFASFI
 DKVRFLEQQNQVLETKWNLLQQLDLNCRKNLEPIYEGYISNLQKQLEMLSGDGVRLDSELNMQDLVED
 YKKRYEVEINRRRTAAENEFVVLKDVDAAYMNKVELQAKVDSL TDEIKFFKCLYEGETIQISHISDTSI
 VLSMDNRRDLDLDSIIAEVRAQYEEIALKSKAEETLYQTKCADLETAIADAEQRGDCALKDARAKLDEL
 EGALHQAKEELARMLREYQELVSLKLALDMEIATYRKLLESEECRMSGEYPNSVSVISSTNAGAGGAG
 FSMGFGASSYSYKTAADVKTGSCGSELKDPLAKTSGSSCATKKASR

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_001146226

ORF Size: 1407 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001146226.2](#)

RefSeq Size: 1904 bp

RefSeq ORF: 1410 bp


Locus ID: 140807

UniProt ID: [Q14CN4](#)

Cytogenetics: 12q13.13

Gene Summary: Keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells. The type II keratins consist of basic or neutral proteins which are arranged in pairs of heterotypic keratin chains coexpressed during differentiation of simple and stratified epithelial tissues. This gene encodes a type II keratin that is specifically expressed in the inner root sheath of hair follicles. The type II keratins are clustered in a region of chromosome 12q12-q13. Multiple transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jun 2009]

Product images:



Circular map for RG228389