

## Product datasheet for **RG228402**

### **NEK3 (NM\_001146099) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	NEK3 (NM_001146099) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	NEK3
Synonyms:	HSPK36
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RG228402 representing NM\_001146099  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**

ATGGATGACTACATGGTCTGAGAATGATTGGGGAGGGCTCCTTCGGCAGAGCTCTTTGGTTTCAGCATG  
 AAAGCAGTAATCAGATGTTTGCCATGAAAGAAATAAGGCTTCCCAAGTCTTTCTCTAATACACAGAATTC  
 TAGGAAGGAGGCTGTTCTTTTAGCCAAAATGAAACACCCTAATATTGTTGCCTTCAAAGAATCATTGAA  
 GCTGAAGGACACTTGATATTGTGATGGAATACTGTGATGGAGGGGATCTAATGCAAAAGATTAACAGC  
 AGAAAGGAAAGTTATTTCTGAAGACATGATACTTAATTGGTTTACCCAATGTGCCTTGGAGTAAATCA  
 CATTACAAGAAACGTGTGCTACACAGAGATATCAAGTCCAAGAATATCTTCTCACTCAGAATGGAAAA  
 GTGAAATTGGGAGACTTTGGATCTGCCGCTTCTCTCCAATCCGATGGCATTGCTTGTACCTATGTGG  
 GAACTCCTTATTATGTGCCTCCAGAAATTTGGGAAAACCTGCCTTATAACAATAAAAGTGACATCTGGTC  
 CTTGGGTTGCATCTGTATGAACTCTGTACCCTTAAGCATCCATTCAGGCAAATAGTTGAAAAAATCTT  
 ATCCTCAAAGTATGTCAAGGGTGCATCAGTCCACTGCCGTCTCATTACTCTATGAACTTCAGTTCCTAG  
 TCAAGCAGATGTTTAAAAGGAATCCCTCACATCGCCCTCGGCTACAACGCTTCTCTCGAGGCATCGT  
 AGCTCGGCTTGTCCAGAAGTGCTTACCCCCGAGATCATCATGGAATATGGTGAGGAAGTATTAGAAGAA  
 AAAAAAATTCGAAGCATAACACACCAAGAAAAAACAAGAGGAAGAACAAGATAGAAAGGGTAGCCATA  
 CTGATTTGGAAGCATTAAAGAAATTTAGTTGAAAGTGCAATGAGAAGAGTAAACAGAGAAGAAAAAGG  
 TAATAAGTCAGTCCATCTGAGGAAAGCCAGTTCACCAAATCTTCATAGACGACAGTGGGAGAAAAATGTA  
 CCCAATACAGCTTTACAGCTTTGGAAAATGCATCCATACTCACCTCCAGTTAACAGCAGAGGACGATA  
 GAGGTGGTTCTGTAATAAAGTACAGCAAAAATACTACTCGTAAGCAGTGGCTCAAAGAGACCCCTGACAC  
 TTTGTTGAACATCCTTAAGAATGCTGATCTCAGCTTGGCTTTTCAAACATACACAATATATAGACCAGGT  
 TCAGAAGGTTCTTGAAGGCCCCCTGTCTGAAGAAACAGAAGCATCGGACAGTGTGATGGAGGTCACG  
 ATTCTGCATTTTGGATCCAGAGCGACTTGAGCCTGGGCTAGATGAGGAGGACACGGACTTTGAGGAGGA  
 AGATGACAACCCCGACTGGGTGTCAGAGCTGAAGAAGCGAGCTGGATGGCAAGGCCTGTGCGACAGA

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA

**Protein Sequence:**

>RG228402 representing NM\_001146099  
 Red=Cloning site Green=Tags(s)

MDDYMLRMIGEGSFRALLVQHESNQMFAMKEIRLPKFSNTQNSRKEAVLLAKMKHPNIVAFKESFE  
 AEGHLYIVMEYCDGGDLMQKIKQKGLFPEDMILNWFQMLGVNHIHKKRVLHRDIKSKNIFLTQNGK  
 VKLGDGFSARLLSNPMAFACTYVGTYPYYPPEIWENLPYNNKSDIWSLGCILYELCTLKHPFQANSWKNL  
 ILKVCQGCISPLPSHYSYELQFLVKQMFKRNPShRPSATTL SRGIVARLVQKCLPPEIIMEYGEEVLEE  
 IKNSKHNTPRKKQEEEDRKGSHDLESINENLVESALRRVNREEKGNKSVHLRKASSPNLHRRQWEKNV  
 PNTALTAENASILTSSLTAEDDRGGSVIKYSKNTTRKQWLKETPDTLLNILKNADLSLAFQTYTYIRPG  
 SEGFLKGPLSEETEASDSDVGGHDSVILDPERLEPGLDEEDTDFEEEDNDPVDWSELKKRAGWQGLCDR

**TRTRPLE** - GFP Tag - V

**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**


**ACCN:** NM\_001146099

**ORF Size:** 1467 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001146099.1](#), [NP\\_001139571.1](#)

**RefSeq Size:** 2345 bp

**RefSeq ORF:** 1470 bp

**Locus ID:** 4752


**UniProt ID:** [P51956](#)

**Cytogenetics:** 13q14.3

**Protein Families:** Druggable Genome, Protein Kinase

**Gene Summary:** This gene encodes a member of the NimA (never in mitosis A) family of serine/threonine protein kinases. The encoded protein differs from other NimA family members in that it is not cell cycle regulated and is found primarily in the cytoplasm. The kinase is activated by prolactin stimulation, leading to phosphorylation of VAV2 guanine nucleotide exchange factor, paxillin, and activation of the RAC1 GTPase. Two functional alleles for this gene have been identified in humans. The reference genome assembly (GRCh38) represents a functional allele that is associated with the inclusion of an additional coding exon in protein-coding transcripts, compared to an alternate functional allele that lacks the exon. [provided by RefSeq, Sep 2019]

**Product images:**



Circular map for RG228402