

Product datasheet for **RG228452**

TENR (ADAD1) (NM_001159295) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	TENR
Synonyms:	Tenr
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



ORF Nucleotide Sequence: >RG228452 representing NM_001159295
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGCTGAAAAAGAACCTGCCAGTTCACCAGCGACAAAGACGATAACTACACCCACAGGATGGTCCTCAG
 AAAGTTACGGCCTGTCCAAGATGGCATCCAAGGTTACGCAAGTAACGGGTAATTTCCAGAGCCGTTGCT
 TTCCAAGAATCTTTCATCTATTTCAAATCTGTCTTCTCCAAAAAATACCTAAGGAATTTATAATG
 AAATACAAACGTGGAGAGATAAATCCTGTGTGACGCTTGCACCAAGTTGCACAAATGCAGCGAGTTCAGC
 TTGACCTTAAGGAACTGTGACAACAGGTAATGTTATGGGACCATATTTGCCTTTTGTGCTGTGGTGA
 TGGTATTCAGTACAAGACTGGACTGGGACAAAAATAAAAGGAGTCTAGATCCAATGCAGCAAAATAGCT
 CTTGATGAGCTTCTACAACCTGGATGAACCTGAACCACGAATTTAGAAACATCAGGTCTCTCTCTTTCC
 CTGCAGAACCTGTTGTTTTATCTGAAGTACATATGTTTCAAAAGTACATTATGAAGGAGACATATTCA
 ATATGCAAGATTTACAGATCGTTAAAGAAAGATTTAATCAACTAATTTCTAATCGTTCAGAATACCTG
 AAATATAGCAGTTCATTGGCTGCTTTTATAATTGAAAGAGCTGGACAACATGAGGTTGTAGCTATAGGCA
 CAGGTGAATACAATTACAGCCAGGACATTAAGCCAGATGGAAGAGTATTGCATGACACTCATGCTGTTGT
 TACAGCAAGAAGGTCTCTTCTTAGGTACTTTTATAGACAACCTCTGCTCTTCTACAGCAAAAATCCTGCT
 ATGATGGAAAAATCAATATTTGTACAGAACCAACTTCTAATCTACTACTCTTAAACAGAAATATCAACA
 TTTGCCCTTTACATGAACCAAGTTGCCTAAAGGATCAGCCAGATTAAGTCACAGTTACGCTCTTAATCCACA
 TTCTATATCTGCATTTGAAGCAATGAAGAACTCTGTCTCCACGTGGCTGTTGAAGGCAAAATTTATCTG
 ACTGTTTACTGCTTAAAGATGGTGTAAAGAATAAGCAGTATGCTCAAGTGACAAATTGACCAGAT
 GGAAGTGTGTTGTTACAGGAGCATTGCTGAGTCATTTTATACAGCCAGTTTACATCAGCAGTATACT
 TATTGGTGATGGGAATTGACAGTGATACCAGAGGCTTAGAAATCGCTATAAAGCAACGTGTTGATGATGCA
 CTCACCTTCAAACTTCCAATGTTTTACTTAGTCAACAGACCTCATATTAGTTTAGTACCCTCTGCATATC
 CCCTTCAAATGAAGTTGGAATATAAATTTCTGAGTCTGAATGGGCACAAGGAGATGTTTCTTTGGAGAT
 TGTGGATGGCCTGAGTGGGAAGTCACTGAAAGTTCCCATTTAAAGTGGTATGTCAATGGCAAGTCGG
 CTTTGTAAAGCTGCAATGTTAAGTCGGTTTAACTGCTTGCCAAAGAAGCTAAAAAGAATTACTTGAAG
 CTGGTACATATCATGCAGCTAAGTGTATGTCTGCCTCTATCAAGAAGCTAAATGTAAGTTGAAATCCTA
 CTTACAACAACATGGCTATGGATCTGGATTGTGAAATCTCCCTGCATAGAGCAATTTAACATG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG228452 representing NM_001159295
 Red=Cloning site Green=Tags(s)

MLKKNLPVQPATKTIITPTGWSSSESYGLSKMASKVTQVTGNFPEPLL SKNLSSISNPVLPKKIPKEFIM
 KYKRGIEINPVSAHQFAQMQRVQLDLKETVTTGNVMGPYFAFCVVDGIQYKTLGLGQNKESRSNAAKLA
 LDELLQLDEPEPRILETSGPPFPAPVVLSELAYVSKVHYEGRHIQYAKISQIVKERFNQLISNRSEYL
 KYSSSLAAFIIERAGQHEVVAIGTGEYNYSQDIKPDGRVLHDTHAVVTARRSLLRYFYRQLLLFYSKNPA
 MMEKSIFCTEPTSNLLTLKQININICLYMNQLPKGSAQIKSQLRLNPHSISAFEANEELCLHVAVEGKIYL
 TVYCPKDGVNRISSMSSDKLTRWEVLGVQGALLSHFIQPVYISSILIGDGNCSDTRGLEIAIKQRVDDA
 LTSKLPMPFYL VNRPHISLVPSAYPLQMNLEYKFLSLNWAQGDVSL EIVDGLSGKITESSPFKSGMSMASR
 LCKAAML SRFNLLAKEAKKELLEAGTYHAAKCMSASYEAKCKLKS YLQQHGYGSWIVKSPCIEQFNM

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001159295

ORF Size: 1674 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001159295.2](#)

RefSeq Size: 2012 bp

RefSeq ORF:	1677 bp
Locus ID:	132612
UniProt ID:	Q96M93
Cytogenetics:	4q27
Protein Families:	Druggable Genome
Gene Summary:	Plays a role in spermatogenesis. Binds to RNA but not to DNA (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG228452