

Product datasheet for **RG228522**

CCDC120 (NM_001163321) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC120 (NM_001163321) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CCDC120
Synonyms:	JM11
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RG228522 representing NM_001163321
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

ATGGAAGTCAAAGGTCAGCTGATCAGCTCTCCTACCTTCAATGCCCCAGCTGCCCTGTTCCGAGAGGCTG
 CCCCCAGGTGAAGTCAGAGCGTCTCGGGGGCTGCTTGACCGGCAGCGACCCTGCAGGAGGCCCTGAG
 CCTGAAACTTCAGGAGCTCCGCAAAGTGTGTCTCCAGGAGCGGAGCTGACTGGCCAGCTGCCCCCTGAG
 TGCCCACTAGAGCCTGGTGAACGGCCCCAGTTGGTCCGCCGGCGGCCCCACAGCCCGCGCCTACCCTC
 CACCGCACCCCAACCAAGCACACCCTCCTTGTGCCCTGCTGAGGAGCTGGCTCTTGAGGCCCTGGAACG
 CGAGGTGTCAAGTCAACAGCAGATCGCGCGGCCGCCCGCCGCTGGCCTTGCCCTGATCTGAGCACC
 GAGCAGCGCCGGCGCGGCCAGGTCCAGGCAGATGACTGAGGAGGCTGCATGAGCTAGAGGAGCAGC
 TCAGGGATGTCCGGGCCCGCCTTGGCCTCCAGTGTCTCCGCTGCCCCAGCCACTGCCACTGTCCACGGG
 GTCAGTGATCACCAACCGAGGAGTCTGCCTGGGCATGCGTCTTGCTCAGCTCAGCAAGAGGACGTAGTT
 CTGCACTCAGAGAGCAGCTCCCTCTCAGAGTCTGGGGCCAGCCATGACAACGAGGAGCCCATGGCTGCT
 TCTCTCTGGCCGAGCGCCCTCACCAACCAAGGCTTGGGACCAGCTGCGGGCAGTATCTGGGGGAGCCC
 TGAGCGGCGAACCCCATGGAACCACTCCATCAGATCTTTATGGGGATCTGAAGAGCCGGCGAACTCT
 GTGGCCAGCCCCACAGCCCCACAGCTCGCTGCCAGGAGTGCTCCAGTTTGAAGGGCGAAGTGTGC
 CTGCCACCCCTGTCTCACCCGGGGCGCTGGCCCCAGCTCTGCAAACTGAAGGCCTTCATTCTCGTCA
 GTGGTCCGGCAGCCAGGACTCCAGATGGGCTTCCCCGGGGGAGCCCTGCCTCCGATCGCGCTCCCTC
 TTCGTAGCTCGCACCCGGCCAGCAACAGTTCTGAGGCCCTGCTGGTGGACCGGGCCGCTGGTGGGGAG
 CTGGCTCCCCGCTGCCCTCTGGCTCCCTCTGCCTCTGGCCCCCAGTCTGCAAGAGCAGTGAGGTGCT
 GTATGAGCGCCCCAACCAACCCCTGCCTTCTCTCCCGCACAGCAGGCCCCCAAGACCCTCCCCGGGCC
 GCGCGGCTAGCTCAGCTGCCCTGCCTCCCGAGGTGCCCGCGGCTCCACCTGTGTGGAGACTTCC
 TCTTGACTATTCTTGAGCCGGGGCTGCCCGCAGTGCGGCTGGAACAGGCTGGGGGAGCTGCCGCC
 TGCAGCTGAGGTCCAGGACCCCTCTCCCGCGGGATGGGCTCCTACCATGCTCCCGGCCACCACT
 GTGTATGCAGCTGACAGCAACAGCCCCCTCTCCGCACCAAGGACCCCAACCCGTGCCACCCGCACTA
 AGCCCTGTGGCTGCCCCAGAGGCTGCCGAAGGCCCTGAGGTGCATCCAAACCTCTGCTGTGGATGCC
 CCCACCCACCGTATCCCTCGGCTGGTGAACGCAGTGCCACAAGAACCTGGCTCTGGAGGGGCTGCGG
 GACTGGTACATCCGGAACCTCGGACTGGCTGCGGGGCCCCAGCGCCGGCCTGTGCTCCCTCCGTGGGCC
 CGCCACACCCACCTTCTCCATGCCCGCTGCTATGAGGTGGGCCAGGCGCTGTACGGGGCCCCAGCCA
 GGGCCACTCCCACTCGAGGAGTTTACGGCGCCCCCTGTCTCTGGCAGATACTACGGGACTTCCTG
 TATCCCCCGAGCTGAGCGCTCGTTAAGTGACCTGACGCTAGAGGGGGAGCAGTCTCCAGTTCTGACA
 CCCAGACCCCGGGGACACTGGTC

ACGCGTACGCGGGCGCTCGAG – GFP Tag – GTTTAA


Protein Sequence: >RG228522 representing NM_001163321
 Red=Cloning site Green=Tags(s)

```
MEVKGQLISSPTFNAPALFGEAAPQVKSERLRGLLDQRQLQEALSLKLQELRKVCLQEALTGQLPPE
CPLEGERPQLVRRRPPTARAYPPHPNQAHSLCPAEELALEALEREVSQQQIAAAARRLALAPDLST
EQRRRRRQVQADALRRLHELEEQLRDVRARLGLPVLPLPQPLPLSTGSVITTQGVCLGMRLAQLSQEDVV
LHSESSSLSESGASHDNEEPHGCFSLAERPSPPKAWDQLRAVSGGSPERRTPWKPPPSDLYGDLKSRNS
VASPTSPTRSLPRSASSFEGRSVPATPVLTRGAGPQLCKPEGLHSRQWSGSQDSQMGFPADPASDRASL
FVARTRRSNSSEALLVDRAAGGGAGSPPAPLAPSASGPPVCKSSEVL YERPQPTPAFSSRTAGPPDPRA
ARPSSAAPASRGAPRLPPVCGDFLLDYSLDRGLPRSGGGTGWGELPPAAEVPGLSRRDGLLTMLPGPPP
VYAADNSPLLRTKDPHTRATRTKPCGLPPEAAEGPEVHPNPLLWMPPTTRIPSAGERSGHKNLALEGLR
DWYIRNSGLAAGPQRRPVLPSVGPPHPFLHARCYESVQALYGAPSQAPLPHSRSFAPPVSGRYADFL
YPPELSARLSDLTLEGEQSSSDTQTPGTLV
```

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001163321

ORF Size: 1983 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001163321.1](#), [NP_001156793.1](#)

RefSeq Size: 3626 bp

RefSeq ORF: 2091 bp


Locus ID: 90060

UniProt ID: [Q96HB5](#)

Cytogenetics: Xp11.23

Gene Summary: This gene encodes a protein that contains a coiled-coil domain. Several alternatively spliced transcript variants encoding multiple isoforms have been found for this gene. [provided by RefSeq, Nov 2012]

Product images:



Circular map for RG228522