

Product datasheet for **RG228548**

ZNF611 (NM_001161500) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF611 (NM_001161500) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ZNF611
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG228548 ORF sequence, **codon optimized**.
Due to the complexity of NM_001161500, the ORF clone is codon optimized for mammalian Expression.
The nucleotide sequence differs from the reference sequence, yet the amino acid sequence remains identical.

Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCTCCGGGAAGAGGCAGCACAGAAACGCAAGGGCAAGGAACCCGGAATGGCCCTGCCGAGGGTCGCC
TGACCTTCGGGATGTTGCTATCGAGTTTTCTACTGGCAGAGTGGAAGTGTCTCAACCCAGCCAGAGAGC
CCTTTATCGGGAAGTCATGCTGGAGAATACAGGAACCTTGAAGCGGTGGACATTTCCAGTAAGTGCATG
ATGAAGGAGGTTCTCAGTACAGGACAGGAAATACGGAGGTTATACACACCGGTACTCTGCAGAGGCAGC
AGTCCCATCACATTGGCGACTTCTGTTTCCAGGAAATAGAGAAGGAGATACATGACATCGAGTTCAGTG
TCAAGAAGACGAACGCAACGGACTGGAGGCCCCCATGACTAAGATCAAGAAGCTCACTGGATCAACAGAC
CAGCACGACCACCAGACGCGGGCAATAAGCCCATCAAAGACCAGCTCGGGTCCAGCTTTTACTCACATC
TGCTGAACTTCACATCTTTAGATTAAAGGGGAAATTGGGAACCAAGTGGAGAAATCCACAAACGACGC
CCCATCCGTGTCCACCTTCCAGCGGATCTCTGCAGGCCGAGACTCAAATCTCTAATAACTACGGTAAT
AATCCATTGAACAGCTCTCTGCTTCCAGAACGAGGAGTACATATGCGAGAGAAAAGCTTCCAGTGTA
ATAAGTCCGAAAGGCTTTTAACTGCTCCAGCCTGTTGAGGAAACACCAGATACCACATCTCGGAGACAA
GCAGTACAAATGTGATGTGTGCGGAAATTTGTTCAATCATGAACAATATCTGGCTTGTCCAGTGGTGC
CACACAGTGGAAAAGCCCTATAAATGTAAGGAGTGGTAAAGACATTTAGTCAGGAGTCTAGTCTTACAT
GTCATAGACGCTGCACACCGGTGTTAAGAGGTACAATTGCAATGAATGCGGTAAGATCTTCGGACAGAA
TTCAGCCCTCCTTATTGACAAAGCCATAGACACCGGGGAAAATCCATACAAGTGAACGAGTGTGACAAA
GCGTTCAACCAACAGTCACAGCTGTACATCATCGGATCCACACCGGAGAGAAAGCCTTACAAATGTGAGG
AATGCGACAAAGTCTTCTCCAGAAAGTCCACCATCGAAACTCACAAACGAATTCACACAGGAGAGAAACC
CTACCGGTGTAAGTGTGATACAGCCTTTACATGGCACTCTCAATTGGCTAGGCACAGACGGATCCAT
ACGGCAAAGAAGACCTATAAGTGAACGAATGCGGAAAAACATTCTCACACAAAAGCAGTCTCGTGTGTC
ATCACAGGCTTTCATGGCGGAGAGAAGTCATATAAGTGCAAGGTGTGCGATAAGGCTTTTGTGGTCCAG
TCAGCTTGCAAAGCACACCAGGATCGATTGCGGAGAAAAACCTTATAAATGTAACGAGTGTGGAAAGACC
TTTGGCCAAAACAGCGATCTGTTGATACACAAATCAATCCACACCGGAGAGCAACCGTACAAGTGTGATG
AGTGTGAAAAGGTGTTTTCTCGAAAATCTAGCCTGGAGACACACAAAATTTGGACACACCGGTGAGAAGCC
CTATAAATGCAAAGTGTGTGACAAAGCATTGCTTGTCACTCCTATTTGGCCAAACACACGGAATTCAC
TCAGGCGAAAAGCCATACAAATGCAACGAATGCTCTAAGACCTTTTCCACCGATCCTATCTCGTCTGCC
ACCATAGAGTTCACTCCGGAGAAAAGCCATACAAATGTAACGAGTGCAGCAAGACATTTAGCCGGAGGAG
CAGCCTGCACTGCCATCGGAGACTGCACAGCGGGGAAAAGCCCTACAAATGCAACGAGTGTGGCAACACT
TTCAGGCACTGTTCTTCTGTATCTACCACCGCAGACTGCATACTGGGGAAAAGTCTTATAAATGTACGA
TTTGTGATAAAGCCTTTGTGAGAACTCTTCTTTCTCGCCATACTAGGATCCATACGGCAGAGAAAGCC
TTACAAGTGAATGAATGCGGAAAGCCTTCAACCAGCAGTCTCACCTGTCCCGACATCACAGAATTCAT
ACGGGCGAGAAGCCA

ACGGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG228548 representing NM_001161500
 Red=Cloning site Green=Tags(s)


MLREAAQKRKGKEPGMALPQGRLTFRDVAIEFSLAEWKCLNPSQRALYREVMLENYRNLEAVDISSKCM
 MKEVLSTGQGNTEVIHTGTLQRHESHIGDFCFQEIEKEIHDIEFQCQEDERNLEAPMTKIKKLTGSTD
 QHDHRHAGNKPIKDQLGSSFYSHLPELHIFQIKGEIGNQLEKSTNDAPSVSTFQRI SCRPTQISNNYGN
 NPLNSSLLPQKQEVHMREKSFQCNKSGKAFNCSSLLRKHQIPHLGDKQYKCDVCGKLFNHEQYLACHDRC
 HTVEKPYKCKEKGKTFSEQSSLTCHRRLLHTGVKRYNCNECGKIFGQNSALLIDKAIDTGENPYKCECDK
 AFNQSSQLSHHRIHTGEKPYKCECDKVF SRKSTIETHKRIHTGEKPYRCKVCDTAFTWHSQARHRIH
 TAKKTYKNECGKTF SHKSSLVCHHRLHGGEKSYCKVCDKAFVWSSQLAKHTRIDCGEKPYKNECGKT
 FGQNSDLLIHKSIHTGEQPYKCECEKVF SRKSSLETHKIGHTGEKPYKCKVCDKAFACHSYLAKHTRIH
 SGEKPYKNECSKTF SHRSYL VCHHRVHSGEKPYKNECSKTF SRRSSLHCHRRLSGEEKPYKNECGNT
 FRHCSSLIYHRRLLHTGEKSYKCTICDKAFVRNLLSRHTRIHTAEKPYKNECGKAFNQSHLSRHRIH
 TGEKP


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:


ACCN: NM_001161500

ORF Size: 2115 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001161500.1](#), [NP_001154972.1](#)

RefSeq Size: 4417 bp

RefSeq ORF: 2118 bp

Locus ID: 81856

UniProt ID: [Q8N823](#)

Cytogenetics: 19q13.41
Protein Families: Transcription Factors
Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]