

## Product datasheet for **RG228557**

### ZC3H14 (NM\_001160104) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZC3H14 (NM_001160104) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ZC3H14
Synonyms:	MRT56; MSUT-2; NY-REN-37; SUT2; UKp68
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>RG228557 representing NM\_001160104  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGGAGATCGGCACCGAGATCAGCCGAAGATCCGGAGTGCCATTAAGGGGAAATTACAAGAATTAGGAG  
CTTATGTTGATGAAGAACTTCTGATTACATTATGGTGATGGTGGCCAACAAGAAAAGTCAGGACCAAAT  
GACAGAGGATCTGTCCCTGTTTCTAGGGAAACAACAATTCGATTCACCGTATGGCTTCATGGTGTATTA  
GATAAACTTCTGCTCTGTTACAACCTGAACCTCTAGTCTGAAGTCTTCTGATACCAACATCTTTGATAGTA  
ACGTGCCCTTCAAACAAGAGCAATTTTCAGTCGGGGAGATGAGAGGAGGCATGAAGCTGCAGTGCCACCACT  
TGCCATTCCTAGCGCGAGACCTGAAAAAAGAGATTCCAGAGTTTCTACAAGTTCGCAGGAGTCAAAAACC  
ACAAATGTCAGACAGACTTACGATGATGGAGCTGCAACCCGACTAATGTCAACAGTGAAACCTTTGAGGG  
AGCCAGCACCCCTCTGAAGATGTGATTGATTAAGCCAGAACCAGATGATCTCATTGACGAAGACCTCAA  
CTTTGTGCAGGAGAATCCCTTATCTCAGAAAAACCTACAGTGACACTTACATATGGTTCTCTCGCCCT  
TCTATTGAAATTTATCGACCACCTGCAAGTAGAAATGCAGATAGTGGTGTTCATTTAAACAGGTTGCAAT  
TTCAACAGCAGCAGAATAGTATTCATGCTGCCAAGCAGCTTGATATGCAGAGTAGTTGGGTATATGAAAC  
AGGACGTTTGTGTGAACCAGAGGTGCTTAACAGCTTAGAAGAAACGTATAGTCCGTTCTTTGAAACAAC  
TCGGAGAAAAATGAGTATGGAGGATGAAAACCTTCGGAAGAGAAAGTTGCCTGTGGTAAGTTCAAGTTGTTA  
AAGTAAAAAATCAATCATGATGGAGAAGAGGAGGAAGAAGATGATGATTACGGGTCTCGAACAGGAAG  
CATCTCCAGCAGTGTGTCTGTGCCTGCAAAGCCTGAAAGGAGACCTTCTCTTCCACCTTCTAAACAAGCT  
AACAAGAATCTGATTTTGAAGGCTATATCTGAAGCTCAAGAATCCGTAACAAAAACAACCTAACTACTCTA  
CAGTTCACAGAAACAGACACTTCCAGTTGCTCCAGAACTCGAACTTCTCAAGAAGAATTGCTAGCAGA  
AGTGGTCCAGGGACAAAGTAGGACCCCAAGAAATAGTCCCCCATTAAGAAGAGGAAACAAAAGGAGAT  
TCTGTAGAAAAAATCAAGGAACTCAACAGAGGCAATTATTATCCCGACTGCAAAATCGACCCAGTAAATGG  
CAGAAACTCTGCAGATGAGTCAAGATTACTATGACATGGAATCCATGGTCCATGCAGACACAAGATCATT  
TATTCTGAAGAAGCCAAAGCTGTCTGAGGAAGTAGTAGTGGCACCACCAAGAGTCGGGGATGAAGACT  
GCAGATTCCTTCGGGTACTTTCAGGACACCTTATGCAGACACGAGATCTTGTACAACAGATAAACCTG  
CAAGTCCCAAGTTTATAGTGACGCTGGATGGTGTCCCCAGCCCCCAGGATACATGTCAGATCAAGAGGA  
GGACATGTGCTTTGAAGGAATGAAACCCGTAACCAAACTGCAGCCTCAAACAAGGGACTCAGAGGTCTC  
CTCCACCACAGCAGTTGCACTTGTGAGCAGGCAGCTTGAAGACCCAAATGCTGAGATGAGTGAAGTGA  
GTGTGGCACAGAAACCAGAAAACTTTTGGAGCGCTGCAAGTACTGGCCTGCTTGTAAAAATGGGGATGA  
GTGTGCCTACCATCACCCATCTCACCTGCAAAGCCTTCCCAATTGTAATTTGCTGAAAAATGTTTG  
TTTGTTACCCAAATGTAATATGATGCAAAGTGTACTAAACCAGATTGTCCTTCACTCATGTGAGTA  
GAAGAATCCAGTACTGTCTCCAAAACAGTTGCACCACCAGCACCTTCCAGTAGTCAGCTCTGCCG  
TTACTTCCCTGCTTGAAGAAGATGGAATGTCCCTTCTATCATCCAAAACATTGTAGGTTAACACTCAA  
TGTACAAGACCGGACTGCACATTCTACCATCCACCATTAATGTCCACCACGACATGCCTTGAATGGA  
TTCGACCTCAAACCAGCGAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >RG228557 representing NM\_001160104  
 Red=Cloning site Green=Tags(s)

```

MEIGTEISRKIRSAIKGKLQELGAYVDEELPDYIMVMVANKKSQDQMTEDLSLFLGNNTIRFTVWLHGVL
DKLRSVTTEPSSLKSSDTNIFDSNVPSNKSNSFRGDERRHEAAVPLAIP SARPEKRDSRVSTSSQESKT
TNVRQTYDDGAATRLMSTVKPLREPAPSEVDIDIKPEPDDLIDEDLNFVQENPLSQKKPTVTLTYGSSRP
SIEIYRPPASRNADSGVHLNRLQFQQQNSIHAAKQLDMQSSWVYETGRLCEPEVLNSLEETYSPPFRNN
SEKMSMEDENFRKRKLPVVSSVVKVKKFNHDGEEEEEDDDYGSRTGSISSVSVPAPKERRPSLPPSKQA
NKNLILKAISEAQESVTKTTNYSTVPQKQTLPVAPRTRTSQEELLAEVVQGQSRTPRISPPIKEEETKGD
SVEKNQGTQQRQLL SRLQIDPVMAETLQMSQDYDMESMVHADTRSFILKKPKLSEEVVAVPNQESGMKT
ADSLRVL SGHLMQTRDLVQDPKPASP KFIVTLDGVPSPPGYMSDQEEDMCFEGMKPVNQTAASNKGLRGL
LHPQQLHLLSRQLED PNAEMSEL SVAQKPEKLLERCKYWPACKNGDECA YHHPISPCKAFPNC KFAEKCL
FVHPNCKYDAKCTKPD CPTHVSRRI PVLSPKPVAPPAPSSSQLCRYFPACKKMECPFYHPKHCRFNTQ
CTRPDCTFYHPTINVPPRHALKWIRPQTSE
  
```


TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**Plasmid Map:**


**ACCN:** NM\_001160104

**ORF Size:** 2190 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001160104.2](#)

**RefSeq Size:** 4190 bp

**RefSeq ORF:** 2193 bp

**Locus ID:** 79882

UniProt ID: [Q6PJT7](#)

Cytogenetics: 14q31.3

Gene Summary: The protein encoded by this gene is a poly(A)-binding protein that can affect gene expression and poly(A) tail length. The encoded protein may influence mRNA stability, nuclear export, and translation. [provided by RefSeq, May 2016]