

Product datasheet for **RG228751**

SSX2 (SSX2B) (NM_001164417) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: SSX2 (SSX2B) (NM_001164417) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: SSX2B
Synonyms: CT5.2; CT5.2b; HOM-MEL-40; SSX
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG228751 representing NM_001164417
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAACGGAGACGACGCCTTTGCAAGGAGACCCACGGTTGGTGCTCAAATACCAGAGAAGATCCAAAAGG
 CCTTCGATGATTTGCCAAATACTTCTCTAAGGAAGAGTGGGAAAAGATGAAAGCCTCAGAGAAAATCTT
 CTATGTGTATATGAAGAGAAAGTATGAGGCTATGACTAACTAGGTTTCAAGGCCACCCTCCACCTTTC
 ATGTGTAATAAACGGGCCGAAGACTTCCAGGGGAATGATTTGGATAATGACCCTAACCGTGGGAATCAGG
 TTGAACGTCCTCAGATGACTTTTCGCGAGGCTCCAGGGAATCTCCCGAAGATCATGCCCAAGAAGCCAGC
 AGAGGAAGGAAAATGATTCGGAGGAAGTCCAGAAAGCATCTGGCCACAAAATGATGGGAAAGAGCTGTGC
 CCCCCGGGAAAACCAACTACCTCTGAGAAGATTACGAGAGATCTGGACCCAAAAGGGGGGAACATGCCT
 GGACCCACAGACTGCGTGAGAGAAAACAGCTGGTGATTTATGAAGAGATCAGCGACCCTGAGGAAGATGA
 CGAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG228751 representing NM_001164417
 Red=Cloning site Green=Tags(s)

MNGDDAFARRPTVGAQIPEKIQKAFDDIAKYFSKEEWEKMKASEKIFYVYMKRKYEAMTKLGFKATLPPF
 MCNKRAEDFQGNLDNDPNRGNQVERPQMTFGRLQGISPIMPKKPAEEGNDSEEVPEASGPQNDGKELC
 PPGKPTTSEKIHRSRGPKRGEHAWTHRLRERKQLVIYEEISDPEEDE

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001164417.1](#), [NP_001157889.1](#)

RefSeq Size: 1306 bp

RefSeq ORF: 567 bp

Locus ID: 727837

UniProt ID: [Q16385](#)

Cytogenetics: Xp11.22

Gene Summary: The product of this gene belongs to the family of highly homologous synovial sarcoma X (SSX) breakpoint proteins. These proteins may function as transcriptional repressors. They are also capable of eliciting spontaneous humoral and cellular immune responses in cancer patients, and are potentially useful targets in cancer vaccine-based immunotherapy. This gene, and also the SSX1 and SSX4 family members, have been involved in t(X;18)(p11.2;q11.2) translocations that are characteristically found in all synovial sarcomas. This translocation results in the fusion of the synovial sarcoma translocation gene on chromosome 18 to one of the SSX genes on chromosome X. The encoded hybrid proteins are likely responsible for transforming activity. Alternative splicing of this gene results in multiple transcript variants. This gene also has an identical duplicate, GeneID: 6757, located about 45 kb upstream in the opposite orientation on chromosome X. [provided by RefSeq, Jul 2013]