

Product datasheet for **RG228762**

APBB2 (NM_001166054) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	APBB2 (NM_001166054) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	APBB2
Synonyms:	FE65L; FE65L1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG228762 representing NM_001166054
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGTCAGAAGTACTTCCAGCTGACTCAGGTGTTGACACCTTGGCAGTGTTTATGGCCAGCAGCGGAATA
CAGACGTACAAAATCGGAACAGCCAGCCACACCACAAAACCCCTTAACCTCCGATCCTCCACAATGA
ACTGTTGAACGCTGAAATAAAACACACAGAAAACCAAGAACAGCACACCTCCCAAATGCAGGAAAAATAT
GCACTAACTAACATCCAGGCGGCCATGGCCTCTCGGATCCAGCTGCACAGCCCTGCTGGGAAATGGCT
CTGCCAACATCAAGCTGGTGAATAATGGGGAGAACCAGCTCCGTAAGGCTGCAGAGCAAGGGCAGCAGGA
CCCCAACAAAACCTGAGCCCACTGCAGTCATCAACATAACTTCTGAGAAGTTAGAGGGTAAAGAGCCC
CACCCACAGGATTCCTCGAGCTGTGAGATTTTACCCTCCAGCCAGGAGAACTAAGAGCTTCCTAAATT
ACTATGCAGATCTGGAACCTCAGCCAGAGAAGTACAGCAGAACCAAGGCAATCACCATGGGACTGCGGA
AGAGAAATCCCAGCCAGTCCAGGGCCAGGCCTCCACCATCATTGGGAATGGTGATTTGCTGCTGCAGAAA
CCAAACAGACCCAGTCCAGCCCTGAAGACGGCCAAGTAGCCACAGTGTCATCCAGCCAGAAAACCAAGA
AGGATCATCCGAAAACAGGGGCCAAAACCGACTGTGCACTGCACCCGATCCAGAACCTGGCACCGAGCGA
TGAGGAGTCCAGCTGGACAACGTTGTCCCAAGACAGTGCCTCACCCAGCTCCCGGATGAAACAGATATA
TGGAGTGATCACTCATTTCAGACTGATCCAGATTTGCCGCTGGCTGGAAAAGAGTCAGTGACATTGCCG
GGACCTATTATTGGCACATCCCAACAGGAACGACTCAGTGGGAACGGCCCGTCTCCATCCCAGCAGATCT
CCAGGGTCTAGGAAAGGGTCACTTAGTTCTGTAAACGCCATCTCCACCCAGAGAACGAGGATTTGCAT
GCAGCCACTGTTAACCCGACCCAGTTTAAAAGAGTTTGAAGGAGCAACCCTACGCTATGCATCTTTGA
AACTCAGAAATGCCCCACACCCTGATGATGATTTCTGTAGTATCAACAGTGACCCAGAACCAAGGAGT
TTTTGCTGTGCGTTCTCTGGGATGGGTAGAGATGGCAGAAGAGGACCTCGCCCCGGTAAAAGTAGTGTT
GCGGTCAACAACCTGCATCAGGCAACTTCTACTGCAAAAATGACATCCGAGACACAGTCGGGATTTGGG
GAGAGGGGAAAGACATGTACCTGATCCTGGAGAATGACATGCTCAGCCTGGTGGACCCCATGGACCGCAG
CGTGCTGCACTCGCAGCCATCGTCAGCATCCGCGTGTGGGGCGTGGCCCGCACAATGGCCGGGATTTT
GCTTATGTAGCAAGAGATAAAGATAACAAGATTTTGAATGTCATGTATTTGATGTGACACACCAGCAA
AAGCCATTGCCACAAGTCTCCACGAGATCTGCTCCAAGATTATGGCTGAACGGAAGAATGCCAAAGCGCT
GGCCTGCAGCTCCTTACAGGAAAGGGCCAATGTGAACCTCGATGTCCTTTGCAAGATTTTCCAACACCA
AAGACTGAGCTGGTCCAGAAGTCCACGTGCACTACTTGGGCATGTTACCTGTAGACAAACCAGTCGGAA
TGATATTTTGAACAGTGCCATAGAAAATCTTATGACCTCATCCAACAAGGAGGACTGGCTGTGAGTGA
CATGAACGTGGCTGATGCCACTGTGACTGTCATCAGTGAAGAATGAAGAGGAAGTCTTAGTGGAAATGT
CGTGTCGATTCTGTCTTTCATGGGTGTTGGGAAGGACGTCCACACATTTGCCTTCATCATGGACACGG
GGAACCAGCGCTTTGAGTGCCACGTTTTCTGGTGCAGCCTAATGCTGGTAACGTGTCTGAGGCGGTGCA
GGCCGCTGCATGTTACGATATCAGAAGTGCTTGGTAGCCAGGCCGCTTCTCAGAAAAGTTCGACCACCT
CCACCGCCAGCAGATTAGTAACCAGAAGAGTCAACCAATGTAACGAGGGGTCTTATCCCTCATTG
ACACTTTGAAACAGAAAACGCCCTGTACCGAAATGCCA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG228762 representing NM_001166054
 Red=Cloning site Green=Tags(s)

```
MSEVLPADSGVDTLAVFMASSTGTTDVTNRNSPATPPNTLNLRSSHNELLNAEIKHTETKNSTPPKCRKKY
ALTNIIQAAMGLSDPAAQPLLGNLSANIKLVKNGENQLRKAEEQGGQDPNKNLSPTAVINITSEKLEGKEP
HPQDSSSCEILPSQPRRTKSFNLNYADLETSARELEQNQGNHHTAEEKSQPVQGGASTIIGNDLLLLQK
PNRPQSSPEDGGVATVSSSPETKKDHPKTGAKTDCALHRIQNLAPSDEESSWTTLSQDSASPSPDETDI
WSDHSFQTDPDLPPGWKRVSDIAGTYWHIPTGTTQWERPVSIADLQGSRKGSLSVTPSPPTPENEDLH
AATVNPDPPLKEFEGATLRYASLKLRNAPHDDDDSCSINSNDPEAKCFVRSGLGWVEMAEDLAPGKSSV
AVNNCIRQLSYCKNDIRDVTGVIWEGEKDMLYLENDMLSLVDPMDRSVLHSQPIVSIRVWVGRDNRDF
AYVARDKDTRILKCHVFRCDTPAKAIATSLHEICSKIMAERKNAKALACSSLQERANVLDVPLQDFPTP
KTEL VQKFHVQYLGMLPVDPKPVGMDILNSAIENLMTSSNKEDWLSVNMNVADATVTVISEKNEEVLVEC
RVRFLSFMGVGKDVHTFAFIMDTGNQRFECVFWCEPNAGNVSEAVQAACMLRYQKCLVARPPSQKVRPP
PPPADSVTRRVTTNVKRGVLSLIDTLKQKRPVTEMP
```


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:


ACCN: NM_001166054

ORF Size: 630 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001166054.1](#), [NP_001159526.1](#)

RefSeq Size: 7012 bp

RefSeq ORF: 633 bp

Locus ID: 323

UniProt ID: [Q92870](#)

Cytogenetics: 4p14-p13

Protein Families: Transcription Factors

Gene Summary: The protein encoded by this gene interacts with the cytoplasmic domains of amyloid beta (A4) precursor protein and amyloid beta (A4) precursor-like protein 2. This protein contains two phosphotyrosine binding (PTB) domains, which are thought to function in signal transduction. Polymorphisms in this gene have been associated with Alzheimer's disease. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2009]