

Product datasheet for **RG228808**

NAT1 (NM_001160170) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NAT1 (NM_001160170) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	NAT1
Synonyms:	AAC1; MNAT; NAT-1; NATI
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG228808 representing NM_001160170 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACATTGAAGCATATCTTGAAAGAATTGGCTATAAGAAGTCTAGGAACAAATTGGACTTGGAAACAT
TAACTGACATTCTTCAACACCAGATCCGAGCTGTTCCCTTTGAGAACCTAACATCCATTGTGGGGATGC
CATGGACTTAGGCTTAGAGGCCATTTTTGATCAAGTTGTGAGAAGAAATCGGGTGGATGGTGTCTCCAG
GTCAATCATCTTCTGACTGGCTCTGACCACTATTGGTTTTGAGACCACGATGTTGGGAGGGTATGTTT
ACAGCACTCCAGCCAAAAATACAGCACTGGCATGATTCACCTTCTCCTGCAGGTGACCATTGATGGCAG
GAACTACATTGTCGATGCTGGGTTTGGACGCTCATACCAGATGTGGCAGCCTCTGGAGTTAATTTCTGGG
AAGGATCAGCCTCAGGTGCCTTGTGCTTCCGTTTACGGAAGAGAATGGATTCTGGTATCTAGACCAAA
TCAGAAGGGAACAGTACATTCCAATGAAGAATTTCTTCAATCTGATCTCCTAGAAGACAGCAAATACCG
AAAAATCTACTCTTTACTCTTAAGCCTCGAACAATTGAAGATTTTGGAGTCTATGAATACATACCTGCAG
ACATCTCCATCATCTGTGTTTACTAGTAAATCATTTTGTTCCTTGCAGACCCAGATGGGGTTCACTGTT
TGGTGGGCTTACCCTCACCCATAGGAGATTCAATTATAAGGACAATACAGATCTAATAGAGTTCAAGAC
TCTGAGTGAGGAAGAAATAGAAAAAGTCTGAAAAATATATTTAATATTTCCCTTGCAGAGAAAGCTTGTG
CCCAAACATGGTGATAGATTTTTACTATT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG228808 representing NM_001160170
 Red=Cloning site Green=Tags(s)

MDIEAYLERIGYKKS RNKLDLETLDILQHQIRAVPFENLNIHCGDAMD LGLEAIFDQVVRNRN RGGWCLQ
 VNHLLYWALTTIGFETTMLGGYVYSTPAKKYSTGMIHLLQVTIDGRNYIVDAGFGRSYQMWQPLELISG
 KDQPQVPCVFRLTEENGFWYLDQIRREQYIPNEEFLHSDLLED SKYRKIYSFTLKPRTIEDFESMNTYLQ
 TSPSSVFTSKSFCSLQTPDGVHCLVGF TLT HRRFNKYKDN TDLIEFKTLSEEEIEKVLKNI FNISLQRKLV
 PKHGDRFFTI

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_001160170

ORF Size: 870 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001160170.4
RefSeq Size:	2349 bp
RefSeq ORF:	873 bp
Locus ID:	9
UniProt ID:	P18440
Cytogenetics:	8p22
Protein Pathways:	Caffeine metabolism, Drug metabolism - other enzymes, Metabolic pathways
Gene Summary:	This gene is one of two arylamine N-acetyltransferase (NAT) genes in the human genome, and is orthologous to the mouse and rat Nat2 genes. The enzyme encoded by this gene catalyzes the transfer of an acetyl group from acetyl-CoA to various arylamine and hydrazine substrates. This enzyme helps metabolize drugs and other xenobiotics, and functions in folate catabolism. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011]