

## Product datasheet for **RG228929**

### TCP1 eta (CCT7) (NM\_001166285) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	TCP1 eta (CCT7) (NM_001166285) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	TCP1 eta
Synonyms:	CCTETA; CCTH; NIP7-1; TCP1ETA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RG228929 representing NM\_001166285  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**

ATGATGCCACACCAGTTATCCTATTGAAAGAGGGGACTGATAGCTCCCAAGGCATCCCCAGCTTGTGA  
 GTAACATCAGTGCCTGCCAGGTGATTGCTGAGGCTGTAAGAAGTACCCTGGGTCCCGTGGCATGGACAA  
 GCTTATTGTAGATGGCAGAGGCAAGCAACAATTTCTAATGATGGGGCCACAATTCTGAAACTTCTTGAT  
 GTTGTCATCCTGCAGCAAAGACTTTGGTAGACATTGCCAAATCCCAAGATGCTGAGGTGGGTGATGGCA  
 CCACCTCAGTGACCTTGTGGCTGCAGAGTTTCTGAAGCAGGTGAAACCCATGTGGAGGAAGTTTACA  
 CCCCAGATCATCATTGAGCTTTCCGCACAGCCACCCAGCTGGCAGTTAACAAGATCAAAGAGATTGCT  
 GTGACCGTGAAGAAGGCAGATAAAGTGGAGCAGAGGAAGCTGCTGAAAAGTGTCCATGACCGCTCTGA  
 GCTCCAAGCTGATCTCCAGCAGAAAGCTTTCTTTGCTAAGATGGTGGTGGATGCAGTATGATGCTCGA  
 TGATTTGCTGCAGCTTAAATGATTGGAATCAAGAAGGTACAGGGTGGAGCCCTCGAGGATTCTCAGCTG  
 GTAGCTGGTGTTCATTCAAGAAGACTTTCTCTTACGCTGGGTTTGAAATGCAACCCAAAAAGTACCACA  
 ATCCCAAGATTGCCCTTTTGAATGTCGAGCTCGAGTTGAAAGCTGAGAAAAGACAATGCTGAGATAAGAGT  
 CCACACAGTTGAGGATTACAGGCAATTGTTGATGCTGAGTGGAAACATTCTATGACAAGTTAGAGAAG  
 ATCCATCATTCTGGAGCCAAAGTTGCTTTGTCCAAACTCCCCATTGGGGATGTGGCCACCCAGTACTTTG  
 CTGACAGGGACATGTTCTGTGCTGGCCGAGTACCTGAGGAGGATCTGAAGAGGACAATGATGGCCTGTGG  
 AGGCTCAATCCAGACCAGTGTGAATGCTCTGTGACGAGATGTGCTGGTTCGATGCCAGGTGTTTGAAGAG  
 ACCCAGATTGGAGCGGAGAGGTACAATTTTTTACTGGCTGCCCAAGGCCAAGACATGCACCTTCATTG  
 TCCGTGGCGCGCCGAGCAGTTTATGGAGGAGACAGAGCGGTCCCTGCATGATCCATCATGATCGTCAG  
 GAGGGCCATCAAGAATGATTCAAGTGGTGGTGGTGGCGGGCCATTGAGATGGAACCTCTCCAAGTACCTG  
 CGGGATTACTCAAGGACTATTCCAGGAAAACAGCAGCTGTTGATTGGGGCTTATGCCAAGCCTTGAGAG  
 TTATCCCACGCCAGCTGTGTGACAATGCTGGCTTTGATGCCACAAACATTCTCAACAAGCTGCGGGCTCG  
 GCATGCCAGGGGGTACATGGTATGGAGTAGACATCAACAACGAGGACATTGCTGACAACCTTTGAAGCT  
 TTCGTGTGGGAGCCAGCTATGGTGGGATCAATGCGCTGACAGCAGCCTCTGAGGCTGCGTGCCTGATCG  
 TGTCTGTAGATGAAACCATCAAGAACCCCGCTCGACTGTGGATGCTCCACAGCAGCAGGCCGGGGCCG  
 TGGTCGTGGCCGCCCCAC

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA

**Protein Sequence:**

>RG228929 representing NM\_001166285  
 Red=Cloning site Green=Tags(s)

MMPTPVILLKEGTDSSQIPQLVSNISACQVIAEAVRRTLGPRGMDKLIIVDGRGKATISNDGATILKLLD  
 VVHPAAKTLVDIAKSQDAEVDGTTSVTLAAEFKQVKPYVEEGLHPQIIIRAFRTATQLAVNKIKEIA  
 VTVKKADKVEQRKLEKCAMTALSSKLSISQKFAFFAKMVVDVAVMMLDLDLQKMIKIKKVVQGALEDSQL  
 VAGVAFKKTFSYAGFEMQPKYHNPKIALLNVELELKAEKDNAEIRVHTVEDYQAIVDAEWNILYDKLEK  
 IHHSAGVVL SKLPIGDVATQYFADRFMFCAGRVPEEDLKRMMACGGSIQTSVNL SADVLGRCQVFEE  
 TQIGGERYNFFTGCPKAKTCTF ILRGAEQFMEETERSLHDAIMIVRRAIKNDSVVAGGGAIEMELSKYL  
 RDYSRTIPGKQQLLIGAYAKALEIIPRQLCDNAGFDATNINLKLRRARHAQGGTWTYGVVDINNEDIADNFEA  
 FVWEPAMVRINALTAASEAACLIIVSVDETIKNPRSTVDAPTAAGRGRGRGRPH

**TRTRPLE** - GFP Tag - V

**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**


**ACCN:** NM\_001166285

**ORF Size:** 1632 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001166285.1](#), [NP\\_001159757.1](#)

**RefSeq Size:** 2044 bp

**RefSeq ORF:** 1500 bp


**Locus ID:** 10574

**UniProt ID:** [Q99832](#)

**Cytogenetics:** 2p13.2

**Gene Summary:** This gene encodes a molecular chaperone that is a member of the chaperonin containing TCP1 complex (CCT), also known as the TCP1 ring complex (TRiC). This complex consists of two identical stacked rings, each containing eight different proteins. Unfolded polypeptides enter the central cavity of the complex and are folded in an ATP-dependent manner. The complex folds various proteins, including actin and tubulin. Alternative splicing results in multiple transcript variants. Related pseudogenes have been identified on chromosomes 5 and 6. [provided by RefSeq, Oct 2009]

**Product images:**



Circular map for RG228929