

Product datasheet for **RG228992**

ZBTB20 (NM_001164347) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZBTB20 (NM_001164347) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ZBTB20
Synonyms:	DPZF; HOF; ODA-8S; PRIMS; ZNF288
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG228992 representing NM_001164347
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGACCGAGCGCATTACAGCATCAACCTTCACAACCTTCAGCAATTCGGTCTCGAGACCTCAACGAGC
AGCGCAACCGTGGCCACTTCTGTGACGTAACGGTGCGCATCCACGGGAGCATGCTGCGCGCACACCGCTG
CGTGCTGGCAGCCGCGAGCCCTTCTTCCAGGACAACTGCTGCTTGGCTACAGCGACATCGAGATCCCG
TCGGTGGTGTCAAGTGCAGTCAAGGCTCATTGACTTCAATGTACAGCGCGGTGCTACGGGTCTCGC
AGTCGGAAGCTCTGCAGATCCTCACGGCCGCGCAGCATCTGCAGATCAAAACAGTCATCGACGAGTGCAC
GCGCATCGTGTACAGAACGTGGGCGATGTGTTCCCGGGATCCAGGACTCGGGCCAGGACACGCCCGCG
GGCACTCCCGAGTCAGGCAGTCAAGCCAGAGCAGCGACACGGAGTCGGGCTACCTGCAGAGCCACCCAC
AGCACAGCGTGGACAGGATCTACTCGGCACTCTACGCGTCTCCATGCAGAATGGCAGCGCGGAGCGCTC
TTTTTACAGCGGCGCAGTGGTCAAGCCACGAGACTGCGCTCGGCTGCCCGCGACCACCATGGAA
GACCCAGCTGGATCACACGCATCCATGAGCGCTCGCAGCAGATGGAGCGCTACCTGTCCACCACCCCG
AGACCACGCACTGCCGAAGCAGCCCGGCCCTGTGCGCATCCAGACCCTAGTGGGCAACATCCACATCAA
GCAGGAGATGGAGGACGATTACGACTACTACGGGCAGCAAAGGGTGCAGATCCTGGAACGCAACGAATCC
GAGGAGTGCACGGAAGACACAGACCAGGCCGAGGGCACCAGAGTGGAGCCAAAGGTGAAAGCTTCGACT
CGGGCGTCAAGTCTCCATAGGCACCGAGCCTGACTCGGTGGAGCAGCAGTTTGGGCTGGGGCGGCGCG
GGACAGCCAGGCTGAACCCACCAACCCGAGCAGGCTGCAGAAGCCCCGCTGAGGGTGGTCCGACAGACA
AACCAGCTAGAAACAGGTGCTTCTCTCCGGAGAGAAGCAATGAAGTGGAGATGGACAGCACTGTTATCA
CTGTCAGCAACAGCTCCGACAAGAGCGTCTACAACAGCCTTCGGTCAACACGTCCATCGGGCAGCCATT
GCCAAGTACCCAGCTCTACTTACGCCAGACAGAAACCCTCACAGCAACCTGAGGATGCCTCTGACCTTG
ACCAGCAACACGCAAGTCAATTGGCACAGCTGGCAACACCTACCTGCCAGCCCTTCTCACTACCCAGCCCG
CGGGCAGTGGCCCAAGCCTTCTCTTTCAGCCTGCCACAGCCCTGGCAGGCCAGCAGACCCAGTTTGT
GACAGTGTTCAGCCCGTCTGTGACCTTTACTGCACAGCTGCCAGCGCCACAGCCCTGGCCTCATCC
GCAGGCCACAGCACAGCCAGTGGGCAAGGCGAAAAAAGCCTTATGAGTGCACCTCTGCAACAAGACTT
TCACCGCAAACAGAACTACGTCAAGCACATGTTCTGTACACACAGGTGAGAAGCCCCACCAATGCAGCAT
CTGTTGGCGCTCCTTCTCCTAAAGGATTACCTTATCAAGCACATGGTGACACACACAGGAGTGAGGGCA
TACCAGTGTAGTATCTGCAACAAGCGCTTCACCCAGAAGAGCTCCCTCAACGTGCACATGCGCCTCCACC
GGGAGAGAAGTCTACGAGTGCTACATCTGCAAAAAGAAGTTCTCTCACAAAGACCCTCCTGGAGCGACA
CGTGGCCCTGCACAGTGGCAGCAATGGGACCCCCCTGCAGGCACACCCCGAGTGGCCCGCTGGCCCC
CCAGGCGTGGTGGCTGCACGGAGGGGACCACTTACGTCTGCTCCGTCTGCCAGCAAAGTTTGACCAAA
TCGAGCAGTTCAACGACCACATGAGGATGCATGTGTCTGACGGA

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG228992 representing NM_001164347
 Red=Cloning site Green=Tags(s)


MTERIHSINLHNFNSVLETLNEQRNRGHFCDVTVRIHGSMMLRAHRCVLAAGSPFFQDKLLLGYSIDIEIP
 SVVSVQSVQKLIDFMYSGLRVLSQSEALQILTAASILQIKTVIDECTRIVSQNVGDVFPGIQDSGQDTPR
 GTPESGTSQGSSDTESGYLQSHPQHSVDRIYSALYACSMQNGSGERSFYSGAVVSHHETALGLPRDHME
 DPSWITRIHERSQMERYLSTTPETTHCRKQPRPVRIQTLVGNIIHKQEMEDDYDYYGQQRVQILERNES
 EECTEDTDQAEGTESEPKGESFDSGVSSSIGTEPDSVEQQFGPGAARDSQAEPQPEQAAEAPAEGGPQT
 NQLETGASSPERSNEVEMDSTVITVSNSSDKSVLQQPSVNTSIGQPLPSTQLYLRTETLTSNLRMPLTL
 TSNTQVIGTAGNTYLPALFTTQPAGSGPKPFLFSLPQPLAGQQTQFVTVFQPLSTFTAQLPAPQPLASS
 AGHSTASGQGEKKPYECTLCNKFTAKQNYVKHMFVHTGEKPHQCSICWRSFSLKDYLKHKMVTHTGVRA
 YQCSICNKRFTQKSSLNVHMRLHRGEKSYECYICKKKFHKTLLEHVALHSASNGTPPAGTPPPGARAGP
 PGVVACTEGTTYVCSVCPAKFDQIEQFNDHMRMHVSDG

SGP TRRRLE - GFP Tag - V


Restriction Sites:

SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_001164347

ORF Size: 2004 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001164347.1](#), [NP_001157819.1](#)

RefSeq Size: 3358 bp

RefSeq ORF: 2007 bp

Locus ID: 26137

UniProt ID: [Q9HC78](#)

Cytogenetics: 3q13.31

Protein Families: Transcription Factors

Gene Summary: This gene, which was initially designated as dendritic cell-derived BTB/POZ zinc finger (DPZF), belongs to a family of transcription factors with an N-terminal BTB/POZ domain and a C-terminal DNA-binding zinc finger domain. The BTB/POZ domain is a hydrophobic region of approximately 120 aa which mediates association with other BTB/POZ domain-containing proteins. This gene acts as a transcriptional repressor and plays a role in many processes including neurogenesis, glucose homeostasis, and postnatal growth. Mutations in this gene have been associated with Primrose syndrome as well as the 3q13.31 microdeletion syndrome. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Feb 2017]