

Product datasheet for **RG229129**

OAZ3 (NM_016178) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OAZ3 (NM_016178) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	OAZ3
Synonyms:	AZ3; OAZ-t; TISP15
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG229129 representing NM_016178 Red=Cloning site Blue=ORF Green=Tags(s)

```
CTGCCTTGTAAGAGGTGTCGCCCTCTGTCTACTCCCTTTCTTATATCAAGAGGGGAAAAACACGTAAC
ACCTCTACCCGATCTGGTCACCATACGCCTATTACCTTTACTGTTACAAGTACCGGATCACTCTCCGGGA
GAAGATGCTGCCTCGTTGTTATAAAAGCATCACTTATAAGGAAGAGGAGGACTTGACACTCCAGCCCGT
TCCTGCCTCCAGTGCTCCTGAGTCCCTAGTAGGCCTCCAGGAGGGCAAAGCACCAGCAGGTAACCAC
GACCAGCTAAAGAACTGTATTCGGCTGGGAACCTGACGGTGCTGGCTACTGACCCCTGCTCCACCAGG
ACCCAGTACAGTTAGACTTTCACTTCCGCCTTACCTCCCAGACCTCTGCCATTGGCAGGCCTTCTCTG
TGACCGTCGACTTTCCTGGATATCCCATATCAGGCCTTGGATCAAGGCAACCGGAAAAGTTTGACTGCA
ACCTGGAGTACGTGGAAGAGAAGACAAATGTGGACTCTGTGTTTGTGAACCTCCAGAATGATCGGAACG
ACAGAGGTGCCCTGCTGCGGGCCTTCAGCTACATGGGCTTTGAGGTGGTCAGACCAGATCACCTGCCCT
CCCTCCCTGGACAATGTCATCTTTATGGTGTATCCCCTTCAAAGGGATGTTGGCCACCTGCCAGTGAG
CCTCCT
```


Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_016178
 ORF Size: 453 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_016178.2 , NP_057262.2
RefSeq Size:	901 bp
RefSeq ORF:	709 bp
Locus ID:	51686
UniProt ID:	Q9UMX2
Cytogenetics:	1q21.3
Gene Summary:	The protein encoded by this gene belongs to the ornithine decarboxylase antizyme family, which plays a role in cell growth and proliferation by regulating intracellular polyamine levels. Expression of antizymes requires +1 ribosomal frameshifting, which is enhanced by high levels of polyamines. Antizymes in turn bind to and inhibit ornithine decarboxylase (ODC), the key enzyme in polyamine biosynthesis; thus, completing the auto-regulatory circuit. This gene encodes antizyme 3, the third member of the antizyme family. Like antizymes 1 and 2, antizyme 3 inhibits ODC activity and polyamine uptake; however, it does not stimulate ODC degradation. Also, while antizymes 1 and 2 have broad tissue distribution, expression of antizyme 3 is restricted to haploid germ cells in testis, suggesting a distinct role for this antizyme in spermiogenesis. Antizyme 3 gene knockout studies showed that homozygous mutant male mice were infertile, and indicated the likely role of this antizyme in the formation of a rigid connection between the sperm head and tail during spermatogenesis. Alternatively spliced transcript variants encoding different isoforms, including one resulting from the use of non-AUG (CUG) translation initiation codon, have been found for this gene. [provided by RefSeq, Dec 2014]