

Product datasheet for **RG229574**

HDAC8 (NM_001166448) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: HDAC8 (NM_001166448) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: HDAC8
Synonyms: CDA07; CDLS5; HD8; HDACL1; KDAC8; MRXS6; RPD3; WTS
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG229574 representing NM_001166448
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGGAGGAGCCGGAGGAACCGCGGACAGTGGCAGTCGCTGGTCCCGTTTATATCTATAGTCCCAGT
ATGTCAGTATGTGTACTCCCTGGCCAAGATCCCCAACGGGCCAGTATGGTGCATTCTTTGATTGAAGC
ATATGCACTGCATAAGCAGATGAGAGATGAAGCATCTGGTTTTTGTATCTCAATGATGCTGTCTGGGA
ATATTACGATTGCGACGGAATTTGAGCGTATTCTCTACGTGGATTGGATCTGCACCATGGAGATGGAA
CAGGTGACGTGTCTGATGTTGGCTAGGGAAGGACGGTACTACAGTGTAATGTGCCCATTCAGGATGG
CATACAAGATGAAAAATATTACCAGATCTGTGAAAGGTACGAACCTCTGCCCAAATCCAGGCCTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG229574 representing NM_001166448
Red=Cloning site Green=Tags(s)

MEEPEEPADSGQSLVPVYIYSPEYVSMCDSLAKIPKRASMVHSLIEAYALHKQMRDEASGFCYLNDVGLG
 ILRLRRKFERILYVDLDLHHGDGTGDVSDVGLGKGRYYSVNVPIQDGIQDEKYYQICERYEPPAPNPGL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001166448

ORF Size: 417 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001166448.2
RefSeq Size:	2071 bp
RefSeq ORF:	420 bp
Locus ID:	55869
UniProt ID:	Q9BY41
Cytogenetics:	Xq13.1
Protein Families:	Druggable Genome, Transcription Factors
Gene Summary:	Histones play a critical role in transcriptional regulation, cell cycle progression, and developmental events. Histone acetylation/deacetylation alters chromosome structure and affects transcription factor access to DNA. The protein encoded by this gene belongs to class I of the histone deacetylase family. It catalyzes the deacetylation of lysine residues in the histone N-terminal tails and represses transcription in large multiprotein complexes with transcriptional co-repressors. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2009]