

Product datasheet for **RG229600**

PAFAH1B2 (NM_001184747) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: PAFAH1B2 (NM_001184747) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: PAFAH1B2
Synonyms: HEL-S-303
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG229600 representing NM_001184747
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGCCAAGGAGACTCAAACCCAGCAGCTATTCGCATGCAGCAGAAGATATCAAGGAGATGACCGAT
GGATGTCTCAGCACAACAGATTTGTTTTGGACTGTAAAGACAAAGAGCCTGATGTACTGTTCTGGGAGA
CTCCATGGTGCAGTTAATGCAGCAATATGAGATATGGCGAGAGCTTTTTCCCACTTCATGCAGTGAAT
TTTGAATTGGGGGAGATACAACAAGACATGTTTTGTGGAGACTAAAGAATGGAGAAGTGGAGAATTA
AGCCTAAGGTCATTGTTGTCTGGGTAGGAACAAATAACCACGAAAATACAGCAGAAGAAGTAGCAGGTGG
GATCGAGGCCATTGTACAACTTATCAACACAAGGCAGCCACAGGCCAAAATCATTGTATTGATTATTAC
TGCAAGATGAGCAGGATTACCATGAAAGAAAAGTCCAGATGGAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG229600 representing NM_001184747
Red=Cloning site Green=Tags(s)

MSQGDSNPAAIPHAIEDIQGDDRWMSQHNRFLVLDCKDKEPDVLFVGDMSVQLMQQYEWRELFSPLHALN
FGIGDTRHVLWRLKNGELENIKPKVIVVWGTNNHENTAEVAGGIEAIVQLINTRQPQAKIIVLIY
WQDEQDYHERKVQMD

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001184747

ORF Size: 465 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001184747.2</u>
RefSeq Size:	1070 bp
RefSeq ORF:	468 bp
Locus ID:	5049
UniProt ID:	<u>P68402</u>
Cytogenetics:	11q23.3
Protein Pathways:	Ether lipid metabolism, Metabolic pathways
Gene Summary:	Platelet-activating factor acetylhydrolase (PAFAH) inactivates platelet-activating factor (PAF) into acetate and LYSO-PAF. This gene encodes the beta subunit of PAFAH, the other subunits are alpha and gamma. Multiple alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Jan 2014]