

Product datasheet for **RG229634**

PLEKHG7 (NM_001178097) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: PLEKHG7 (NM_001178097) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: PLEKHG7
Synonyms: C12orf74
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG229634 representing NM_001178097
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGAAAACAGAGTCATTCTGTCCAGAGGTGCCACCCCAAGACTGTGGAGCCTCTCCTCGGCCCTCGC
TGAGGAGCCTGCCAAAGAACCAGGGGAGTCTCCTCCAGTTTGACCGCAAGCCCCAGGCCGATCTCCAC
CTCGCCACTTTGAGGAGATTAAGGACCCGTGGCTGTGGACAAGGCAGGATGCCTGGCAGGTGACCACC
TGGGAAGCTGGGAGCTCCTGTGGCTTCCCATGTTACCTTTGAAGAGCCTGCCAGGAAGCCCAAGG
ATTCTTACACTTGCTGTACCCTTGAGACTCCACTCAAGATTGACCTCTGAACCTGAAAGGCCCTGAA
TGCAGCTGACTACTGGAGCCCCAAACCCGGCCCACTGACAAGTATCTCCCTCCTGAGCTTCAGCCTGTC
AATGAAGGGTCCCTTACCAGGCCTCTCTCGGCAGCAAGAAGGCCACTTCTGCCAGCCCCACCCTAC
GACACCCTAGTCTCAGGTAACACAGCTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG229634 representing NM_001178097
Red=Cloning site Green=Tags(s)

MEKTESFCPEVPPQDCGASPRPSLRSLPKNQGSLLQFDRQAPGRISTSPTLRRLRTRGCGTRQDAWQVTT
WGSWGAPVGFPCYLKSLPGSPKDSHLLSPLRLHSRLTSEPERALNAADSLEPQTRPTDKYLPPELQPV
NEGSLHQASLRQQEGHFLPSPTLRHPSQVTQL

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001178097

ORF Size: 519 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001178097.3</u>
RefSeq Size:	2259 bp
RefSeq ORF:	522 bp
Locus ID:	338809
UniProt ID:	<u>Q32Q52</u>
Cytogenetics:	12q22