

## Product datasheet for **RG229708**

### Cadherin like 23 (CDH23) (NM\_001171936) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Cadherin like 23 (CDH23) (NM_001171936) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CDH23
Synonyms:	CDHR23; PITA5; USH1D
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG229708 representing NM_001171936 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCGGTCCTGGTTCCAGCAGGATCCTATGGTGGGAGCATGCACCACAGGCACCAGGGCCTCACACCCCA  
AAGCCAACCTGTGTGGCTGGATCCCTTCTGTCCGAACCTGGAGCTGGCCGCCAGGGCGGAGCATGAGGA  
TGACCTACCGGAGAACCTGAGTGAGATCGCCGACCTGTGGAACAGCCCCACGCGCACCCATGGAACTTTT  
GGCGGTGAGCCAGCAGCTGTCAAGCCTGATGATGACCGATACCTGCGGGCTGCCATCCAGGAGTATGACA  
ACATTGCCAAGCTGGGCCAGATCATTCTGTGAGGGCCAATCAAGCTGATACAGACTGAGCTGGACGAGGA  
GCCAGGAGACCACAGCCAGGGCAGGGTAGCCTGCGCTTCCGCCACAAGCCACCAGTGGAGCTCAAGGGG  
CCCGATGGGATCCATGTGGTGCACGGCAGCACGGGCACGCTGCTGGCCACCGACCTCAACAGCCTGCCCCG  
AGGAAGACCAGAAGGGCCTGGGCCGCTCGCTGGAGACGCTGACCCTGCCGAGGCCACTGCCTTCGAGCG  
CAACGCCCGCACAGAATCCGCCAAATCCACACCCTGCACAACTTCGCGACGTGATCATGGAGACCCCC  
CTGGAGATCACAGAGCTG

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG229708 representing NM\_001171936  
Red=Cloning site Green=Tags(s)

MRSWFQQDPMVGACTTGRASHPKANPVWLDPFRCRNLELAAQAEHDDLLENLSEIADLWNSPTRTHGTG  
 GREPAAVKPD DRYLRAAIQEYDNI AKLGQI IREGPIKLIQTELDEEPGDHSPGQGS LFRFRHKPPVELKG  
 PDGIHVHVGSTG TLLATDLNLSLPEEDQKGLGRSLET LTAAEATAFERNARTESAKSTPLHLKLRD VIMETP  
 LEITEL


TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001171936

**ORF Size:** 648 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001171936.1</a> , <a href="#">NP_001165407.1</a>
<b>RefSeq Size:</b>	1739 bp
<b>RefSeq ORF:</b>	651 bp
<b>Locus ID:</b>	64072
<b>UniProt ID:</b>	<a href="#">Q9H251</a>
<b>Cytogenetics:</b>	10q22.1
<b>Protein Families:</b>	Transmembrane
<b>Gene Summary:</b>	This gene is a member of the cadherin superfamily, whose genes encode calcium dependent cell-cell adhesion glycoproteins. The encoded protein is thought to be involved in stereocilia organization and hair bundle formation. The gene is located in a region containing the human deafness loci DFNB12 and USH1D. Usher syndrome 1D and nonsyndromic autosomal recessive deafness DFNB12 are caused by allelic mutations of this cadherin-like gene. Upregulation of this gene may also be associated with breast cancer. Alternative splice variants encoding different isoforms have been described. [provided by RefSeq, May 2013]