

## Product datasheet for **RG229743**

### Tissue Factor (F3) (NM\_001178096) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tissue Factor (F3) (NM\_001178096) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** F3  
**Synonyms:** CD142; TF; TFA  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG229743 representing NM\_001178096  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGAGACCCCTGCCTGGCCCCGGGTCCCAGCCCGAGACCGCCGTCGCTCGGACGCTCCTGCTCGGCT  
GGGTCTTCGCCCAGGTGGCCGGCGCTTCAGGCACTACAAATACTGTGGCAGCATATAATTTAACTGGAA  
ATCAACTAATTTCAAGACAATTTGGAGTGGGAACCCAAACCCGTCATCAAGTCTACACTGTTCAAATA  
AGCACTAAGTCAAGGATTTGAAAAGCAAATGCTTTTACACAACAGACACAGAGTGTGACCTCACCGACG  
AGATTGTGAAGGATGTGAAGCAGACGACTTGGCAGGGTCTTCTCTACCCGGCAGGGAATGTGGAGAG  
CACCGGTTCTGCTGGGGAGCCTCTGTATGAGAATCCCCAGAGTTCACACCTTACCTGGAGACAAACCTC  
GGACAGCCAACAATTCAGAGTTTGAACAGGTGGGAACAAAAGTGAATGTGACCGTAGAAGATGAACGGA  
CTTTAGTCAGAAGGAACAACACTTTCTAAGCCTCCGGGATGTTTTTGGCAAGGACTTAATTTATACACT  
TTATTATTGAAATCTTCAAGTTCAGGAAAGAAATATTCTACATCATTGGAGCTGTGGTATTTGTGGTCA  
TCATCCTGTATCATCCTGGCTATATCTCTACACAAGTGTAGAAAGGCAGGAGTGGGCGAGAGCTGGAA  
GGAGAACTCCCCAC

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG229743 representing NM\_001178096  
 Red=Cloning site Green=Tags(s)

METPAWPRVPRPETAVARTLLLGWVFAQVAGASGTTNTVAAYNL TWKSTNFKTILEWEPKPVNQVYTVQI  
 STKSGDWKSKCFYTTDTECDL TDEIVKDVKQTYLARVFSYPAGNVESTGSAGEPL YENSPEFTPYLETNL  
 GQPTIQSFEQVGTKVNVTVEDERTLVRRNNTFLSLRDVFGKDLIYTLYYWKSSSSGKKYSTLELWYLWS  
 SSLSSSWLYLYTSVERQEWGRGRRTPH

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM\_001178096

ORF Size: 714 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001178096.1</a> , <a href="#">NP_001171567.1</a>
<b>RefSeq Size:</b>	2233 bp
<b>RefSeq ORF:</b>	717 bp
<b>Locus ID:</b>	2152
<b>UniProt ID:</b>	<a href="#">P13726</a>
<b>Cytogenetics:</b>	1p21.3
<b>Protein Families:</b>	Druggable Genome, Transmembrane
<b>Protein Pathways:</b>	Complement and coagulation cascades
<b>Gene Summary:</b>	This gene encodes coagulation factor III which is a cell surface glycoprotein. This factor enables cells to initiate the blood coagulation cascades, and it functions as the high-affinity receptor for the coagulation factor VII. The resulting complex provides a catalytic event that is responsible for initiation of the coagulation protease cascades by specific limited proteolysis. Unlike the other cofactors of these protease cascades, which circulate as nonfunctional precursors, this factor is a potent initiator that is fully functional when expressed on cell surfaces, for example, on monocytes. There are 3 distinct domains of this factor: extracellular, transmembrane, and cytoplasmic. Platelets and monocytes have been shown to express this coagulation factor under procoagulatory and proinflammatory stimuli, and a major role in HIV-associated coagulopathy has been described. Platelet-dependent monocyte expression of coagulation factor III has been described to be associated with Coronavirus Disease 2019 (COVID-19) severity and mortality. This protein is the only one in the coagulation pathway for which a congenital deficiency has not been described. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Aug 2020]