

Product datasheet for **RG229861**

Lunatic Fringe (LFNG) (NM_001166355) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lunatic Fringe (LFNG) (NM_001166355) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	LFNG
Synonyms:	SCDO3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG229861 representing NM_001166355 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGATGAACAGACAGGAAGGCTCAGGCTGGACACGTATTGTATGAGTGCCAAGCAGATCTGGGCATGGA
GCAAATGCTCAGGAAGGCTGTGGGATGAGCACATGAAATGGATGGAAGGATGGACGGACAGATGGACAGA
TGGATGGATGGATGGATGGATGAGTGGAGCCCAACACCAGCTCTCAGGTCTACGGAGGTGGCCTC
TCTCAGCAGACGTTTCATCTTCACTGACGGGGAAGATGAGGCCCTGGCCAGGCACACGGGCAACGTGGTCA
TCACAAACTGCTCGGCCGCCACAGCCGCCAGGCGCTGTCTGCAAGATGGCCGTGGAGTATGACCGCTT
CATCGAGTCCGGCAGGAAGTGGTTCTGCCACGTGGACGATGACAACCTACGTCAACCTGCGGGCCCTGCTG
CGGCTGTGGCCAGCTACCCGCACACGCGGGACGTCTACGTGGCAAGCCAGCCTGGACAGGCCATCC
AGGCCATGGAGCGGGTCAGCGAGAACAAGGTGCGTCTGTCCACTTCTGGTTTGCCACGGGCGGCGCTGG
CTTCTGCATCAGCCGTGGGCTGGCTCTGAAGATGAGCCGTGGGCCAGCGGGGTCACCTCATGAATACG
GCTGAGCGGATCCGGCTGCCTGATGACTGCACCATCGGCTACATCGTGGAGGCCCTGCTGGGTGTGCCCC
TCATCCGACGCGCCTCTTCCACTCCACCTGGAGAACCTGCAGCAGGTGCCACCTCGGAGCTCCACGA
GCAGGTGACGCTGAGCTACGGTATGTTTTGAAAACAAGCGGAACGCCGTCCACGTGAAGGGGCCCTTCTCG
GTGGAGGCCGACCCATCCAGTTCGGCTCCATCCACTGCCACCTGTACCCGGACACACCTGGTGTCCCC
GCACTGCCATCTTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG229861 representing NM_001166355
 Red=Cloning site Green=Tags(s)


MDEQTGRLRLDITYCMSAKQIWAWSKCSGRLWDEHMKWMEGWTDRWTDGWMGMDDEWSPTPALRSYGGGL
 SQQTFIFTDGEDEALARHTGNVITNCSAAHSRQALSCKMAVEYDRFIESGRKWFCHVDDDNVNLRAL
 RLLASYPHTRDVYVGKPSLDRPIQAMERVS ENKVRPVHFWFATGGAGFCISRGLALKMSPWASGGHFMNT
 AERIRLPDDCTIGYIVEALLGVPLIRSLFHSHLENLQQVPTSELHEQVTL SYGMFENKRNAVHVKGPF S
 VEADPSRFRSIIHCHLYPDTWPCRTAIF

TRTRPLE - GFP Tag - V


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001166355

ORF Size: 924 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001166355.1 , NP_001159827.1
RefSeq Size:	2280 bp
RefSeq ORF:	927 bp
Locus ID:	3955
UniProt ID:	Q8NES3
Cytogenetics:	7p22.3
Protein Families:	Transmembrane
Protein Pathways:	Notch signaling pathway
Gene Summary:	<p>This gene is a member of the glycosyltransferase 31 gene family. Members of this gene family, which also includes the MFNG (GeneID: 4242) and RFNG (GeneID: 5986) genes, encode evolutionarily conserved glycosyltransferases that act in the Notch signaling pathway to define boundaries during embryonic development. While their genomic structure is distinct from other glycosyltransferases, these proteins have a fucose-specific beta-1,3-N-acetylglucosaminyltransferase activity that leads to elongation of O-linked fucose residues on Notch, which alters Notch signaling. The protein encoded by this gene is predicted to be a single-pass type II Golgi membrane protein but it may also be secreted and proteolytically processed like the related proteins in mouse and Drosophila (PMID: 9187150). Mutations in this gene have been associated with autosomal recessive spondylocostal dysostosis 3. [provided by RefSeq, May 2018]</p>