

Product datasheet for **RG229873**

MYD88 (NM_001172567) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MYD88 (NM_001172567) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MYD88
Synonyms:	IMD68; MYD88D
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG229873 representing NM_001172567 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCGACCCGACCGCGCTGAGGCTCCAGGACCGCCGCCATGGCTGCAGGAGGTCCCGCGCGGGTCTG
CGGCCCGGTCTCCTCCACATCCTCCCTCCCTGGCTGCTCAACATGCGAGTGCGGCCCGCCTGTC
TCTGTTCTTGAACGTGCGGACACAGGTGGCGGCCGACTGGACCGCTGGCGGAGGAGATGGACTTTGAG
TACTTGGAGATCCGGCAACTGGAGACACAAGCGGACCCCACTGGCAGGCTGCTGGACGCTGGCAGGGAC
GCCCTGGCGCCTCTGTAGCCGACTGCTCGAGCTGCTTACCAAGCTGGGCCGCGACGACGTGCTGCTGGA
GCTGGGACCCAGCATTGAGGAGGATTGCCAAAAGTATATCTTGAAGCAGCAGCAGGAGGAGCTGAGAAG
CCTTTACAGGTGGCCGCTGTAGACAGCAGTGTCCACGGACAGCAGAGCTGGCGGGCATCACCACACTTG
ATGACCCCTGGGGCATAATGCCTGAGCGTTTCGATGCCTTCATCTGCTATTGCCCCAGCGACATCCAGTT
TGTGCAGGAGATGATCCGGCAACTGGAACAGACAACTATCGACTGAAGTTGTGTGTCTGACCGCGAT
GTCCTGCCTGGCACCTGTGTCTGGTCTATTGCTAGTGAGCTCATCGAAAAGAGGTTGGCTAGAAGCCAC
GGGGTGGTGCCCGGATGGTGGTGGTGTCTCTGATGATTACCTGCAGAGCAAGGAATGTGACTTCCA
GACCAAAATTTGCACTCAGCCTCTCTCCAGGTGCCCATCAGAAGCGACTGATCCCCATCAAGTACAAGGCA
ATGAAGAAAGAGTTCCCCAGCATCTGAGGTTTCATCACTGTCTGGACTACACCAACCCCTGCACCAAAAT
CTTGTTCTGACTCGCCTTGCCAAGCCTTGTCCCTGCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG229873 representing NM_001172567
 Red=Cloning site Green=Tags(s)

MRPDRAEAPGPPAMAAGGPGAGSAAVVSSTSSPLAALNMRVRRRLSLFLNVRTQVAADWTALAEEMDFE
 YLEIRQLETQADPTGRLLDAWQGRPGASVGRLLLELLTKLGRDDVLELGPSTIEEDCQKYILKQQQEEAEK
 PLQVAAVDSSVPRTAELAGITTLDDPLGHMPERFDFICYCPSDIQFVQEMIRQLEQTNYRLKLCVSDRD
 VLPGTCVWSIASELIEKRLARRPRGGCRMVVVVSDDYLSQKEDCFQTKFALSLSPGAHQKRLIPIKYKA
 MKKEFPSILRFITVCDYTNPCTKSWFWTRLAKALSPL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001172567

ORF Size: 951 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001172567.1](#), [NP_001166038.1](#)

RefSeq Size: 2886 bp

RefSeq ORF: 915 bp

Locus ID: 4615


Cytogenetics: 3p22.2

Protein Families: Druggable Genome

Protein Pathways: Apoptosis, Toll-like receptor signaling pathway

Gene Summary: This gene encodes a cytosolic adapter protein that plays a central role in the innate and adaptive immune response. This protein functions as an essential signal transducer in the interleukin-1 and Toll-like receptor signaling pathways. These pathways regulate that activation of numerous proinflammatory genes. The encoded protein consists of an N-terminal death domain and a C-terminal Toll-interleukin1 receptor domain. Patients with defects in this gene have an increased susceptibility to pyogenic bacterial infections. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Feb 2010]

Product images:



Circular map for RG229873