

Product datasheet for **RG229954**

CD244 (NM_001166663) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: CD244 (NM_001166663) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: CD244
Synonyms: 2B4; NAIL; NKR2B4; Nmrk; SLAMF4
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG229954 representing NM_001166663
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCTGGGGCAAGTGGTCACCCTCATACTCCTCCTGCTCCTCAAGGTGTATCAGGGCAAAGGATGCCAGG
 GATCAGCTGACCATGTGGTTAGCATCTCGGAGTGCCCTTTCAGTTACAACCAAACAGCATAACAGACGAA
 GGTTGACAGCATTGCATGGAAGAAGTTGCTGCCCTCACAAATGGATTTTCATCACATATTGAAGTGGGAG
 AATGGCTCTTTCCTTCCAATACTTCCAATGATAGATTCAAGTTTATAGTCAAGAACTTGAGTCTTCTCA
 TCAGGCAGCTCAGCAGCAGGACAGTGGCCTCTACTGCCTGGAGGTCACCAGTATATCTGGAAAAGTTCA
 GACAGCCACGTTCCAGGTTTTGTATTTGAATCTCTGCTTCCAGATAAAGTTGAGAAACCCCGCCTACAG
 GGGCAGGGGAAGATCCTGGACAGAGGGAGATGCCAAGTGGCTCTGTCTTGTCTCCAGGGATGGCA
 ATGTGTCTATGCTTGGTACAGAGGGAGCAAGCTGATCCAGACAGCAGGGAACTCACCTACCTGGACGA
 GGAGGTTGACATTAATGGCACTCACACATACCTGCAATGTCAGCAATCCTGTTAGCTGGGAAAGCCAC
 ACCCTGAATCTCACTCAGGACTGTCAGAATGCCATCAGGAATTCAGATTTTGGCCGTTTTTGGTGTCA
 TCGTGATTCTAAGCGCACTGTTCTTGGCACCCCTGCCTGCTTCTGTGTGGAGGAGAAAGAGGAAGGA
 GAAGCAGTCAGAGACCAGTCCCAAGGAATTTTGAACAATTTACGAAGATGTCAAGGATCTGAAAACCAGG
 AGAAATCAGAGCAGGAGCAGACTTTTCTGGAGGGGGAGCACCATCTACTCTATGATCCAGTCCCAGT
 CTCTGCTCCACGTCACAAGAAGCTGCATATACATTATTCATTAATTCAGCCTTCCAGGAAGTCTGG
 ATCCAGGAAGAGGAACCCAGCCCTTCTTCAATAGCACTATCTATGAAGTGATTGGAAAGAGTCAACCT
 AAAGCCAGAACCCTGCTCGATTGAGCCGCAAAGAGCTGGAGAACCTTGTGTTTATTCC

ACGGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG229954 representing NM_001166663
Red=Cloning site Green=Tags(s)

MLGQVVTLLLLLLLVYQKGKCGSADHVVSISGVPLQLQPNSIQTKVDSIAWKLLPSQNGFHILKWE
 NGLSPNNTSNDRF SFIVKNLSLLIKAAQQQDGLYCLEVTSISGKVTATFQVVFESLLPDKVEKPRLO
 GQGKILDRGRCQVALSCLVSRDGNVSYAWYRGSKLIQTAGNLTYLDEEVDINGHTYTCNVSNPVSWEH
 TLNLTQDCQNAHQEFRFWPFLVIVILSALFLGLTACFCVWRRKRKEKQSETSPKEFLTIVYEDVKDLKTR
 RNHEQEQTFFPGGSTIYSMIQSQSSAPTSQEPAYTLYSLIQPSRKSRSRKRNHSPSFNSTIYEVIGKSQP
 KAQNPARLSRKELENFDVYS


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001166663

ORF Size:	1110 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001166663.2
RefSeq Size:	2543 bp
RefSeq ORF:	1113 bp
Locus ID:	51744
UniProt ID:	Q9BZW8
Cytogenetics:	1q23.3
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane
Protein Pathways:	Natural killer cell mediated cytotoxicity
Gene Summary:	<p>This gene encodes a cell surface receptor expressed on natural killer (NK) cells (and some T cells) that mediate non-major histocompatibility complex (MHC) restricted killing. The interaction between NK-cell and target cells via this receptor is thought to modulate NK-cell cytolytic activity. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Oct 2009]</p>