

Product datasheet for **RG230047**

Pregnancy Specific Glycoprotein 1 (PSG1) (NM_001184825) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pregnancy Specific Glycoprotein 1 (PSG1) (NM_001184825) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	PSG1
Synonyms:	B1G1; CD66f; DHFRP2; FL-NCA-1/2; PBG1; PS-beta-C/D; PS-beta-G-1; PSBG-1; PSBG1; PSG95; PSGGA; PSGIIA; SP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG230047 representing NM_001184825
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGGAACCCCTCTCAGCCCTCCCTGCACACAGCGCATCAAATGGAAGGGGCTCCTGCTCACAGCATCAC
 TTTTAAACTTCTGGAACCTGCCACCCTGCCAAGTCACGATTGAAGCCGAGCCAACAAAGTTCCGA
 GGGGAAGGATGTTCTTCTACTTGTCCACAATTTGCCCCAGAATCTTACCGCTACATCTGGTACAAAGG
 CAAATGAGGGACCTCTACCATTACATTACATCATATGTAGTAGACGGTAAATAATTATATATGGCCTG
 CATATAGTGGACGAGAAACAGCATATTTCAATGCATCCCTGCTGATCCAGAATGTCACCCGGGAGGACGC
 AGGATCTACACCTTACACATCATAAAGGGAGATGATGGGACTAGAGGAGTAAGTGGACGTTTCACCTTC
 ACCTTACACCTGGAGACTCCTAAGCCCTCCATCTCCAGCAGCAACTTAAATCCCAGGGAGACCATGGAGG
 CTGTGAGCTTAACCTGTGACCCTGAGACTCCAGACGCAAGCTACCTGTGGTGGATGAATGGTCAGAGCCT
 CCCTATGACTCACAGCTTGAAGCTGTCCGAAACCAACAGGACCCTCTTTCTATTGGGTGTCACAAAGTAT
 ACTGCAGGACCCTATGAATGTGAAATACGGAACCCAGTGAGTGCCAGCCGAGTACCCAGTACCCCTGA
 ATCTCCTCCCGAAGCTGCCAAGCCCTACATCACCATCAACAACCTAAACCCAGGGAGAATAAGGATGT
 CTTAAACTTACCTGTGAACCTAAGAGTGAGAACTACACCTACATTTGGTGGCTAAATGGTCAGAGCCTC
 CCGGTCAGTCCCAGGGTAAAGCGACCCATTGAAAACAGGATCCTCATTCTACCCAGTGTCCAGAGAAATG
 AAACAGGACCCTATCAATGTGAAATACGGGACCGATATGGTGGCATCCGAGTGACCCAGTACCCCTGAA
 TGTCTCTATGGTCCAGACCTCCCAGAATTTACCCTTACCTATTACCGTTCAGGAGAAGTCCCTC
 TACTTGTCTGTCTGCGGACTTAACCCACGGCACAGTATTCTTGGACAATTAATGAAAAGTTTCAGC
 TACCAGGACAAAAGCTCTTTATCCGCCATATTACTACAAAGCATAGCGGGCTCTATGTTTGTCTGTTCG
 TAACTCAGCCACTGGCAAGGAAAGCTCCAATCCATGACAGTCGAAGTCTCTGACTGGACAGTTC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG230047 representing NM_001184825
 Red=Cloning site Green=Tags(s)

MGTL SAPPCTQRIKWKGLLLTASLLNFWNLPTTAQVTIEAEPK VSEKDVLLLVHNL PQLNGYIWKY
 QMRDLYHYITSYVVDGEIIIYGPAYSGRETAYSNASLLIQNVTR EDAGSYTLHIKGGDGRGVTGRFTF
 TLHLETPKPSISSNLPRETMEAVSLTCDPETPDASYLWWMNGQSLPMT HSLKLSETRTLFLLGVTKY
 TAGPYECEIRNPVSASRSDPVTLNLLPKLPKYITINNLNPRENKDVLNFTCEPKSENYTYIWWLNGQSL
 PVSPRVKRP IENRILILPSVTRNETGPYQCEIRDYGGIRSDPVTLNVL YGPDLPRIYPSFTYYRSGEVL
 YLSCSADSNPPAQYSWTINEKFQLPGQLFIRHITTKHSGLYVCSVRNSATGKESKSMTEVSDWTVP

TRTRPLE – GFP Tag – V


Restriction Sites:

Sgfl-Mlul

Cloning Scheme:



Plasmid Map:



ACCN: NM_001184825

ORF Size: 1257 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001184825.2](#)

RefSeq Size: 2071 bp

RefSeq ORF: 1260 bp

Locus ID: 5669

UniProt ID: [P11464](#)

Cytogenetics: 19q13.2

Protein Families: Secreted Protein

Gene Summary: The human placenta is a multihormonal endocrine organ that produces hormones, enzymes, and other molecules that support fetal survival and development. Pregnancy-specific beta-1-glycoprotein (PSBG, PSG) is a major product of the syncytiotrophoblast, reaching concentrations of 100 to 290 mg/l at term in the serum of pregnant women (Horne et al., 1976 [PubMed 971765]). PSG is a member of the immunoglobulin (Ig) superfamily (Watanabe and Chou, 1988 [PubMed 3257488]; Streydio et al., 1988 [PubMed 3260773]).[supplied by OMIM, Oct 2009]