

Product datasheet for **RG230084**

MPP1 (NM_001166462) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MPP1 (NM_001166462) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MPP1
Synonyms:	AAG12; DXS552E; EMP55; MRG1; PEMP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG230084 representing NM_001166462 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGTCCTGGGCTGTATCGCATCCATTGAATACTGTGACCGAGGACATGTACACCAACGGGTCTCCTG
CCCCAGGTAGCCCTGCCAGGTCAAGGGACAGGAGGTGCGGAAAGTGCAGCTCATACAGTTTGAGAAGGT
CACAGAAGAGCCCATGGGAATCACGCTGAAGCTGAATGAAAAACAGTCCTGTACGGTGGCCAGAATTCTT
CATGGTGGCATGATCCATAGACAAGGCTCCCTCACGTGGGGATGAGATCCTAGAAATCAATGGCACAA
ATGTGACAAATCATTAGTGGATCAGCTGCAGAAGGCGATGAAAGAAACCAAGGAATGATCTCATTAAA
AGTAATCCCAACCAGCAAAGCCGCTTCTCCTGACTACAGATGTTTCATGAGAGCGCAGTTTGACTATGAT
CCCCAAAAGGACAATCTGATCCCTTGAAGGAGGCGGACTGAAGTTTGCTACTGGGACATTATCCAGA
TTATCAACAAGGATGACAGCAATTGGTGGCAGGGACGGGTGGAAGGCTCCTCAAGGAGTCAGCAGGATT
GATCCCTTCCCCTGAGCTGCAGGAATGGCGAGTGGCAAGTATGGCTCAGTCAGCTCCTAGCGAAGCCCCG
AGCTGCAGTCCCTTTGGGAAGAAGAAGTACAAAGACAAATATCTGGCCAAGCACAGCTCGATTTTGT
ATCAGTTGGATGTTGTTTCTACGAGGAAGTCGTTCCGCTCCCTGCATTCAGAGGAAGACCCTGGTGT
GATCGGAGCCAGTGGGGTGGTTCGACGCCACATTAAGAATGCCCTGCTCAGCCAGAATCCGGAGAAGTTT
GTGTACCTGTCCATATACAACACGGCCGCAAGGAAGAGTGAAGGAAGTGGGAAGGATACCACCTTAA
TCTCAACGGAGGAGATGACGAGGAACATCTCTGCCAATGAGTTCTTGAGTTTGGCAGCTACCAAGGCAA
CATGTTTGGCACCAAAATTTGAAACAGTGCACCAGATCCATAAGCAGAACAAGATTGCCATCCTTGACATT
GAGCCCCAGACCTGAAAATTGTTTCGGACAGCAGAACTTTCGCTTTCATTGTGTTTCATTGCACCTACTG
ACCAGGGCACTCAGACAGAAGCCCTGCAGCAGCTGCAGAAGGACTCTGAGGCCATCCGCAGCCAGTACGC
TCACTACTTTGACCTCTCACTGGTCAATAATGGTGTGATGAAACCCCTTAAGAAATTACAAGAAGCCTTC
GACCAAGCGTGCAGTTCTCCACAGTGGTGCCTGTCTCCTGGTTTAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG230084 representing NM_001166462
Red=Cloning site Green=Tags(s)

MESWAVSHPLNTVTEDMYTNNGSPAPGSPAQVKGQEVKRVRLIQFEKVTEEPMGITLKLNEKQSCSTVARIL
 HGGMIHRQGLSHVGDEILEINGTNVTNHSVDQLQKAMKETKGMISLKVIPNQSRPALQFMRAQFDYD
 PKKDNLIPCKEAGLK FATGDIIQIINKDDSNWWQGRVEGSSKESAGLIPSPELQEWVRSMAQSAPEAP
 SCSPFSGKKKKYKDKYLAKHSSIFDQLDVVSYEEVVRLPAFKRRTLVLIGASGVGRSHIKNALLSQNPEKF
 VYPVYTRPPRKSEEDGKEYHFISTEEMTRNISANEFLEFGSYQGNMFGTKFETVHQIHKQNKIAILDI
 EPQTLKIVRTAELSPFIVFIAPTDQGTQTEALQQLQKDSEAIRSQYAHYFDLSLVNNGVDETLKKLQEAF
 DQACSSPQWVPVSWVY

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001166462

ORF Size: 1308 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001166462.1](#), [NP_001159934.1](#)

RefSeq Size: 2216 bp

RefSeq ORF: 1311 bp

Locus ID: 4354


UniProt ID: [Q00013](#)

Cytogenetics: Xq28

Protein Families: Druggable Genome

Gene Summary: This gene encodes the prototype of the membrane-associated guanylate kinase (MAGUK) family proteins. MAGUKs interact with the cytoskeleton and regulate cell proliferation, signaling pathways, and intercellular junctions. The encoded protein is an extensively palmitoylated membrane phosphoprotein containing a PDZ domain, a Src homology 3 (SH3) motif, and a guanylate kinase domain. This gene product interacts with various cytoskeletal proteins and cell junctional proteins in different tissue and cell types, and may be involved in the regulation of cell shape, hair cell development, neural patterning of the retina, and apico-basal polarity and tumor suppression pathways in non-erythroid cells. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2009]

Product images:



Circular map for RG230084