

## Product datasheet for **RG230141**

### TFEB (NM\_001167827) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	TFEB (NM_001167827) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	TFEB
Synonyms:	ALPHATFEB; BHLHE35; TCFEB
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RG230141 representing NM\_001167827  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCGCGATCGCC

ATGACAGCAAGCTCAGGCTGGGAGCCAGCGCCGCGAGCCACCATGGCGTCACGCATAGGGTTGCGCATGC  
 AGCTCATGCGGGAGCAGGCGCAGCAGGAGGAGCAGCGGGAGCGCATGCAGCAACAGGCTGTGCATGATTA  
 CATGCAGCAGCAGCAGCAGCAGCAACAGCAGCAGCTCGGAGGGCCGCCACCCCGCCATCAATACCCCC  
 GTCCACTTCCAGTCGCCACCACCTGTGCCTGGGGAGGTGTTGAAGGTGCAGTCCTACCTGGAGAATCCCA  
 CATCTACCATCTGCAGCAGTCGCAGCATCAGAAGGTGCGGGAGTACCTGTCCGAGACCTATGGGAACA  
 GTTTGCTGCCACATCAGCCCAGCCAGGGCTCTCGAAACCCACCAGCCGCTCCCCAGGGGTGCGA  
 GCTGGACAGTGTCTCCTCCGCTGGCAACAGTGTCCCAATAGCCCATGGCCATGCTGCACATTG  
 GCTCCAACCTGAGAGGGAGTTGGATGATGTATTGACAACATTATGCGTCTGGACGATGCTCTGGCTA  
 CATCAATCTGAAATGCAGATGCCAACACGCTACCCCTGTCCAGCAGCCACCTGAATGTGTACAGCAGC  
 GACCCCAAGTACAGCCTCCCTGGTGGCGTCAACCAGCAGCTCCTGCCCTGCGGACCTGACCCAGAAGC  
 GAGAGCTCACAGATGCTGAGAGCAGGGCCCTGGCCAAGGAGCGGCAAGAAAGACAATCACAATTAAT  
 TGAAAGGAGACGAAGTTCAACATCAATGACCGCATCAAGGAGTTGGAAATGCTGATCCCAAGGCCAAT  
 GACCTGGACGTGCGCTGGAAACAAGGCCACCATCCTCAAGGCCCTGTGGATTACATCCGGAGGATGCAGA  
 AGGACCTGCAAAAGTCCAGGGAGCTGGAGAACCACTCTCGCCGCTGGAGATGACCAACAAGCAGCTCTG  
 GCTCCGTATCCAGGAGCTGGAGATGCAGGCTCGAGTGCACGGCTCCCTACCACCTCCCGTCCGGCATG  
 AACATGGCTGAGCTGGCCAGCAGGTGGTGAAGCAGGAGCTGCCTAGCGAAGAGGGCCAGGGGAGGCC  
 TGATGCTGGGGGCTGAGTCCCTGACCTGACCTGCCAGCTCTGCCCGCAAGCCCGCTGCCCT  
 GCCACCCAGCCACCATCCCATCCATCACCTGGACTTCAGCCACAGCCTGAGCTTTGGGGCAGGGAG  
 GACGAGGGTCCCGGGCTACCCGAACCCCTGGCGCGGGGATGGCTCCCATCCAGCCTGTCCA  
 AGAAGGATCTGGACCTCATGCTCCTGGACGACTACTGCTACCGTGGCCTCTGATCCACTTCTGTCCAC  
 CATGTCCCGGAGGCTCCAAGGCCAGCAGCCCGGAGCAGTTTTCAGCATGGAGGAGGGCGATGTGCTG

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

**Protein Sequence:**

>RG230141 representing NM\_001167827  
 Red=Cloning site Green=Tags(s)

MASRIGLRMQLMREQAQQEEQRERMQQQAVMHYMQQQQQQQQLGGPPTPAINTPVHFQSPPPVPGEVL  
 KVQSYLENPTSYHLQSQHQKVVREYLSEYGNKFAAHISPAQGSPPKPPAASPGVRAGHVLSSSAGNSAP  
 NSPMAMLHIGSNPERELDDVIDNIMRLDDVLYINPEMQMPNTLPLSSSHLNVYSSDPQVTASLVGVTSS  
 SCPADLTQKREL TDAESRALAKERQKDNHNL IERRRRFNINDRIKELGMLIPKANDLDVRWNGKILKA  
 SVDYIRRMQKDLQKSRELENHSRREMTNKQLWLRIQELEMQARVHGLPTTSPSGMNAELAQQVVKQEL  
 PSEEGPGEALMLGAEVPDPEPLPALPPQAPLPLPTQPPSPFHHLDFSHLSFGGREDEGPPGYPEPLAPG  
 HGSPFPLSKKDLDLMLLDDSLPLASDPLLSTMSPEASKASSRRSSFMEEGDVL

TRTRPLE – GFP Tag – V

**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**


**ACCN:** NM\_001167827

**ORF Size:** 1428 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001167827.3](#)

RefSeq Size: 2198 bp

RefSeq ORF: 1473 bp

Locus ID: 7942


UniProt ID: [P19484](#)

Cytogenetics: 6p21.1

Protein Families: Druggable Genome, Transcription Factors

**Gene Summary:** Transcription factor that specifically recognizes and binds E-box sequences (5'-CANNTG-3'). Efficient DNA-binding requires dimerization with itself or with another MiT/TFE family member such as TFE3 or MTF. In association with TFE3, activates the expression of CD40L in T-cells, thereby playing a role in T-cell-dependent antibody responses in activated CD4(+) T-cells and thymus-dependent humoral immunity. Specifically recognizes and binds the CLEAR-box sequence (5'-GTCACGTGAC-3') present in the regulatory region of many lysosomal genes, leading to activate their expression. It thereby plays a central role in expression of lysosomal genes. Acts as a positive regulator of autophagy by promoting expression of genes involved in autophagy. Specifically recognizes the gamma-E3 box, a subset of E-boxes, present in the heavy-chain immunoglobulin enhancer. Plays a role in the signal transduction processes required for normal vascularization of the placenta. Regulates lysosomal positioning in response to nutrient deprivation by promoting the expression of PIP4P1 (PubMed:29146937). [UniProtKB/Swiss-Prot Function]

### Product images:



Circular map for RG230141