

Product datasheet for **RG230208**

CUG BP1 (CELF1) (NM_001172639) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CUG BP1 (CELF1) (NM_001172639) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CUG BP1
Synonyms:	BRUNOL2; CUG-BP; CUGBP; CUGBP1; EDEN-BP; hNab50; NAB50; NAPOR
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG230208 representing NM_001172639
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCTGCGTTTTAAGTTGGATTTCTTCCAGAATGATGGTGGATCATTGCTCTTTGAATCCAGTCCCG
 TCTCAAAGAAAATGAACGCACCCTGGACCACCCAGACCAACCAGATCTTGATGCTATCAAGATGTTTGT
 GGGCCAGTTTCCAAGGACCTGGTCTGAAAAGGACTTGCAGGAACTCTCGAACAGTATGGTGTGTGTAT
 GAAATCAACGCTCCTAAGGGATAGGAGCCAAAACCCGCCTCAGAGCAAAGGGTGTGTTTTTTTACATTTT
 ACACCCGTAAAGCTGCATTAGAAGCTCAGAATGCTCTTCAACATGAAAGTCTCCAGGGATGCATCA
 CCCTATACAGATGAAACCTGCTGACAGTGAGAAGAACAATGTGGAAGACAGGAAGCTGTTTATTGGTATG
 ATTTCCAAGAAGTGCCTGAAAATGACATCCGAGTCATGTTCTCTCGTTTGGACAGATTGAAGAATGCC
 GGATATTGCGGGGACCTGATGGCCTGAGCCGAGGTTGTGCATTTGTGACTTTTACAACAAGGCCATGGC
 ACAGACGGCTATCAAGGCAATGCACCAAGCACAGACCATGGAGGGTTGCTCATCACCCATGGTGGTAAAA
 TTTGCTGATACACAGAAGGACAAAAGAAGAGAATGGCCAGCAGCTCCAGCAGCAGATGCAGCAAA
 TCAGCGCAGCATCTGTGTGGGAAACCTTGTGCTAAATACTCTTGGACCCAGTATTTAGCACTTTA
 TTTGCAGCTCCTTTCAGCAGACTGCCTCCTCTGGGAACCTCAACACCCTGAGCAGCCTCCACCAATGGGA
 GGGTTGAATGCAATGCAGTTACAGAATTTGGCTGCACTAGCTGCTGCAGCTAGTGCAGCTCAGAACACAC
 CAAGTGGTACCAATGCTCTCACTACATCCAGCAGTCCCTCAGCGTGTCTACTAGTTCAGGGTCTCACC
 TAGCTCTAGCAGCAGTAATCTGTCAACCCATAGCCTCACTTGGAGCCCTGCAGACATTAGCTGGAGCA
 ACGGCTGGCCTCAATGTTGGCTCTTTGGCAGGAATGGCTGCTTTAAATGGTGGCCTGGGCAGCAGTGGCC
 TTTCCAATGGCACCCGGGAGCACCATGGAGGCCCTCACTCAGGCCTACTCGGGTATCCAGCAATATGCTGC
 TGCTGCGCTCCCACTCTGTACAACCAGAATCTTCTGACACAGCAGAGTATTGGTGTGCTGGAAGCCAG
 AAGGAAGTCCAGAGGGAGCCAACCTGTTTATCTACCACCTGCCCCAGGAGTTTGGTGTGATCAGGACCTGC
 TGCAGATGTTTATGCCCTTTGGGAATGTCGTGTCTGCCAAGGTTTTTCATAGACAAGCAGACAAAACCTGAG
 CAAGTGTGTTTGGTGTGTAAGTTACGACAATCCTGTTTCGGCCCAAGCTGCCATCCAGTCCATGAACGGC
 TTTAGATTGGCATGAAGCGGCTTAAAGTGCAGCTCAACGTTTGAAGAATGACAGCAAGCCCTAC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG230208 representing NM_001172639
 Red=Cloning site Green=Tags(s)

MAAFKLDLPEMVDHCSLNSSPVSKMNGTLDHPDQPDLDIAIKMFVGVQVPRTWSEKDLRELFEQYGA
 VEINVLDRSRQNPQSKGCCFVTFYTRKAALAEQNALHNMKVLPGMHHPIQMKPADSEKNNVEDRKLFIGM
 ISKKTENDIRVMFSSFGQIEECRILRGPDLRGCAFVTFTRAMAQTAIKAMHQAQTMEGCSSPMVVK
 FADTQKDKQKRMAQQLQQMQQISAAASVWGNLAGLNTLGPQYLALYLQLLQQTASSGNLNTLSSLHPMG
 GLNAMQLQNLAALAAAASAAQNTPSGTNALTSSSPLSVLTSSGSSPSSSSNSVNPISLALQTLAGA
 TAGLNVGSLAGMAALNGGLGSSGLSNGTGSTMEALTAQYSGIQYAAAALPTLYNQNLTTQQSIGAAGSQ
 KEGPEGANLFIYHLPQEFGDQDLLQMFMFVGNVSAKVFIDKQTNLSKCFGVSYDNPVSAQAAIQSMNG
 FQIGMKRLKVQLKRSKNDKPY

TRTRPLE – GFP Tag – V

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:


ACCN: NM_001172639

ORF Size: 1536 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001172639.2](#)

RefSeq Size: 7434 bp

RefSeq ORF: 1539 bp

Locus ID: 10658


UniProt ID: [Q92879](#)

Cytogenetics: 11p11.2

Protein Families: Druggable Genome

Gene Summary: Members of the CELF/BRUNOL protein family contain two N-terminal RNA recognition motif (RRM) domains, one C-terminal RRM domain, and a divergent segment of 160-230 aa between the second and third RRM domains. Members of this protein family regulate pre-mRNA alternative splicing and may also be involved in mRNA editing, and translation. This gene may play a role in myotonic dystrophy type 1 (DM1) via interactions with the dystrophin myotonic-protein kinase (DMPK) gene. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG230208