

Product datasheet for **RG230413**

TTC39B (NM_001168340) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TTC39B (NM_001168340) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	TTC39B
Synonyms:	C9orf52
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG230413 representing NM_001168340
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGATGTCTGTGCTGGCCTGCCGCTCCGGGGCGGGGAAATCGAGTAGCTGCACTCAGGCCCGCCCGA
GACCCGGAGGGAGCGCCGGCCGTCCTTTTGCCTCCTTTGCCTGGGCTGAGCCAGAGCCGAGAGC
AGGAGTCGGCTCTGAGTCCCTGCTTGGTTTTGGGTGGCAGCAGCCAGAGGAGGAATATGGCGCTTTA
GGCAGCCGAGCGGAGCTGGAAGCGGACGAGGACGTTTTCGAAGATGCCTTGAAACCATCTCAATATCTT
CTCACTCAGATATGGCAACAAGCAGCCTTCACTTTGCATCATGTGATACACAACAGGCGCCAGACAGCG
CGGAGCGTCTACTGTCAGCAGTTCATCATCAACCAAGGTGGATCTCAAGAGTGGCCTCGAAGAATGTGCT
GTGGCATTGAACTATTTCTAAGCAACAAATTTACAGACGCTTAGAATTGCTTCGCCCCGGGCTAAGG
AGAGTATGTACCATGCCTTGGGCTACAGTACCATTGTGGTGTGCAGGCTGTCTGACCTTCGAGCAACA
GGACATCCAGAACGGCATTCTGCCATGAAGGACGCTTTACAAACCTGCCAAAAATACAGGAAAAATAC
ACAGTTGTAGAATCTTTCTCAAGTCTTCTTTCTAGAGGATCCCTGGAGCAACTGAGTGAAGAGGAAATGC
ATGCTGAAATCTGTTATGCCGAGTGTCTCTACAGAAAGCTGCACTCACGTTTGTGCAGGATGAAATAT
GATCAACTTCATCAAAGGTGGACTCAAAATTAGAACAAGTTACCAAAATATATAAAGAATGTCTTTCTATC
CTGCATGAAATTCAGAAGAACAACTTCAGCAGGAATCTTCTATGAGTTTGAAGGGGTGTCAAGCTTG
GCAGTGGGGCCTTAATTTGATGCTGTCCCTTTACCAGCAAGGATTATCAGACTACTAGAATTTATTGG
ATTTCTGGGAATAGGGAGTTGGCCCTCTGCAGTTACGTGAAGGTGCTTCAGGAAGGAGCATGAGGTCT
GCACTGTGCTGCTTAACCATCTTAGCTTTTCACTTACATCTCCCTAATACTTGGTACAGGAGAAGTGA
ATGTTGCGGAAGCAGAGCGTCTCCTGGCACCTTCCCTCCAGCAGTTTCCAAATGTACGCTCAGTATTTCA
AAGTGCACAAGAAGTATTTCAAAAATGCATTTTCAGTTCAAGAAGAATGGAAACAGTTTCACCATCTCTGC
TACTGGGAGCTAATGTGGATTAATGTTTTCCAACAAAAGTGGATGCAGGCATATTACTATTCAGATCTGC
TTTGCAAAGAGAGTAAATGGTCCAAGGCACTTATGTGTTCTTGAAAGCAGCAATTTTGAGTATGCTTCC
AGAGGAGGATGTAGTAGCAACGAATGAGAATGTGGTAACTTTATTCAGACAGGTGGACAGCTTGAAGCAG
AGAATTGCCGGGAAATCTATCCCTACTGAGAAGTTTGTGTGAGGAAGGCTCGGCCTTACTCCGCCTCCT
TACCTGCACCTGTGAAGCTCATTTACCTGCCCTGAAATGATGTATGTCTGGAATGGTTTTTCAATAGT
GAGCAAAAAGAAAAGACCTTTCTGAAAATCTGTTAGTAACTGTTGAAAAAGCTGAGGCAGCTTTACAAAGT
CAAAATTTTAAACAGCTTCTCTGTGGATGATGAGTGCTTAGTGAAGTTACTTAAAGGATGTTGCCTCAAGA
ACTTACAGCGGCCCTTGCAAGCTGAACTATGTTACAATCATGTTGTGAAAGTGAAGGCTACTGAAGTA
TGACCACTACCTAGTGCCGTTCACTCTATTTGAATTGGCATCTTTGTATAAAAGCCAGGGGGAAATTGAC
AAGGCCATAAAGTTCTAGAAACTGCAAGGAACAACACTACAAAGATTACTCCCTGGAGTCCAGACTACACT
TCAGAATTCAGGCAGCTCTTACCTCTGGAGAAAACCTTCTCAGAT

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG230413 representing NM_001168340
Red=Cloning site Green=Tags(s)

```
MDAVLACRLRGRGNRVAALRPRPRPGGSAGPSPFALLCAGLSPEPRAGVGSEFPWFLLGGSSQRRNMALL
GSRAELEADEDVFEDALETISISSHSDMATSSSLHFASCDTQQAPRQRGASTVSSSSSTKVDLKSGLLECA
VALNLFSLNKFDTALELLRPWAKESMYHALGYSTIVVLQAVLTFEQQDIQNGISAMKDALQTCQKYRKKY
TVVESFSSLLSRGSLEQLSEEEMHAEICYAECLLQKAALTFVQDENMINFIKGGKIRTSYQIYKECLSI
LHEIQKNKLQEFFYEFEGGVKLGSGAFNLMLSLLPARIIRLLEFIFGSGNRELGLLQLREGASGRSMRS
ALCCLTILAFHTYISLILGTGEVNVAEERLLAPFLQQFPNVRVVFQSAQEVFQKCSVQEWEKQFHHLIC
YWELMWINVFQQNWMQAYYSDLLCKESKWSKATYVFLKAAILSMLPEEDVVATNENVVTLFRQVDSLKQ
RIAGKSIPTEKFAVRKARRYSASLPAPVKLILPALEMMYVWNGFSIVSKRKDLSENLLVTEKAEAAALQS
QNFNSFSVDDECLVKLLKGCCLKNLQRPLQAEALCYNHVVESEKLLKYDHYLVPFTLFELASLYKSQGEID
KAIKFLETARNNYKDYSLESRLHFRIQAALHLWRKPSSD
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001168340

ORF Size: 2007 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001168340.1](#), [NP_001161812.1](#)

RefSeq Size: 3237 bp

RefSeq ORF: 2010 bp


Locus ID: 158219

UniProt ID: [Q5VTQ0](#)

Cytogenetics: 9p22.3

Gene Summary: Regulates high density lipoprotein (HDL) cholesterol metabolism by promoting the ubiquitination and degradation of the oxysterols receptors LXR (NR1H2 and NR1H3). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG230413