

Product datasheet for **RG230423**

SUN1 (NM_001171944) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SUN1 (NM_001171944) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SUN1
Synonyms:	UNC84A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG230423 representing NM_001171944
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGATTTTTCTCGGCTTACATGTACAGTCTCCCCAGTGTGTGCCGGAGAACACGGGCTACACGTATG
CGCTCAGTTCAGCTATTCTTCAGATGCTCTGGATTTTGAGACGGAGCACAAATTGGACCTGTATTGA
TTCTCCACGGATGTCCCCTAGTTTGCCTGGCCACGACAGCATGCACCCTGGGGATGGTGAGGCT
GTGGGTGCCGACAGCGGCACCAGCAGCGTGTCTCCCTGAAGAACCAGCGGCCAGAACAACAAACAGC
GCAGAAGCACAAACAAATCAGCTTTTAGTATCAACCACGTGTCAAGGCAGGTACAGTCTCTGGCGTCA
CCACGGCGGCACTGTGAGCTGCAGGATGCTGTGACTCGACGGCCTCTGTATTGGACGAGTCTTGGATT
CGTGAACAGACCACAGTGGACACTTCTGGGGAAAGGCAGCCTCTGGAGTGTCTGGTGGCTGGGGATTG
GATGGTACCAGTTTGTACTTTGATTTCTGGTGAATGTGTTTCTTACCAGGTGCCTCGAAACAT
CTGCAAGTTTTAGTCTTCTCATCCACTCTTCTTTTACTAGGTCTCTCTTACGGGGCCAGGGCAAT
TTCTTTTCGTTCTTGCCGTGTTGAACTGGGCAAGCATGCATAGAACACAGCGGGTGGATGACCCCCAGG
ACGTGTTTAAACCCACGACTTCTCGCTGAAGCAGCCTCTGCAGGGTGACAGTGAGGCTTTTCCGTGGCA
TTGGATGAGTGGCGTGGAGCAGCAGGTGGCCTCTCTGTCTGGACAGTGCCACCACCATGGTGAAGTCTC
CGAGAGCTGACCACTTTGCTACAGAAGCTGCAGGCTCGGGTGGACCAGATGGAAGCGCGCTGCCGGGC
CGTCAGCTTCGGTCAGAGACGCTGTGGGACAGCCCCGAGGGAGACTGACTTTATGGCCTTTCACCAAGA
ACATGAAGTGGTATGTCACACTTGAAGATATTCTGGGAAAACAGAGAAAAATCTGAGGCCATCCAG
AAGGAAGTAGAACAGACCAAGCAAAAAACAATCAGTGGGTTGGTGAGCAGCTCCTGCCACAGTCGAGC
ACCTCCAGCTGGAGCTGGATCAGCTAAAGTCAGAGCTGTCCAGCTGGCGACACGTGAAGACCGGCTGTGA
GACAGTGGATGCCGTACAAGAAAGAGTGGACGTGCAAGTCAGAGAAATGGTGAAACTCCTGTTTTCCGAA
GATCAGCAAGGCGTTCTCTGGAACAGCTGCTGCAGAGTTTCTCATCAGTTTTGTGAGCAAAGGCGACT
TGCAGACGATGCTGCGAGACCTGCAGCTGCAGATCCTGCGGAACGTACCCACCACGTTTCCGTGACCAA
GCAGCTCCCAACCTCAGAAGCCGTGGTGTCTGCTGTGAGCGAGGCGGGGGCGTCTGGAATAACAGAGGCG
CAAGCACGTGCCATCGTGAACAGCGCCTTGAAGCTGTATTCCCAAGATAAGACCGGGATGGTGGACTTTG
CTCTGGAATCTGGTGGTGGCAGCATCTTGAAGTACTCGTGTCTGAAACTTACGAAACAAAACGGCGCT
GATGAGTCTGTTGGGATCCCGCTGTGGTACTTCTCGCAGTCCCCGCGCTGGTATCCAGCCTGACATT
TACCCCGTAACTGCTGGGCATTTAAAGGCTCCAGGGTACCTGGTGGTGGGCTCTCCATGATGATCC
ACCCAGCCGCTTCACTCTGGAGCACATCCCTAAGACGCTGTGCGCAACAGGCAACATCAGCAGCGCCCC
CAAGGACTTCGCCGTCTATGGATTAGAAAATGAGTATCAGGAAGAAGGGCAGCTTCTGGGACAGTTCAGG
TATGATCAGGATGGGGAGTCGCTCCAGATGTTCCAGGCCCTGAAAAGACCCGACGACACAGCTTTCCAAA
TAGTGGAACTTCGGATTTTTCTAACTGGGGCCATCCTGAGTACCTGTCTGTATCGGTTCCAGAGTTCA
TGGCGAACCTGTCAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG230423 representing NM_001171944
 Red=Cloning site Green=Tags(s)


```
MDFSRLHMYSPPQCVPE NTGYTYALSSSYSSDALDFE TEHKLDPVFDS PRMSRRSLRLATTACTLGDGEA
VGADSGTSSAVSLKNRAARTTKQRRSTNKSAFSINHVS RQVTSSGVSHGGTVSLQDAVTRRPVLD ESWI
REQTVDHFHWGKAASGVFWLGI GWYQFVTLISWLVN FLLTRCLRNICKFLVLLIPLFLLGLSLRGQGN
FFSFLPVLN WASHMRTQRVDDPQDVFKPTTSRLKQPLQGDSEAF PWHWMSGVEQQVASLSGQCHHHGENL
RELTTLLQKLQARVDQMEGGAAGPSASVRDAVGQPPRE TDFMAFHQEHEVRMSHLEDILGKLRKSEAIQ
KELEQTKQKTISAVGEQLLPTVEHLQLELDQLKSELSSWRHVKTGCETVDAVQERVDVQVRE MVKLLFSE
DQQGGSLEQLLQRFSSQFVSKGDLQTMLRDLQLQILRN VTHHVSVKQLPTSEAVVS AVSEAGASGITEA
QARAIVNSALKL YSQDKTGMVDFALESGGGSILSTRCSE TYETKTALMSLFGIPLWYFSQSPRVVIQ PDI
YPGNCWAFKGSQ GYLVVRLSMMIHPAAFTLEHIPKTL SPTGNISSAPKDFAVYGLENEYQE EGQLLGQFT
YDQDGESLQMFQAL KRPDDTAFQIVELRIFSNWGHPEY TCLYRFRVHGEPVK
```

TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001171944

ORF Size: 2046 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001171944.2](#)

RefSeq Size: 3727 bp

RefSeq ORF: 2049 bp

Locus ID: 23353


UniProt ID: [O94901](#)

Cytogenetics: 7p22.3

Protein Families: Transmembrane

Gene Summary: This gene is a member of the unc-84 homolog family and encodes a nuclear envelope protein with an Unc84 (SUN) domain. The protein is involved in nuclear anchorage and migration. Alternatively spliced transcript variants have been described. [provided by RefSeq, Jan 2019]

Product images:



Circular map for RG230423