

## Product datasheet for **RG230432**

### **NKRF (NM\_001173488) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	NKRF (NM_001173488) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	NKRF
Synonyms:	ITBA4; NRF
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>RG230432 representing NM\_001173488  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGAAAAATTCTCAAATGGCTGAAGGTATTGATATTGGGAGATGCCTTCATATGATCTGGTGTCTGT  
CCAAACCTTCCAAAGGTCAAAAACGCCACTCTCAACATGTGATGGTCAAAATCCTCCTAAAAAGCAAGC  
CGTTCCAAATCCATGCGAGACCTCGTTTTGAGCCTGTACATTTGTAGCTAGTAGTTCAAAGATGAA  
AGACAGGAAGATCCTTATGGCCCTCAAACAAAAGAGGTAATGAACAAACACATTTTCCAGCATGCCAA  
GAGACATCTACCAAGATTATACTCAAGACTCTTTCAGTATACAAGATGGGAATTCTCAGTATTGTGATTC  
ATCAGGATTCATTCTCACAAAAGACCAGCCTGTAACAGCCAACATGTATTTTGACAGTGGGAACCTGCC  
CCAAGCACCACATCACAGCAGGCAAACCTCAGTCAACTCCTGAGCCTTACCATCACAGACATTTCCCG  
AGTCTGTGGTAGCCGAGAAGCAGTATTTTATTGAAAAATTAACGGCGACAATCTGGAAGAACCTTTCTAA  
TCCAGAAATGACTTCTGGATCTGATAAAATTAATTATACATATATGTTAACTCGTTGTATTACAGCGTGT  
AAGACAAATCCTGAGTATATATGCTCCTTTAAAGGAAATTCCTCCTGCCGACATCCCCAAAAATAAAA  
AACTTCTAACTGATGGCTATGCTTGTGAAGTTAGATGCCAAAATATCTACTTAACTACAGGTTATGCTGG  
CAGCAAGAATGGGTCCAGGGATCGAGCTACAGAGCTAGCTGTAAGTCTTGCAGAAACGTATTGAAGTT  
AGAGTTGTCCGGCGGAAATTCAGCATAACATTTGGAGAGGACCTCGTGGTGTGTGAGATTGGCATGTCTCT  
CCTATGAATTTCTCCAGCTCTGAAGCCACCAGAAGACCTGGTGGTGTGGGTAAAGATGCTTCCGGGCA  
GCCAATTTTTAATGCTTCTGCCAAACACTGGACCAATTTTGTATTACAGAAAATGCAAAATGATGCAATT  
GGTATCCTTAACAATTCTGCCTCATTCAACAAGATGTCAATTGAATACAAATATGAGATGATGCCAAATC  
GCACATGGCGTTGTGAGTGTTTTTACAAGATCACTGCTTAGCTGAAGGTTATGGAACCAAGAAAACAAG  
TAAACATGCAGCTGCCGACGAGGCTTTGAAAATCTTCAAAAAACACAGCCCACTTATCCATCTGTCAAA  
AGTTCACAATGCCATACAGGCTTTCACCCAGAGGATCTGGAAAGAAGAAAGATATAAAGGATCTTGTAG  
TTTATGAGAATCTTCAAATCCCGTGTGCACGCTGAACGACACAGCTCAGTTTAAACCGAATGACAGTTGA  
GTATGTCTATGAAAGGATGACAGGCCTCCGCTGGAAATGCAAAGTGATTCTAGAGAGTGAAGTAATTGCA  
GAAGCAGTTGGGTGAAGAAAACGTCAAATATGAAGTGTGGGGAAGCTGTGAAAACCTCAAAAAGA  
CCCAGCCAACCTGTCATTAACAACCTGAAGAAAGGAGCTGTTGAAGATGTGATTTCAAGAAATGAAATTCA  
GGGCCGCTCAGCAGAGGAGGCTTACAAACAGCAAATCAAAGAAGATAATATTGAAATCAGCTGCTGAGA  
AAGATGGGTTGGACTGGTGGTGGTTTAGGTAATCTGGTGAGGGCATAACGGGAGCCTATCTCAGTAAAAG  
AGCAGCATAAGCGGGAAGGGCTTGGTCTGGATGTAGAGAGGGTGAATAAAATGCCAAGAGAGATATTGA  
ACAGATCATCAGAAACTACGCCCGCTCCGAGAGCCACACAGATTTGACTTCTCTAGAGAGCTGACTAAT  
GATGAACGGAAGCAAATACATCAGATTGCCAGAAGTATGGTCTTAAGAGTAAGTCTCATGGGGTGGGCC  
ATGATAGGTACCTAGTGGTAGGTAGAAAAGACGGAAGGAAGACCTACTAGATCAGCTCAAACAGGAAGG  
CCAAGTGGGCCATTACGAGCTTGTATGCCTCAAGCAAAT

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >RG230432 representing NM\_001173488  
 Red=Cloning site Green=Tags(s)


MEKILQMAEGIDIGEMPSYDLVLSKPSKGQKRHLSTCDGQNPCKQAGSKFHARPRFEPVHFVASSSKDE  
 RQEDPYGPQTKEVNEQTHFASMPRDIYQDYTDQSF SIQDGN SQYCDSSGFILTKDQPV TANMYFDSGNPA  
 PSTTSQQANSQSTPEPSPSQTFPE SVVAEKQYF IEKLTATIWKNL SNPEMTSGSDKIN YTYMLTRCIQAC  
 KTNPEYIYAPLKEIPPADIPKKNKLLTDGYACEVRCQNIYLTTGYAGSKNGSRDRATELAVKLLQKRIEV  
 RVVRRKFKHTFGEDLVVCQIGMSSYEFPPALKPPEDLVVLGKDASGQPIFNASAKHWTNFVITENANDAI  
 GILNNSASFNKMSIEYKYEMMPNRTWRCRVFLQDHCLAEGYGTKKT SKHAAADEALKILQKTQPTYPSVK  
 SSQCHTGS SPRSGKKKDIKDLVYENSSNPVCTLNDTAQFNRM TVEYVYERMTGLRWKCKVILESEVIA  
 EAVGVKKTVKYEAAGEAVKTLKKTQPTVINLKKGAVEDVISRNEIQGRSAEEAYKQKIKEDNIGNQLLR  
 KMGWTGGGLGKSGEGIREPISVKEQHKREGLDVERVNKIAKRDI EQIRNYARSESHTDLTFSREL TN  
 DERKQIHQIAQKYGLKSKSHGVGHDRYL VVGRKRRKEDLLDQLKQEGQVGHYELVMPQAN


TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



**Plasmid Map:**


**ACCN:** NM\_001173488

**ORF Size:** 2070 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001173488.1](#), [NP\\_001166959.1](#)

**RefSeq Size:** 3747 bp

RefSeq ORF: 2073 bp

Locus ID: 55922

UniProt ID: [O15226](#)

Cytogenetics: Xq24

Protein Families: Transcription Factors

**Gene Summary:** This gene encodes a transcriptional repressor that interacts with specific negative regulatory elements to mediate transcriptional repression of certain nuclear factor kappa B responsive genes. The protein localizes predominantly to the nucleolus with a small fraction found in the nucleoplasm and cytoplasm. Alternate splicing results in multiple transcript variants.  
[provided by RefSeq, Mar 2010]