

## Product datasheet for **RG230435**

### Angiotensin Converting Enzyme 1 (ACE) (NM\_001178057) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Angiotensin Converting Enzyme 1 (ACE) (NM_001178057) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Angiotensin Converting Enzyme 1
Synonyms:	ACE1; CD143; DCP; DCP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RG230435 representing NM\_001178057  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGGCCAGGGTTGGGCTACTGCAGGACTTCCCAGCCTCCTCTTCTGCTGCTCTGCTACGGGCACCCTC  
 TGCTGGTCCCCAGCCAGGAGGCATCCCAACAGGTGACAGTCACCCATGGGACAAGCAGCCAGGCAACAAC  
 CAGCAGCCAGACAACCACCCACCAGGCGACGGCCACCAGACATCAGCCAGAGCCAAAACCTGGTGACT  
 GATGAGGCTGAGGCCAGCAAGTTTGTGGAGGAATATGACCGGACATCCCAGGTGGTGTGGAACGAGTATG  
 CCGAGGCCAACTGGAACAACACCAACATCACCACAGAGACCAGCAAGATTCTGCTGCAGAAGAACAT  
 GCAAATAGCCAACCACACCCTGAAGTACGGCACCCAGGCCAGGAAGTTTGTATGTGAACAGTTGCAGAAC  
 ACCACTATCAAGCGGATCATAAAGAAGGTTTCAGGACCTAGAACGGGCAGCACTGCCTGCCAGGAGCTGG  
 AGGAGTACAACAAGATCCTGTTGGATATGGAACACCTACAGCGTGGCCACTGTGTGCCACCCGAATGG  
 CAGCTGCCTGCAGCTCGAGCCAGATCTGACGAATGTGATGGCCAGTCCCGGAAATATGAAGACCTGTTA  
 TGGGCATGGGAGGGCTGGCGAGACAAGGCGGGGAGAGCCATCCTCCAGTTTTACCCGAAATACGTGGAAC  
 TCATCAACCAGGCTGCCCGGCTCAATGGCTATGTAGATGCAGGGGACTCGTGGAGGTCTATGTACGAGAC  
 ACCATCCCTGGAGCAAGACCTGGAGCGGCTCTCCAGGAGCTGCAGCCACTCTACCTCAACCTGCATGCC  
 TACGTGCGCCGGGCCCTGCACCGTCACTACGGGGCCAGCACATCAACCTGGAGGGGCCATTCTGCTC  
 ACCTGTGGGGAACATGTGGGCGCAGACCTGGTCCAACATCTATGACTTGGTGGTGGCCCTCCCTTCAGC  
 CCCCTCGATGGACACCACAGAGGCTATGCTAAAGCAGGGCTGGACGCCAGGAGGATGTTAAGGAGGCT  
 GATGATTTCTTACCTCCCTGGGGCTGCTGCCCGTGCCTCCTGAGTTCTGGAACAAGTCGATGCTGGAGA  
 AGCCAACCGACGGGCGGAGGTGGTCTGCCACGCCTCGGCCTGGGACTTCTACAACGGCAAGGACTCCG  
 GATCAAGCAGTGCACCACCGTGAACCTGGAGGACCTGGTGGTGGCCACCACGAAATGGGCCACATCCAG  
 TATTTTCATGCAGTACAAAGACTTACCTGTGGCCTTGAGGGAGGGTGCCAACCCCGGCTTCCATGAGGCCA  
 TTGGGGACGTGCTAGCCCTCTCAGTGTCTACGCCAAGCACCTGCACAGTCTCAACCTGCTGAGCAGTGA  
 GGGTGGCAGCGACGAGCATGACATCAACTTTCTGATGAAGATGGCCCTTGACAAGATCGCCTTTATCCCC  
 TTCAGCTACCTCGTCGATCAGTGGCGCTGGAGGGTATTTGATGGAAGCATCACCAGGAGAACTATAACC  
 AGGAGTGGTGGAGCCTCAGGCTGAAGTACCAGGGCCTTGCCCCCAGTGCCAGGACTCAAGGTGACTT  
 TGACCCAGGGGCAAGTCCACATTCCTTCTAGCGTGCCTTACATCAGGACCGCCATGAAGCTGGGCTTC  
 AGTAGGCCGTGGCCGGAAGCCATGCAGCTGATCAGGGCCAGCCAAACATGAGCGCCTCGGCCATGTTGA  
 GCTACTTCAAGCCGCTGCTGGACTGGCTCCGCACGGAGAACGAGCTGCATGGGGAGAAGCTGGGCTGGCC  
 GCAGTACAACGGACGCCGAACCTCCGCTCGCTCAGAAGGGCCCTCCCAGACAGCGGCCGCGTCAGCTTC  
 CTGGGCCCTGGACCTGGATGCGCAGCAGGCCCGCTGGGCCAGTGGCTGCTGCTTCTCTGGGCATCGCCC  
 TGCTGGTAGCCACCCTGGGCCTCAGCCAGCGGCTCTTACGATCCGCCACCAGCCTCCACCGGCCTC  
 CCACGGGCCCCAGTTCGGCTCCGAGGTGGAGCTGAGACTCC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >RG230435 representing NM\_001178057  
 Red=Cloning site Green=Tags(s)


MGQGWATAGLPSLLFLLLCYGHPLLVPSEQEASQVTVTHGTSSQATTSSQTTTHQATAHQ TSAQSPNLVT  
 DEAEASKFVEEYDRTSQVWVNEYAEANWNYNTNITTETSKILLQKNMQIANHTLKYG TQARKFDVNQLQN  
 TTIKRIIKKVQDLERAALPAQELEEYNKILLDMETTYSVATVCHPNGSCLQLEPDL TNVMATSRKYEDLL  
 WAWEGWRDKAGRAILQFYPKYVELINQAARLNGYVDAGDSWRSMYETPSLEQDLERLFQELQPLYLNLHA  
 YVRRALHRHYGAQHINLEGIPIAHL LGNMWAQTWSNIYDLVVPFSPASMDTTEAMLKQGWTPRRMFKEA  
 DDFFTSLGLLPVPPF FWNKSMLEKPTD GREVVCHASAWDFYNGKDFRIKQCTTVNLEDLVVAHHEMGIQ  
 YFMQYKDL PVALREGANPGFHEAIGDVLALSVSTPKHLHSLNLLSSEGG SDEHDINFLMKMALDKIAFIP  
 FSYLVQWRWRVFDGSITKENYQEWWSLR LKYQGLCPPVPRTQGDFDPGAKFHIPSSVPIYRTAMKLG  
 SRPWPEAMQLITGQPNMSASAMLSYFKPLLDWLR TENELHGEKLGWPQYNWTPNSARSEGPLPDSGRV  
 SFGVSRVSLGLDLDAQARVQWLLFLGLIALLVATLGLSQRLF SIRHRS LHRHSHGPFQFGSEVELRHS

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



**ACCN:** NM\_001178057

**ORF Size:** 2073 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001178057.2](#)

**RefSeq Size:** 3141 bp

**RefSeq ORF:** 2076 bp

**Locus ID:** 1636

**UniProt ID:** [P12821](#)


**Cytogenetics:** 17q23.3

**Protein Families:** Druggable Genome, ES Cell Differentiation/IPS, Protease, Secreted Protein, Transmembrane

**Protein Pathways:** Hypertrophic cardiomyopathy (HCM), Renin-angiotensin system

**Gene Summary:** This gene encodes an enzyme involved in blood pressure regulation and electrolyte balance. It catalyzes the conversion of angiotensin I into a physiologically active peptide angiotensin II. Angiotensin II is a potent vasopressor and aldosterone-stimulating peptide that controls blood pressure and fluid-electrolyte balance. This angiotensin converting enzyme (ACE) also inactivates the vasodilator protein, bradykinin. Accordingly, the encoded enzyme increases blood pressure and is a drug target of ACE inhibitors, which are often prescribed to reduce blood pressure. This enzyme additionally plays a role in fertility through its ability to cleave and release GPI-anchored membrane proteins in spermatozoa. Many studies have associated the presence or absence of a 287 bp Alu repeat element in this gene with the levels of circulating enzyme. This polymorphism, as well as mutations in this gene, have been implicated in a wide variety of diseases including cardiovascular pathophysiologies, psoriasis, renal disease, stroke, and Alzheimer's disease. Regulation of the homologous ACE2 gene may be involved in progression of disease caused by several human coronaviruses, including SARS-CoV and SARS-CoV-2. Alternative splicing results in multiple transcript variants encoding both somatic (sACE) and male-specific testicular (tACE) isoforms. [provided by RefSeq, Sep 2020]

## Product images:



Circular map for RG230435