

Product datasheet for **RG230445**

DDX4 (NM_001166533) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DDX4 (NM_001166533) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	DDX4
Synonyms:	VASA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG230445 representing NM_001166533
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGGAGATGAAGATTGGGAAGCAGAAATCAACCCTCATATGTCTTCTATGTTCCCATATTTGAGAAGG
 ATAGGTATTCTGGAGAAAATGGAGACAATTTAACAGGACTCCAGCTTCATCATCAGAAATGGATGATGG
 ACCTTCTCGAAGAGATCATTTCATGAAAAGTGGATTTGCCTCTGGGCGGAATTTGGAAAACAGAGATGCT
 GGTGAGTGAATAAGCGAGATAATACATCCACAATGGGTGGTTTTGGAGTTGAAAAGAGTTTTGGAAACA
 GAGGTTTTTCAAACAGCAGGTTTGAAGATGGTGATAGCTCTGGTTTTCTGGAGAGGCTATCGAGATGGAAA
 TAATTCAGAAGCTTCAGGGCCATACAGAAGAGGTGGAAGAGGTAGTTCCGAGGTTGCCGTGGAGGATTT
 GGTCTAGGAAGTCCAATAATGACTTAGACCCAGACGAATGTATGCAGCGCACTGGTGGCCTTTTTGGTT
 CTAGAAGACCAGTATTAAGTGGCACAGGTAATGGTGATACTTCTCAAAGCAGAAGTGGCAGTGGAAAGTGA
 ACGAGGTGGTTACAAAGTTTTAAATGAAGAAGTAATAACAGGCTCTGAAAAGAATTCTTGGAAAGTCAGAA
 GCAGAAGGAGGAGAAAAGTAGTGATACTCAAGGACCAAAAGTGACCTACATACCCCTCCTCCACCTGAGG
 ATGAGGACTCCATCTTTGCACATTATCAGACAGGCATAAACTTCGACAAAACGACACTATTCTTGTGGA
 AGTGTCTGGACATGATGCACCACCAGCAATTCTGACTTTTGAAGAAGCTAATCTCTGTGACAGACTGAAT
 AACAACATTGCTAAAGCTGGTTATACTAAGCTTACTCCTGTGCAAAAATACAGTATTCTATCATACTTG
 CAGGACGAGATTTGATGGCTTGCCTCAACAGGGTCTGGGAAGACTGCGGCTTTTCTCCTACCAATTTT
 GGCTCATATGATGCATGATGGAATAACTGCCAGTCGTTTTAAAGAGTTGCAGGAACCAGAGTGATATTATT
 GTAGCACCAACTCGAGAATTGGTCAACCAGATTTATTTGGAAGCCAGAAAATTTCTTTTGGGACTTGTG
 TAAGAGCTGTTGTTATATATGGGGGAACCCAGCTGGGACATTCAATTCGACAAAATAGTACAAGGCTGTAA
 TATATTATGTGCTACTCCTGGAAGACTGATGGATATCATAGGCAAGAAAAGATTGGTCTCAAACAGATC
 AAATACTTAGTTTTGGATGAAGCTGATCGCATGTTGGATATGGTTTTGGTCCAGAAATGAAGAAGTTAA
 TTTCTTGGCCAGGAATGCCATCAAAGGAACAGCGCCAAACCTTATGTTCAAGTCAACTTTTCCAGAGGA
 AATTCAAAGGTTGGCTGCAGAGTTTTTAAAGTCAAATTATCTGTTTGTGCTGTTGGACAAGTGGGTGGA
 GCATGTAGAGATGTTGAGCAGACCGTTCTCAAGTTGGCCAGTTCTCAAAAAGAGAAAAGCTCGTTGAAA
 TTCTGCGAAACATAGGGGATGAAAGAATATGGTCTTTGTTGAAACTAAGAAAAAGCAGATTTTATTGC
 AACTTTTCTTTGTCAAGAAAAAATCAACTACAAGTATTCATGGTGATCGGGAACAGAGAGCGGGAG
 CAAGCTCTTGGAGATTTTCGCTTTGAAAAGTGCCAGTCTTGTGCTACTTCAGTAGCTGCCAGAGGGC
 TGGATATTGAAAATGTGCAACATGTTATCAATTTTATGATCTTCTTCTACCATTGATGAATATGTTTATCG
 AATTGGGCGTACTGGTCGTTGTGGGAATCTGGCAGAGCAATTTCTTTTTTGTGATCTTGAATCGGATAAC
 CATTTAGCACAGCCTCTAGTAAAAGTATTGACAGATGCTCAACAGGATGTTCTGCATGGTTGGAAGAAA
 TTGCCCTTAGTACATACATTCTGGCTTCAAGTGGTAGTACAAGAGGAAACGTGTTTGCATCAGTTGATAC
 CAGAAAGGGCAAGAGCACTTTGAACACAGCTGGGTTTTCTTCTTCAAGCTCCCAATCCAGTAGATGAT
 GAGTCATGGGAT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG230445 representing NM_001166533
 Red=Cloning site Green=Tags(s)

MGDEDWEAEINPHMSSYVPIFEKDRYSGENDNFNRTPASSEMDDGPSRRDHFMSGFASGRNFGNRDA
 GECNKRNDTSTMGGFVVGKSFNGRGSNSRFEDGDSSGFWRGYRDGNNSEASGPYRRGGRGSFRGCRGGF
 GLGSPNNDLDPDECMQRTGGLFGSRRPVL SGTNGDTSQSRSVSGSERGGYKGLNEEVITGSGKNSWKSE
 AEGGESDQTQGPVITYIPPPPEDEDSIFAHYQTGINFDKYDTILVEVSGHDAPPAILTFEEANLCQTLN
 NNIKAGYTKLTPVQKYSIPIILAGRDLMACAQTGSGKTA AFLLPILAHMMHDGITASRFKELQEPECII
 VAPTRELVNQIYLEARKFSFGTCVRAVVIYGGTQLGHSIRQIVQGCNILCATPGRMLDIIIGKEIKLQKI
 KYLVLDEADRMLDMFGPEMKKLI SCPGMPSKEQRQTLMSATFPEEIQLAAEFLKSNLYFVAVGQVGG
 ACRDVQQTVLVQVQFSKREKLV EILRNIGDERTMVFVETKKKADF IATFLCQEKISTTSIHGDREQRERE
 QALGDFRFGKCPVLVATSVAARGLDIENVQHVINFDPSTIDEYVHRIGRTGRCGNTGRAISFFDLESND
 HLAQPLVKVLTDAQQDVPALWEEIAFSTYIPGFSGSTRGNVFAVDTRKKGKSTLNTAGFSSSQAPNPVDD
 ESWD


TRTRPLE - GFP Tag - V


Restriction Sites:

SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:


ACCN: NM_001166533

ORF Size: 2112 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001166533.1](#), [NP_001160005.1](#)

RefSeq Size: 2466 bp

RefSeq ORF: 2115 bp

Locus ID: 54514

UniProt ID: [Q9NQI0](#)

Cytogenetics: 5q11.2

Gene Summary: DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases. They are implicated in a number of cellular processes involving alteration of RNA secondary structure such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Based on their distribution patterns, some members of this family are believed to be involved in embryogenesis, spermatogenesis, and cellular growth and division. This gene encodes a DEAD box protein, which is a homolog of VASA proteins in *Drosophila* and several other species. The gene is specifically expressed in the germ cell lineage in both sexes and functions in germ cell development. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2009]