

Product datasheet for **RG230465**

LDB3 (NM_001171610) Human Tagged ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | LDB3 (NM_001171610) Human Tagged ORF Clone |
| Tag: | TurboGFP |
| Symbol: | LDB3 |
| Synonyms: | CMD1C; CMH24; CMPD3; CYPHER; LDB3Z1; LDB3Z4; LVNC3; MFM4; ORACLE; PDLIM6; ZASP |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-AC-GFP (PS100010) |
| E. coli Selection: | Ampicillin (100 ug/mL) |



[View online »](#)

ORF Nucleotide
Sequence:

>RG230465 representing NM_001171610
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTCTTACAGTGTGACCCTGACTGGGCCCGGCCCTGGGGCTTCGCTGCAGGGGGCAAGACTTCA
ACATGCCCTCACTATCTCCCGGATCACACCAGGCAGCAAGGCAGCCAGTCCCAGCTCAGCCAGGGTGA
CCTCGTGGTGGCCATTGACGGCGTCAACACAGACACCATGACCCACCTGGAAGCCAGAACAAGATCAAG
TCTGCCAGCTACAATTGAGCCTCACCTGCAGAAATCAAAGCGTCCCATTCCCATCTCCACGACAGCAC
CTCCAGTCCAGACCCCTCTGCCGGTGATCCCTACCAGAAGGACCCCGCTCTGGACACGAACGGCAGCCT
GGTGGCACCCAGCCAGCCCTGAGGGCAGGGCCAGCCAGGCACCCAGGCACCCCGGAGCTCAGGCC
ACCTTTAGCCCTGCCTTCTCCGGCCCTCCGCCTTCTCCTCACTCGCCGAGGCCTCTGACCTGGCCCTC
CGCGGGCCAGCCTGAGGGCCAAGACCAGCCAGAGGGGGCCGGGACCTACTCGGCCAAAAGCCCTGCC
GGGCTCGAGCCAGCCGAGGCAATAAACAACCCATTGGCTGTACTCGGCAGAGACCCTGAGGGAGATG
GCTCAGATGTACCAGATGAGCCTCCGAGGGAAGGCCTCGGGTGTGGACTCCAGGAGGCCCGACTACC
AGGAACGCTTCAACCCAGTGCCCTGAAGGACTCGGCCCTGTCCACCCACAAGCCATCGAGGTGAAGGG
GCTGGGCGGAAGGCCACCATCATCCATGCGCAGTACAACACGCCCATCAGCATGTATCCAGGATGCC
ATCATGGATGCCATCGCTGGGCAGGCCAAGCCAAAGGCAGTGACTTCAGTGGGAGCCTCCCTATTAAGG
ACCTTGCCGTAGACAGCGCCTCTCCCGTCTACCAGGCTGTGATTAAGAGCCAGAACAAGCCAGAAGATGA
GGCTGACGAGTGGGCACGCCGTTCTCCAACTGCAGTCTCGTCTCTCCGCATCTGGCCAGATGACG
GGGACAGAATTCATGCAAGACCCGTGATGAAGAAGCTCTGCGAAGGTCAAGGCCCCAGGCCTTCTCTACA
GCCCGCAGTGGCCGCTTCTCAGCACCTGCCACCCACACCAGCTACAGTGAGGGGCCCGCCCTTCTTACA
ACCCAAGCCCGGGTGTCAACACTGCCAGCATCCGGCCTTCTGTCTACCAGCCAGTGCCTGCATCTACC
TACAGCCCGTCCCGAGGGCCAATTACAGTCCCACTCCCTACACCCCTCCCTGCCCTGCCTACACCC
CCTCCCTGCCCTGCCTACACCCCTCACCTGTCCACCTACACTCCATCCCGAGCACCGCCTATAC
CCCCTCACCTGCCCCAACTATAACCTGCACCCTCGGTGGCCTACAGCGGGGGCCCTGCGGAGCCTGCC
AGCCGTCCACCCTGGGTGACAGATGATAGCTTCTCCAGAAGTTTGCCCCGGGCAAGAGCACCACTCCA
TCAGCAAGCAGACCCTGCCCGGGGAGGCCAGCCTACACCCAGCGGGTCTCAGGTGCCACCACTTGC
CAGGGGACCGTCCAGAGGGCTGAGCGATTCCAGCCAGCAGCCGACTCCACTCTGCGGTCACTGCAAC
AATGTCATCCGGGCCCCATTTCTGGTAGCCATGGGCCGTTCTTGGCACCTGAAGAGTTCACCTGTGCCT
ACTGCAAGACTTCCCTGGCAGATGTGTGCTTTGTGGAAGAGCAGAACAAGTTTACTGTGAGCGATGTTA
TGAGCAATTCTTTGCCCGCTGTGTGCCAAGTGAACACCAAAATTATGGGGGAAGTAATGCATGCCTTG
AGACAGACATGGCACACCACCTGCTTCTGTGTGCGGCCTGCAAGAAGCCTTTTGGGAACAGCCTTTCC
ACATGGAAGACGGGGAGCCCTACTGCGAGAAAGACTACATCAATCTGTTACAGACCAAGTGCCATGGCTG
CGATTTCCCGTGGAGGCTGGCGACAAGTTTATCGAAGCCCTGGGCCACACTTGGCACGACACCTGCTTC
ATTTGCGCAGTCTGCCATGTGAATCTGGAGGGCAGCCGTTCTACTCCAAGAAGGACAGACCCCTGTGCA
AGAAGCACGCACACCATCAACTTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG230465 representing NM_001171610
 Red=Cloning site Green=Tags(s)

```

MSYSVTLTGPWPWFRLQGGKDFNMPLTISRITPGSKAAQSQLSQGDLVVAIDGVNTDTMTHLEAQNKIK
SASYNLSLTLQSKRPIPISTTAPPVQTPLPVIHPHQDPALDTNGSLVAPSPSPEARASPGTPTPELRP
TFSPAFSRPSAFSSLAESDPGPPRASLRKTSPEGARDLLGPKALPGSSQPRQYNNPIGLYSAETLREM
AQMYQMSLRGKASGVGLPGGADYQERFNPSALKDSALSTHKPIEVKGLGGKATIIHAQYNTPIISMYSQDA
IMDAIAGQAQAQGSDFSGSLPIKDLAVDSASPVYQAVIKSQNKPEDEADEWARRSSNLQSRFRILAQMT
GTEFMQDPDEEALRRSRPQASSYSPAVAASSAPATHTSYSEGPAAPAPKPRVVTASIRPSVYQVPVAST
YSPSPGANYSPTPYTPSPAPAYTPSPAPAYTPSPVPTYTPSPAPAYTPSPAPNYNPAPSVAYSGGPAEPA
SRPPWVTDDSFQKFAPGKSTTSISKQTLPRGGPAYTPAGQVPLARGTVQRAERFPASSRTPLCGHNC
NVIRGPFVLVAMGRSWHPEEFTCA YCKTSLADVCFVEEQNNVYCERCYEQFFAPLCAKCN TKIMGEVMHAL
RQTWHHTCFVCAACKKPF GNSLFHMEDGEPYCEKDYINL FSTKCHGCDFPVEAGDKFIEALGHTWHHTCF
ICAVCHVNLEGQPFYSKKDRPLCKKHAHTINL
  
```

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_001171610

ORF Size: 2196 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001171610.2](#)

RefSeq Size: 5436 bp

RefSeq ORF: 2199 bp

Locus ID: 11155

UniProt ID: [O75112](#)

Cytogenetics: 10q23.2

Gene Summary: This gene encodes a PDZ domain-containing protein. PDZ motifs are modular protein-protein interaction domains consisting of 80-120 amino acid residues. PDZ domain-containing proteins interact with each other in cytoskeletal assembly or with other proteins involved in targeting and clustering of membrane proteins. The protein encoded by this gene interacts with alpha-actinin-2 through its N-terminal PDZ domain and with protein kinase C via its C-terminal LIM domains. The LIM domain is a cysteine-rich motif defined by 50-60 amino acids containing two zinc-binding modules. This protein also interacts with all three members of the myozenin family. Mutations in this gene have been associated with myofibrillar myopathy and dilated cardiomyopathy. Alternatively spliced transcript variants encoding different isoforms have been identified; all isoforms have N-terminal PDZ domains while only longer isoforms (1, 2 and 5) have C-terminal LIM domains. [provided by RefSeq, Jan 2010]