

Product datasheet for **RG230489**

AMPD1 (NM_001172626) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AMPD1 (NM_001172626) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	AMPD1
Synonyms:	MAD; MADA; MMDD
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG230489 representing NM_001172626
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAATGTCAGGATATTTTATAGTGTCACTCAGTCACTCACCCACAGTCTCTCTCTCTCTTTTCTACTGTG
CTATCCTAGAATCAAGGATTTTCAGCAACAATGCCTCTGTTCAAACCTCCAGAAATTGATGATGCAATGCG
CAACTTTGCTGAAAAAGTGTTCCTCTGAAGTCAAAGATGAAGGAGGTCGTCAGGAGATTTCCCCCTTT
GATGTGGATGAGATCTGTCGATTTCTCATCATGAGATGCAAGCACACATATTCCATCTGGAGACTCTGT
CCACCTCCACAGAAGCCAGGAGAAAAAGCGTTTCAAGGACGGAAGACTGTTAATTTGTCCATTCCACT
AAGTGAACATCTTCCACCAAACGTCCACATTGATGAATACATTTCTCATCTCCAACCTACCAGACC
GTGCTGATTTTTCAGAGAGTGCAGATTACTGGTGACTATGCCTCTGGGTTACAGTTGAAGATTTTGAAA
TTGTTTCAAAGGTCGTATCGGCCACTATGCATACGTGAGAAAACATGCAGAAGTCGTTTCAGAGGTT
CCCTAAAACCCCTTCAAATACTTGCGGAACATTGATGGTGAGGCTTGGGTAGCAAATGAGAGCTTCTAT
CCAGTCTTTACTCTCTGTGAAGAAGGAGAGGACCCCTTCGAACAGACAACCTTCTGAAAACCTGG
GCTATCACCTCAAATGAAGGACGGTGTAGTTACGTCTATCCTAATGAAGCAGCAGTCAGCAAAGATGA
GCCTAAGCCACTTCTTACCCAACTCTGGACACCTTCTTAGACGATATGAATTTTTACTTGCTTTAATT
GCTCAAGGACCTGTTAAGACCTATACCCACCGCGCCTGAAGTTCCTCTCCTCCAAGTCCAGGTCATC
AGATGCTTAACGAGATGGACGAGTTAAAGGAGCTGAAAAACAACCCCAACCGAGATTTTTATAACTGCAG
GAAGGTGGACACCCATATCCATGCAGCCGCTTGCATGAACCAGAAAACATCTGCTGCGTTTTATAAGAAA
TCTTACCAAATGATGCTGACAGAGTGGTCTATAGCACAAAGAGAAGAATCTGACCCTAAAGGAACCTT
TTGCTAAATTAATAATGCATCCTTATGACCTGACTGTTGATTCTCTGGATGTTGATGCTGACGCCAGAC
CTTCCAGCGTTTTGATAAGTTCAATGACAAATATAATCCTGTAGGAGCAAGTGAGCTACGGGACCTCTAC
TTGAAGACAGACAATTACATTAATGGGGAATTTTTGCCACTATCATCAAGGAGGTAGGTGCGGACCTGG
TGGAGGCCAAGTACCAGCATGCTGAGCCCCGCTGTCCATCTATGGCCGAGTCCTGATGAGTGGAGCAA
ACTCTCTCTGCTGCTGCAATCGCATCCACTGCCCAACATGACATGGATGATCCAGGTTCCAGG
ATCTATGATGTGTTCCGTTCCAAGAAATTCCTTCCACATTTTGGAAAAATGCTGGAGAATATTTTCATGC
CAGTGTTTGGAGCCACCATCAACCCCAAGGCTGACCCAGAAGTCAAGTGTCTTCTCAAGCATATCACTGG
CTTTGACAGTGTGGATGATGAGTCCAAACACAGTGGCCACATGTTCTCCTCCAAGAGTCCAAGCCCCAG
GAGTGGACATTGAAAAGAATCCATCTTACACTTACTATGCCTACTACATGTATGCAAACATCATGGTGC
TCAACAGCCTGAGAAAGAACGAGGCATGAATACGTTTTCTGTTCCGACCTCACTGTGGAGAAGCTGGAGC
CCTCACCCATCTCATGACAGATTTCATGATAGCAGATGATATCTCTCATGGCCTAAATTTAAAAAAGAGT
CCCGTGCTACAGTACTGTTTTCTTAGCCAAATTTCCATCGCCATGTCACCACTAAGTAACAATAGCC
TATTTCTAGAGTATGCCAAAAATCCTTTTTTGATTTCCTTCAGAAAGGGCTAATGATCTCACTGTCTAC
AGATGACCCAATGCAATTCACCTTTACCAAGGAGCCCCTAATGGAAGAATATGCTATTGCTGCACAAGTC
TTCAAGCTGAGCACCTGTGATATGTGCGAAGTGGCAAGGAACAGTGTCTTGCAGTGTGGAATTTCTCATG
AGGAGAAAGTAAAGTTTCTGGGCGACAATTACCTTGAGGAAGGCCCTGCTGAAAATGATATCCGGAGGAC
AAATGTAGCCCAAATCCGCATGGCCTATCGCTATGAAACCTGGTGTATGAACTCAATTTAATTGCTGAG
GGTCTTAAATCAACAGAA

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG230489 representing NM_001172626
 Red=Cloning site Green=Tags(s)

MNVRIFYSVSQSPHLSLLLFYCAILESRI SATMPLFKLPEIDDAMRNFAEKVFASEVKDEGGRQEI SPF
 DVDEICPISHHEMQAHIFHLET LSTSTEARRKKRFQGRKTVNLSIPLSETSSSTKLSHIDEYISSSPTYQT
 VPDFQRVQITGDYASGVTVEDFEIVCKGLYRALCIREKYMQKSFQRFPKTPSKYLRNIDGEAWVANESFY
 PVFTPPVKKGEDPFRDNL PENLGYHLKMKDGVVYVYPNEAAVSKDEPKPLPYPNLDTFLLDMNFL LALI
 AQGPVKTYTHRRLLKFLSSKFQVHQMLNEMDELKELKNNPHRDFYNCRKVDTHIAAACMNQKHL LRFIKK
 SYQIDADR VVYSTKEKNLTLKELFAKLKMPYDLTVDSL DVHAGRQTFQRFDKFNDKYNPVGASELRDLY
 LKTDNYINGEYFATI IKEVGADLVEAKYQHAEPRLSIYGRSPDEWSKLSWFVCNRIHCPNMTWMIQVPR
 IYDVFRSKNFLPHFGKMLENIFMPVFEATINPQADPELSVFLKHITGFDSVDDSKHSGHMFSSKSPKPQ
 EWTLEKNPSYTYAYMYANIMVLSLNRKERGMNTFLFRPHCGEAGAL THLMTAFMIADDISHGLNLKKS
 PVLQYLFFLAQIPIAMSPLSNNSLFLEYAKNPFDFLQGLMISLSTDDPMQHFHTKEPLMEEYIAAAQV
 FKLSTCDMCEVARNSVLQCGISHEEKVKFLGDNYLEEGPAGNDIRRTNVAQIRMAYRYETWCYELN LIAE
 GLKSTE

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001172626


ORF Size: 2328 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001172626.1</u> , <u>NP_001166097.1</u>
RefSeq Size:	2414 bp
RefSeq ORF:	2232 bp
Locus ID:	270
UniProt ID:	<u>P23109</u>
Cytogenetics:	1p13.2
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Purine metabolism
Gene Summary:	Adenosine monophosphate deaminase 1 catalyzes the deamination of AMP to IMP in skeletal muscle and plays an important role in the purine nucleotide cycle. Two other genes have been identified, AMPD2 and AMPD3, for the liver- and erythrocyte-specific isoforms, respectively. Deficiency of the muscle-specific enzyme is apparently a common cause of exercise-induced myopathy and probably the most common cause of metabolic myopathy in the human. Alternatively spliced transcript variants encoding different isoforms have been identified in this gene.[provided by RefSeq, Feb 2010]

Product images:



Circular map for RG230489