

Product datasheet for **RG230504**

ZFX (NM_001178085) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZFX (NM_001178085) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ZFX
Synonyms:	ZNF926
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG230504 representing NM_001178085
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGATGAAGATGGGCTTGAATTACAACAAGAGCCAAACTCATTTTTTGTGCAACAGGAGCTGATGGTA
 CACACATGGATGGTGATCAAATTGTTGTGGAAGTACAAGAACTGTTTTGTTTCAGATGTTGTGGATTC
 AGACATAACTGTGCATAACTTTGTTCTGATGACCCAGATTCAGTTGTAATCCAAGATGTTATTGAGGAC
 GTTGTTATAGAAGATGTTCAAGTCCAGATATCATGGAAGAAGCAGATGTGTCTGAAACGGTCATCATT
 CTGAGCAAGTGGTGGACTCAGATGTAAGTGAAGAAGTTTCTTTAGCACATTGCACAGTCCCAGATGATGT
 TTTAGCTTCTGACATTACTTCAGCCTCAATGTCTATGCCAGAACACGCTTTGACGGGTGATTCTATACAT
 GTGTCTGACGTTGGACATGTTGGACATGTTGGACATGTTGAACATGTGGTTCATGATAGTGTAGTGAAG
 CAGAAATTGTCACTGATCCTCTGACTACCGACGTAGTTTCAGAAGAAGTATTGGTAGCAGACTGTGCCTC
 TGAAGCAGTCATAGATGCCAATGGGATCCCTGTGGACCAGCAGGATGATGACAAAGGCAACTGTGAGGAC
 TACCTTATGATTTCTTGGATGATGCTGGCAAAATAGAACACGATGGTTCCTTCTGGAATGACCATGGACA
 CAGAGTCGGAAATTGATCCTTGTAAAGTGGATGGCACTTGCCCTGAGGTCATCAAGGTGTACATTTTTAA
 AGCTGACCCTGGAGAAGTACTTAGGTGGAAGTGTAGACATTGTGGAGAGTGAGCCTGAGAATGATCAT
 GGAGTTGAACTGCTTGTATCAGAACAGCAGTATTCGTGTTCCAGGGAAAAGATGGTTTATGACTGTCA
 ATGACTCTCAGCCAGAAGATGAAGATTTAAATGTTGCTGAAATCGCTGACGAAGTTTATATGGAAGTGT
 CGTAGGAGAGGAGGATGCTGCAGCAGCAGCGCCAGCCGCCGCTGCACGAGCAGCAATGGATGACAA
 TAAATGCAAAACCTTCATGCCGATTGCATGGCAGCAGCTTATGGTAATAATTCTGATGGAATTGAAAA
 GGAAATGGCACTGCAAGTGCCTCTTGCACATAGATGAGTCTGCTGGCCTCGCAGACTGGCTAAACAAAA
 ACCAAAGAAAAGGAGAAGACCTGATCCAGCAGTACCAACAGCAATAATTATTGGCCCTGATGGACAT
 CCTTTGACTGTCTATCCTTGCATGATTTGTGGGAAGAAGTTTAAAGTCGAGAGTTTTTTGAAAAGGCACA
 TAAAAACCATCCGAACACCTTGCCAAGAAGAAATACCGCTGTACTGACTGTGATTACACTACCAACAA
 GAAGATAAGTTTACACAACCACCTGGAGAGCCACAAGCTGACCAGCAAGGCAGAGAAGGCCATTGAATGC
 GATGAGTGTGGGAAGCATTCTCTCATGCAGGGGCTTTGTTTACTCACAAAATGGTGCATAAAGAAAAAG
 GAGCCAACAAAATGCACAAGTGTAAATCTGTGAATACGAGACAGCTGAACAAGGTTATTGAATCGCCA
 CCTTTGGCAGTCCACAGCAAGAACTTTCTCATATTTGTGTGGAGTGTGGTAAGGGTTTTCTGCACCCG
 TCAGAGCTCAAAAAGCACATGAGAAATCCATACTGGGGAGAAGCCGTACCAATGCCAGTACTGCGAATATA
 GGTCTGCAGACTTTCTAACTTGAACACGATGTCAAAACTAAGCATAGTAAAGAGATGCCATTCAGTGT
 TGACATTTGTCTTCTGACTTTCTCGGATACCAAGAGGTGCAGCAACATGCTCTTATCCACCAAGAAAGC
 AAAACACACCAGTGTTCATTGCGACCACAAGAGTTCGAACTCAAGTATTTGAAACGACACATAATTT
 CAGTTCACACGAAAGACTACCCCATAAAGTGTGACATGTGTGATAAAGGCTTTCACAGGCCTTCAGAACT
 CAAGAAACACGTGGCTGCCACAAAGGGCAAAAAATGCACCAGTGTAGACATTGTGACTTTAAGATTGCA
 GATCCATTTGTTCTAAGTCGCCATTTCTCTCAGTTCACACAAAGGATCTTCCATTTAGGTGCAAGAGAT
 GTAGAAAGGGATTTAGGCAACAGAGTGTGCTTAAAAAGCATATGAAGACACACAGTGGCAGGAAAGTGT
 TCAGTGTGAGTACTGTGAGTATAGCACTACAGATGCCTCAGGCTTTAAACGGCACGTTATTTCCATTAC
 ACGAAAGACTATCCTCACCGGTGTGAGTACTGCAAGAAAGGCTTCCGAAGACCTTCAGAAAAGAACCCAGC
 ACATAATGCGACATCATAAAGAAGTTGGCTGCC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG230504 representing NM_001178085
 Red=Cloning site Green=Tags(s)


```

MDEDGLELQQEPNSFFDATGADGTHMDGDQIVVEVQETVVFVSDVVSDITVHNFVDDPDSVVIQDVIED
VVIEDVQCPDIMEEADVSETVIIPEQVLDSVTEEVSLAHCTVPDDVLASDITSASMSMPEHVLTGDSIH
VSDVGHVGHVGHVEHVHDSVVEAEIVTDPLTTDVVSEEVLVADCASEAVIDANGIPVDQDDDKGNCE
YLMISLDDAGKIEHDGSSGMTDTESEIDPCKVDGTCPEVIKVIYIFKADPGEDDLGGTVDIVSEPENDH
GVELLDQNSSIRVPREKMYMTVNDSPQPEDEDLNVAEIADEVYMEVIVGEEDAAAAAAAAAVHEQQMDDN
EIKTFMPIAWAAAYGNNSDGIENRNGTASALLHIDESAGLGRLAKQKPKKRRRPDSRQYQTAIIGPDGH
PLTVYPCMICGKKFKSRGFLKRHMKNHPEHLAKKKYRCTDCDYTTNKKISLHNHLESHKLTSKAEKAI
DECCKHFSHAGALFTHKMVHKEKGANKMHKCKFCYETAEQGLLRHLLAVHKNFPHICVECGKGRHP
SELKKHMRHTGEKPYQCQYCEYRSADSSNLKTHVKTKHSKEMPFKCDICLLTFSDTKEVQHALIHQES
KTHQCLHCDHKSNSDLKRHIISVHTKDYPHKCDMCDKGFHRPSELKKHVAHKGKKMHQCRHCDFKIA
DPFVLSRHILSVHTKDLPRCKRCRKGFRQQSELKKHMKTHSGRKVYQCEYCEYSTTDASGFKRHVISIH
TKDYPHRCEYCKGFRPSEKNQHIMRHHKEVGLP
  
```

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI

Plasmid Map:



ACCN: NM_001178085

ORF Size: 2415 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001178085.1](#), [NP_001171556.1](#)

RefSeq Size: 7740 bp

RefSeq ORF: 2418 bp

Locus ID: 7543

UniProt ID: [P17010](#)

Cytogenetics: Xp22.11

Protein Families: Transcription Factors

Gene Summary: This gene on the X chromosome is structurally similar to a related gene on the Y chromosome. It encodes a member of the krueppel C2H2-type zinc-finger protein family. The full-length protein contains an acidic transcriptional activation domain (AD), a nuclear localization sequence (NLS) and a DNA binding domain (DBD) consisting of 13 C2H2-type zinc fingers. Studies in mouse embryonic and adult hematopoietic stem cells showed that this gene was required as a transcriptional regulator for self-renewal of both stem cell types, but it was dispensable for growth and differentiation of their progeny. Multiple alternatively spliced transcript variants encoding different isoforms have been identified, but the full-length nature of some variants has not been determined. [provided by RefSeq, May 2010]