

Product datasheet for **RG230686**

PARD3 (NM_001184790) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PARD3 (NM_001184790) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	PARD3
Synonyms:	ASIP; Baz; PAR3; PAR3alpha; PARD-3; PARD3A; PPP1R118; SE2-5L16; SE2-5L1T1; SE2-5T2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG230686 representing NM_001184790 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAAGTGACCGTGTGCTTCGGACGGACCCGGGTGGTCGTGCCGTGCGGGGACGGCCACATGAAAGTTT
TCAGCCTCATCCAGCAGGCGGTGACCCGCTACCGGAAGGCCATCGCCAAGGATCCAACTACTGGATACA
GGTGCATCGTTGGAACATGGAGATGGAGGAATACTAGACCTTGATGACATTCTTTGTGATGTAGCAGAC
GATAAAGACAGACTGGTAGCAGTGTGGATGAGCAGGATCCACATCACGGAGGTGATGGCACCAGTGCCA
GTTCCACGGGTACCCAGAGCCCAGAGATATTTGGTAGTGAGCTTGGCACCAACAATGTCTCAGCCTTTCA
GCCTTACCAAGCAACAAGTAAAATTGAGGTCACACCTTCAGTCCTTCGAGCAAATATGCCTTTCATGTT
CGACGCAGTAGTGACCCAGCTCTAATTGGCCTCTCCACTTCTGTGAGTATAGTAATTTTTCTCTGAAG
AGCCTTCAAGGAAAAATCCCACACGCTGGTCAACAACAGCTGGCTTCTCAAGCAGAACACTGCTGGGAG
TCCTAAAACCTGCGACAGGAAGGATGAGGATGGGACAGAAGAGGATAACAGTCGTGTTGAACCTGTTGGA
CATGCTGACACGGGTTGGAGCATATACCAACTTTTCTCTGGATGATATGGTAAAGCTCGTAGAAGTCC
CCAACGATGGAGGGCCTCTGGGAATCCATGTAGTGCCCTTCAGTGCTCGAGGCGGCAGAACCCCTGGGTT
ATTAGTAAAACGATTGGAGAAAGGTGGTAAAGCTGAACATGAAAACTTTTTCTGAGAAATGATTGCATT
GTCAGGATTAATGATGGCGACCTTCGAAATAGAAGATTTGAACAAGCACAAACATATGTTTCGCCAAGCCA
TGCGTACACCCATCATTGGTTCCATGTGGTTCCCTGCAGCAAATAAAGAGCAGTATGAACAACATATCCCA
AAGTGAGAAGAACAATTACTATTCAAGCCGTTTTAGCCCTGCAGCCAGTATATTGACAACAGGAGTGTG
AACAGTGCAGGGCTTACACCGTGCAGAGAGCACCCCGACTGAACCACCCGCTGAGCAGATAGACTCTC
ACTCAAGACTACCTCATAGCGCACACCCCTCGGGAAAACCACCATCCGCTCCAGCCTCGGCACCTCAGAA
TGTATTTAGTACGACTGTAAGCAGTGGTTATAACACCAAAAAAATAGGCAAGAGGCTTAATATCCAGCTT
AAGAAAGTACAGAAGGTTGGGATTGAGCATCACTTCCAGAGATGTAACAATAGGTGGCTCAGCTCCAA
TCTATGTGAAAAACATTCTCCCCGGGGGGCGCCATTAGGATGGCCGACTTAAGGCAGGAGACAGACT
TATAGAGTAAATGGAGTAGATTTAGTGGCAAATCCCAAGGAAGTTGTTTCGCTGTTGAGAAGCACC



[View online »](#)

AAGATGGAAGGAAGTGTGAGCCTTCTGGTCTTTCGCCAGGAAGACGCCTTCCACCCAAGGGAAGTAAAAG
CAGAAGATGAGGATATTGTTCTTACACCTGATGGCACCAGGGAATTTCTGACATTTGAAGTCCCCTTAA
TGATTCAGGATCTGCAGGCCTTGGTGTCAAGTGTCAAAGGTAACCGGTCAAAGAGAACCCAGCAGATTTG
GGAATCTTTGTCAAGTCCATTATTAATGGAGGAGCAGCATCTAAAGATGGAAGGCTTCGGGTGAATGATC
AACTGATAGCAGTAAATGGAGAATCCCTGTTGGCAAGACAAACCAAGATGCCATGGAACCCCTAAGAAG
GTCTATGTCTACTGAAGCAATAAACGAGGAATGATCCAGCTTATTGTTGCAAGGAGAATAAGCAAGTGC
AATGAGCTGAAGTCACTGGGAGCCCCCTGGACCTGAGCTGCCATTGAAACAGCGTTGGATGATAGAG
AACGAAGAATTTCCATTCCCTCTACAGTGGGATTGAGGGGCTTGATGAATCGCCAGCAGAAATGCTGC
CCTCAGTAGGATAATGGGTAATACCAGCTGTCCCTACAGTGAATATGCCCAAGATGACACTGTCATT
ATAGAAGATGACAGGTTGCCAGTCTTCTCCACATCTCTGACCAGTCTTCCAGCTCCCATGATG
ATGTGGGTTTGTGACGGCAGATGCTGGTACTTGGCCAAAGGCTGCAATCAGTGATTAGCCGACTGCTC
TTTGAGTCCAGATGTTGATCCAGTCTTCTTCAACGAGAAGGATTTGGACGTGAGATAGCTGACGAG
ACTAACTCAATACAGTGGATGACCAGAAAGCAGGTTCTCCAGCAGAGATGTGGTCTTCCCTGGGTC
TGAAGAAGTCAAGCTCGTTGGAGAGTCTGCAGACCGCAGTTGCCGAGGTGACTTTGAATGGGGATATTCC
TTTCCATCGTCCACGGCCGCGGATAATCAGAGGCAGGGGATGCAATGAGAGCTTCAGAGCTGCCATCGAC
AAATCTTATGATAAACCCGCGGTAGATGATGATGAAGGCATGGAGACCTTGAAGAAGACACAGAAG
AAAGTTCAAGATCAGGGAGAGAGTCTGTATCCACAGCCAGTGATCAGCCTTCCCCTCTCTGGAGAGACA
AATGAATGGAACCAAGAGAAAGGTGATAAGACTGATAGAAAAAGGATAAACTGGAAAAGAAAAGAG
AAAGATAGAGATAAGGAGAAGGATAAAATGAAAGCCAAGAAGGGAATGCTGAAGGGCTTGGGAGACATGT
TCAGGTTTGGCAAACATCGAAAAGATGACAAGATTGAGAAAACGGGTAATAAAAAACAGGAATCCTT
TACATCAGAAGAGGAGAGGATACGAATGAAGCAGGAGCAGGAGAGGATTCAGCCAAAACTCGAGAATTT
AGGGAACGACAAGCTCGAGAGCGTGACTATGCTGAAATTCAGATTTTCATCGGACATTTGGCTGTGATG
ATGAGTTAATGTATGGGGAGTTTCTTCTTATGAAGTTCCATGGCTCTCAACGCTAGACCTCAGAGCCC
ACGAGAAGGGCATATGATGGATGCTTTGTATGCCCAAGTCAAGAAGCCGCGGAATTCAAAACCTCACCT
GTAGACAGTAACAGATCAACTCCTAGCAATCATGATCGGATACAGCGTCTGAGGCAAGAATTTAGCAAG
CAAAGCAAGATGAAGATGTAGAAGATCGTCGGCGGACCTATAGTTTTGAGCAACCCTGGCCGAACGCACG
GCCGGCGACGCAGAGCGGGCGACACTCGGTGTCCGTGGAGGTGCAGATGCAGCGGCAGCGGAGGAGGAG
CGCGAGAGCTCCCAGCAGGCCAGCGCCAGTACAGCTCTCTGCCTCGGCAAAGCAGGAAAAATGCCAGCT
CGGTCTCCAGGACTCTTGGGAGCAGAATACTCCCTGGGAAGGCTTCCAGAGTGCCAAAGAGAAACCC
CAGGTACTCCAGTACCAAGGCTCCAGGAACGGTACCTGGGAGGACATGGCTTCAACGCCAGGGTCATG
CTGAAACTCAGGAGCTCCTTCGCCAGGAACAGAGGCGGAAGGAGCAGCAGATGAAGAAGCAGCCTCCTT
CCGAGGGGCCAGCAACTATGACTCGTATAAGAAAGTCCAGGACCCAGTTACGCCCTCCCAAGGGGCC
CTTCCGGCAAGATGTGCCCCCTCCCCTTCTCAGGTTGCGAGGCTGAACAGACTTCAGACTCCTGAGAAA
GGGAGGCCCTTCTATTCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG230686 representing NM_001184790
 Red=Cloning site Green=Tags(s)

MKVTVCFGRTRVVVPCGDGHMKVFSLIQQAVTRYKAIKDPNYWIQVHRLEHGDGGILDLDLDDILCDVAD
 DKDRLVAVFDEQDPHHGGDGTASSTGTQSP E I F G S E L G T N N V S A F Q P Y Q A T S E I E V T P S V L R A N M P L H V
 R R S S D P A L I G L S T S V S D S N F S S E E P S R K N P T R W S T T A G F L K Q N T A G S P K T C D R K D E D G T E E D N S R V E P V G
 H A D T G L E H I P N F S L D D M V K L V E V P N D G G P L G I H V V P F S A R G G R T L G L L V K R L E K G G K A E H E N L F R E N D C I
 V R I N D G D L R N R R F E Q A Q H M F R Q A M R T P I I W F H V V P A A N K E Q Y E Q L S Q S E K N N Y S S R F S P D S Q Y I D N R S V
 N S A G L H T V Q R A P R L N H P P E Q I D S H S R L P H S A H P S G K P P S A P A S A P Q N V F S T T V S S G Y N T K K I G K R L N I Q L
 K K G T E G L G F S I T S R D V T I G G S A P I Y V K N I L P R G A A I Q D G R L K A G D R L I E V N G V D L V G K S Q E E V V S L L R S T
 K M E G T V S L L V F R Q E D A F H P R E L K A E D E D I V L T P D G T R E F L T F E V P L N D S G S A G L G V S V K G N R S K E N H A D L
 G I F V K S I I N G G A A S K D G R L R V N D Q L I A V N G E S L L G K T N Q D A M E T L R R S M S T E G N K R G M I Q L I V A R R I S K C
 N E L K S P G S P P G P E L P I E T A L D D R R E R I S H S L Y S G I E G L D E S P S R N A A L S R I M G K Y Q L S P T V N M P Q D D T V I
 I E D D R L P V L P P H L S D Q S S S S H D D V G F V T A D A G T W A K A A I S D S A D C S L S P D V D P V L A F Q R E G F G R Q I A D E
 T K L N T V D D Q K A G S P S R D V G P S L G L K K S S L E S L Q T A V A E V T L N G D I P F H R P R P R I I R G R G C N E S F R A A I D
 K S Y D K P A V D D D E G M E T L E E D T E E S S R S G R E S V S T A S D Q P S H S L E R Q M N G N Q E K G D K T D R K K D K T G K E K K
 K D R D K E K D K M K A K K G M L K G L G D M F R F G K H R K D D K I E K T G K I K I Q E S F T S E E E R I R M K Q E Q E R I Q A K T R E F
 R E R Q A R E R D Y A E I Q D F H R T F G C D D E L M Y G G V S S Y E G S M A L N A R P Q S P R E G H M M D A L Y A Q V K K P R N S K P S P
 V D S N R S T P S N H D R I Q R L R Q E F Q Q A K Q D E D V E D R R R T Y S F E Q P W P N A R P A T Q S G R H S V S V E V Q M Q R Q R Q E E
 R E S S Q Q A Q R Q Y S S L P R Q S R K N A S S V S Q D S W E Q N Y S P G E G F Q S A K E N P R Y S S Y Q G S R N G Y L G G H G F N A R V M
 L E T Q E L L R Q E Q R R K E Q M K K Q P P S E G P S N Y D S Y K K V Q D P S Y A P P K G P F R Q D V P P S P S Q V A R L N R L Q T P E K
 G R P F Y S


TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_001184790

ORF Size: 3798 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001184790.1](#), [NP_001171719.1](#)

RefSeq Size: 5743 bp

RefSeq ORF: 3801 bp

Locus ID: 56288

UniProt ID: [Q8TEW0](#)

Cytogenetics: 10p11.22-p11.21

Protein Pathways: Adherens junction, Chemokine signaling pathway, Endocytosis, Neuroactive ligand-receptor interaction, Tight junction

Gene Summary: This gene encodes a member of the PARD protein family. PARD family members interact with other PARD family members and other proteins; they affect asymmetrical cell division and direct polarized cell growth. Multiple alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Oct 2011]