

Product datasheet for **RG230699**

PARD3 (NM_001184787) Human Tagged ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | PARD3 (NM_001184787) Human Tagged ORF Clone |
| Tag: | TurboGFP |
| Symbol: | PARD3 |
| Synonyms: | ASIP; Baz; PAR3; PAR3alpha; PARD-3; PARD3A; PPP1R118; SE2-5L16; SE2-5L1T1; SE2-5T2 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-AC-GFP (PS100010) |
| E. coli Selection: | Ampicillin (100 ug/mL) |
| ORF Nucleotide Sequence: | >RG230699 representing NM_001184787 Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAAGTGACCGTGTGCTTCGGACGGACCCGGGTGGTCGTGCCGTGCGGGGACGGCCACATGAAAGTTT
TCAGCCTCATCCAGCAGGCGGTGACCCGCTACCGGAAGGCCATCGCCAAGGATCCAACTACTGGATACA
GGTGCATCGTTGGAACATGGAGATGGAGGAATACTAGACCTTGATGACATTCTTTGTGATGTAGCAGAC
GATAAAGACAGACTGGTAGCAGTGTGGATGAGCAGGATCCACATCACGGAGGTGATGGCACCAGTGCCA
GTTCCACGGGTACCCAGAGCCCAGAGATATTTGGTAGTGAGCTTGGCACCAACAATGTCTCAGCCTTTCA
GCCTTACCAAGCAACAAGTAAAATTGAGGTCACACCTTCAGTCCTTCGAGCAAATATGCCTTTCATGTT
CGACGCAGTAGTGACCCAGCTCTAATTGGCCTCTCCACTTCTGTGAGTATAGTAATTTTTCTCTGAAAG
AGCCTTCAAGGAAAAATCCCACACGCTGGTCAACAACAGCTGGCTTCTCAAGCAGAACACTGCTGGGAG
TCCTAAAACCTGCGACAGGAAGAAAGATGAAAACACAGAACCTCCCGCGGGATACTAGTAAGTGGTCT
AACCAATTTAGAGAGACAATGCTCGCTCGTCTGAGTGCCAGTCAACCAATGGTGGGCAAGTGGCTGG
AGAAACAAGAACAGGATGAGGATGGGACAGAAGAGGATAACAGTCGTGTTGAACCTGTTGGACATGCTGA
CACGGGTTTGGAGCATATACCAACTTTTCTCTGGATGATATGGTAAAGCTCGTAGAAGTCCCAACGAT
GGAGGGCTCTGGGAATCCATGTAGTGCCTTTTCAGTGCTCGAGGGCGCAGAACCTGGGGTTATTAGTAA
AACGATTGGAGAAAGGTGGTAAAGCTGAACATGAAAATCTTTTTCTGAGAAATGATTGCATTGTCAGGAT
TAATGATGGCGACCTTCGAAATAGAAGATTTGAACAAGCACAAACATATGTTTCGCAAGCCATGCGTACA
CCCATCATTTGGTTCCATGTGGTTCCTGCAGCAAATAAAGAGCAGTATGAACAACATCCCAAAGTGA
AGAACAATTACTATTCAAGCCGTTTTAGCCCTGACAGCCAGTATATTGACAACAGGAGTGTGAACAGTGC
AGGGCTTACACGGTGCAGAGAGCACCCCGACTGAACCACCCGCTGAGCAGATAGACTCTCACTCAAGA
CTACCTCATAGCGCACACCCCTCGGAAAACCACCATCCGCTCCAGCCTCGGCACCTCAGAATGTATTTA
GTACGACTGTAAGCAGTGGTTATAACACCAAAAAAATAGGCAAGAGGCTTAATATCCAGCTTAAGAAAGG
TACAGAAGGTTTGGGATTCAGCATCACTCCAGAGATGTAACAATAGGTGGCTCAGTCCAATCTATGTG



[View online >](#)

AAAAACATTCTCCCCGGGGGGCGGCCATTTCAGGATGGCCGACTTAAGGCAGGAGACAGACTTATAGAGG
TAAATGGAGTAGATTTAGTGGGCAAATCCCAAGAGGAAGTTGTTTCGCTGTTGAGAAGCACCAAGATGGA
AGGAACTGTGAGCCTTCTGGTCTTTCGCCAGGAAGACGCCTTCCACCCAAGGGAAGTGAATGCAGAGCCA
AGCCAGATGCAGATTCAAAAAGAAACGAAAGCAGAAGATGAGGATATTGTTCTTACACCTGATGGCACCA
GGGAATTTCTGACATTTGAAGTCCCCTTAATGATTCAGGATCTGCAGGCCTTGGTGTGAGTGTCAAAGG
TAACCGGTCAAAGAGAACCACGCAGATTTGGGAATCTTGTCAAGTCCATTATTAATGGAGGAGCAGCA
TCTAAAGATGGAAGGCTTCGGGTGAATGATCAACTGATAGCAGTAAATGGAGAATCCCTGTTGGCAAGA
CAAACCAAGATGCCATGGAAACCCTAAGAAGGTCTATGTCTACTGAAGGCAATAAACGAGGAATGATCCA
GCTTATTGTTGCAAGGAGAATAAGCAAGTGCAATGAGCTGAAGTCACTGGGAGCCCCCTGGACCTGAG
CTGCCATTGAAACAGCGTTGGATGATAGAGAACGAAGAATTTCCATTCCCTCTACAGTGGGATTGAGG
GGCTTGATGAATCGCCAGCAGAAATGCTGCCCTCAGTAGGATAATGGGTGAGTCAGGTAATACCAGCT
GTCCCCTACAGTGAATATGCCCAAGATGACACTGTCATTATAGAAGATGACAGGTTGCCAGTGTCTCT
CCACATCTCTGACCAGTCTCTCCAGCTCCCATGATGATGTGGGTTTGTGACGGCAGATGCTGGTA
CTTGGGCCAAGGCTGCAATCAGTGATTCAGCCGACTGCTCTTGGTCCAGATGTTGATCCAGTCTTTGC
TTTTCAACGAGAAGGATTTGGACGTCAGAGTATGTCAGAAAAACGCACAAAGCAATTTTCAGATGCCAGT
CAATTGGATTTTCGTTAAACACGAAAATCAAAAAGCATGGATTTAGGTATAGCTGACGAGACTAACTCA
ATACAGTGGATGACCAGAAAGCAGGTTCTCCAGCAGAGATGTGGGTCTTCCCTGGGTCTGAAGAAGTC
AAGCTCGTTGGAGAGTCTGCAGACCGCAGTTGCCGAGGTGACTTTGAATGGGGATATTCCTTTCCATCGT
CCACGGCCCGGATAATCAGAGGCAGGGGATGCAATGAGAGCTTCAGAGCTGCCATCGACAAATCTTATG
ATAAACCCGCGGTAGATGATGATGATGAAGGATGGAGACCTTGAAGAAGACACAGAAGAAAGTTCAAG
ATCAGGGAGAGAGTCTGTATCCACAGCCAGTATCAGCCTTCCCCTCTCTGGAGAGACAAATGAATGGA
AACCAAGAGAAGGTGATAAGACTGATAGAAAAAGGATAAACTGAAAAAGAAAGAAGAAAGATAGAG
ATAAGGAGAAGGATAAAATGAAAGCCAAGAAGGGAATGCTGAAGGGCTTGGGAGACATGTTCAGGATTC
AGCCAAAACCTCGAGAATTTAGGGAACGACAAGCTCGAGAGCGTGACTATGCTGAAATTCAGATTTTCAT
CGGACATTTGGCTGTGATGATGAGTAAATGTATGGGGAGTTTCTTCTTATGAAGTTCCATGGCTCTCA
ACGCTAGACCTCAGAGCCACGAGAAGGCATATGATGGATGCTTTGTATGCCCAAGTCAAGAAGCCGCG
GAATTCAAAACCTCACCTGTAGACAGTAACAGATCAACTCCTAGCAATCATGATCGGATACAGCGTCTG
AGGCAAGAATTTAGCAAGCAAAGCAAGATGAAGATGTAGAAGATCGTCGGCGGACCTATAGTTTGGAGC
AACCTGGCCGAACGCACGGCCGGCAGCAGAGCGGGCGACACTCGGTGTCCGTGGAGGTGCAGATGCA
GCGGCAGCGGCAGGAGGAGCGGAGAGCTCCCAGCAGGCCAGCGCCAGTACAGCTCTCTGCCTCGGCAA
AGCAGGAAAAATGCCAGCTCGGTCTCCCAGGACTCTTGGGAGCAGAACTACTCCCCTGGGAAGGCTTCC
AGAGTGCCAAAGAGAACCCAGGTACTCCAGCTACCAAGGCTCCAGGAACGGCTACCTGGGAGGACATGG
CTTCAACGCCAGGTCATGCTGGAACTCAGGAGCTCCTTCGCCAGGAACAGAGGCGGAAGGAGCAGCAG
ATGAAGAAGCAGCCTCCTTCCGAGGGGCCAGCAACTATGACTCGTATAAGAAAGTCCAGGACCCAGT
ACGCCCCCTCCAAGGGGCCCTTCCGGCAAGATGTGCCCCCTCCCCTTCTCAGGTTGCGAGGCTGAACAG
ACTTCAGACTCCTGAGAAAGGGAGGCCCTTCTATTCC

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG230699 representing NM_001184787
 Red=Cloning site Green=Tags(s)

MKVTVCFGRTRVVVPCGDGHMKVFSLIQQAVTRYRKAIAKDPNYWIQVHRLHEHGGGILDLDLDDILCDVAD
 DKDRLVAVFDEQDPHHGGDGTASSTGTQSP E I F G S E L G T N N V S A F Q P Y Q A T S E I E V T P S V L R A N M P L H V
 R R S S D P A L I G L S T S V S D S N F S S E E P S R K N P T R W S T T A G F L K Q N T A G S P K T C D R K K D E N Y R S L P R D T S N W S
 N Q F Q R D N A R S S L S A S H P M V G K W L E K Q E Q D E D G T E E D N S R V E P V G H A D T G L E H I P N F S L D D M V K L V E V P N D
 G G P L G I H V V P F S A R G G R T L G L L V K R L E K G K A E H E N L F R E N D C I V R I N D G D L R N R R F E Q A Q H M F R Q A M R T
 P I I W F H V V P A A N K E Q Y E Q L S Q S E K N N Y S S R F S P D S Q Y I D N R S V N S A G L H T V Q R A P R L N H P P E Q I D S H S R
 L P H S A H P S G K P P S A P A S A P Q N V F S T T V S S G Y N T K K I G K R L N I Q L K K G T E G L G F S I T S R D V T I G G S A P I Y V
 K N I L P R G A A I Q D G R L K A G D R L I E V N G V D L V G K S Q E E V S L L R S T K M E G T V S L L V F R Q E D A F H P R E L N A E P
 S Q M Q I P K E T K A E D E D I V L T P D G T R E F L T F E V P L N D S G S A G L G V S V K G N R S K E N H A D L G I F V K S I I N G G A A
 S K D G R L R V N D Q L I A V N G E S L L G K T N Q D A M E T L R R S M S T E G N K R G M I Q L I V A R R I S K C N E L K S P G S P P G P E
 L P I E T A L D D R R R I S H S L Y S G I E G L D E S P S R N A A L S R I M G E S G K Y Q L S P T V N M P Q D D T V I I E D D R L P V L P
 P H L S D Q S S S S H D D V G F V T A D A G T W A K A A I S D S A D C S L S P D V D P V L A F Q R E G F G R Q S M S E K R T K Q F S D A S
 Q L D F V K T R K S K S M D L G I A D E T K L N T V D D Q K A G S P S R D V G P S L G L K S S S L E S L Q T A V A E V T L N G D I P F H R
 P R P R I I R G R C N E S F R A A I D K S Y D K P A V D D D E G M E T L E E D T E E S S R S G R E S V S T A S D Q P S H S L E R Q M N G
 N Q E K G D K T D R K K D K T G K E K K K D R D K E K D K M A K K G M L K G L G D M F R I Q A K T R E F R E R Q A R E R D Y A E I Q D F H
 R T F G C D D E L M Y G G V S S Y E G S M A L N A R P Q S P R E G H M M D A L Y A Q V K K P R N S K P S P V D S N R S T P S N H D R I Q R L
 R Q E F Q Q A K Q D E D V E D R R R T Y S F E Q P W P N A R P A T Q S G R H S V S V E V Q M Q R Q R Q E E R E S S Q Q A Q R Q Y S S L P R Q
 S R K N A S S V S Q D S W E Q N Y S P G E G F Q S A K E N P R Y S S Y Q G S R N G Y L G G H G F N A R V M L E T Q E L L R Q E Q R R K E Q Q
 M K K Q P P S E G P S N Y D S Y K K V Q D P S Y A P P K G P F R Q D V P P S P S Q V A R L N R L Q T P E K G R P F Y S

TRTRPLE - GFP Tag - V


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001184787

ORF Size: 3957 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001184787.1](#), [NP_001171716.1](#)

RefSeq Size: 5902 bp

RefSeq ORF: 3960 bp

Locus ID: 56288

UniProt ID: [Q8TEW0](#)

Cytogenetics: 10p11.22-p11.21

- Protein Pathways:** Adherens junction, Chemokine signaling pathway, Endocytosis, Neuroactive ligand-receptor interaction, Tight junction
- Gene Summary:** This gene encodes a member of the PARD protein family. PARD family members interact with other PARD family members and other proteins; they affect asymmetrical cell division and direct polarized cell growth. Multiple alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Oct 2011]