

Product datasheet for **RG230797**

TRAF1 (NM_001190947) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: TRAF1 (NM_001190947) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: TRAF1
Synonyms: EBI6; MGC:10353
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG230797 representing NM_001190947
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCTCCAGCTCAGGCAGCAGTCTCGCCCGGCCCTGATGAGAATGAGTTTCCCTTTGGGTGCCCTC
 CCACCGTCTGCCAGGACCAAAGGAGCCCAGGGCTCTCTGCTGTGCAGGCTGTCTCTGAGAACCCGAG
 GAATGGCGAGGATCAGATCTGCCCAAATGCAGAGGGGAAGACCTCCAGTCTATAAGCCAGGAAGCCGT
 CTTGCAACTCAGGAGAAGGCTCACCCGAGGTGGCTGAGGCTGGAATTGGGTGCCCTTTGCAGGTGTCG
 GCTGCTCCTTCAAGGGAAGCCACAGTCTGTGCAAGAGCATGAGGTCACCTCCCAGACCTCCCACCTAAA
 CCTGCTGTTGGGTTTCATGAAACAGTGAAGGCCCGGCTGGGCTGTGGCTGGAGTCTGGGCCCATGGCC
 CTGGAGCAGAACCTGTCAGACCTGCAGCTGCAGGCAGCCGTGGAAGTGGCGGGGACCTGGAGGTCGATT
 GCTACCGGGCACCTGCTCCGAGAGCCAGGAGGAGCTGGCCCTGCAGCACTTCATGAAGGAGAAGCTTCT
 GGCTGAGCTGGAGGGGAAGCTGCGTGTGTTGAGAACAATTGTTGCTGTCTCAACAAGGAGGTGGAGGCC
 TCCCACCTGGCCCTGGCCACCTATCCACCAGAGCCAGCTGGACCGTGAGCGCATCCTGAGCTTGGAGC
 AGAGGGTGGTGGAGCTTCAGCAGACCCTGGCCAGAAAGACCAGGCCCTGGGCAAGCTGGAGCAGAGCTT
 GCGCCTCATGGAGGAGGCTCCTTCGATGGCACTTTCTGTGGAAGATCACCAATGTCACCAGGCGGTGC
 CATGAGTCGGCCTGTGGCAGGACCGTCAGCCTTTCTCCCAGCCTTCTACACTGCCAAGTATGGCTACA
 AGTTGTGCTGCGGCTGTACCTGAATGGAGATGGCACTGGAAAGAGAACCATCTGTGCTCTTTCATCGT
 GATCATGAGAGGGAGTATGATGCGCTGCTGCCGTGGCCTTTCCGGAACAAGGTACCTTCATGCTGCTG
 GACCAGAACAAACCGTGAGCACGCCATTGACGCTTCCGGCCTGACCTAAGCTCAGCGTCTTCCAGAGGC
 CCCAGAGTGAAACCAACGTGGCCAGTGGATGCCCACTTCTTCCCCTCAGCAAAGTGCAGTACCCCAA
 GCACGCTACGTGAAGGACGACACAATGTTCTCAAGTGCATTGTGGAGACCAGCACT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG230797 representing NM_001190947
Red=Cloning site Green=Tags(s)

MASSSSGSSPRPAPDENEFPPGCPPTVCQDPKEPRALCCAGCLSENPRNGEDQICPKCRGEDLQSI SPGSR
 LRTQEKAHPEVAEAGIGCPFAGVGC SFKGS PQSVQEHEVTSQTSHLNL LLGFMKQWKARLGCGLESGPMA
 LEQNLSDLQLQA AVEVAGDLEVDCYRAPCSE SQEELALQHFMKEKLLAELEGKLRVFENI VAVLNKEVEA
 SHLALATSIHQSQLDRERILSLEQRVVELQQT LAQKDQALGKLEQSLRLMEEASF DGTFLWKITNVTRRC
 HESACGRTVSLFSPA FYTAKYGYKLC LRLYLNGDGTGKRTHLSL FIVIMRGEYDALLPWPFRNKVTFMLL
 DQNNREHAIDAFRPDLSSASFQRPQSETN VASGCPLFFPLSKLQSPKHAYVKDDTMFLKCIVETST


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001190947

ORF Size:	882 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001190947.1 , NP_001177876.1
RefSeq Size:	3715 bp
RefSeq ORF:	885 bp
Locus ID:	7185
UniProt ID:	Q13077
Cytogenetics:	9q33.2
Protein Families:	Druggable Genome
Protein Pathways:	Pathways in cancer, Small cell lung cancer
Gene Summary:	The protein encoded by this gene is a member of the TNF receptor (TNFR) associated factor (TRAF) protein family. TRAF proteins associate with, and mediate the signal transduction from various receptors of the TNFR superfamily. This protein and TRAF2 form a heterodimeric complex, which is required for TNF-alpha-mediated activation of MAPK8/JNK and NF-kappaB. The protein complex formed by this protein and TRAF2 also interacts with inhibitor-of-apoptosis proteins (IAPs), and thus mediates the anti-apoptotic signals from TNF receptors. The expression of this protein can be induced by Epstein-Barr virus (EBV). EBV infection membrane protein 1 (LMP1) is found to interact with this and other TRAF proteins; this interaction is thought to link LMP1-mediated B lymphocyte transformation to the signal transduction from TNFR family receptors. Three transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Jul 2010]