

## Product datasheet for **RG230921**

### CCDC179 (NM\_001195637) Human Tagged ORF Clone

#### Product data:

Product Type: Expression Plasmids  
 Product Name: CCDC179 (NM\_001195637) Human Tagged ORF Clone  
 Tag: TurboGFP  
 Symbol: CCDC179  
 Mammalian Cell Selection: Neomycin  
 Vector: pCMV6-AC-GFP (PS100010)  
 E. coli Selection: Ampicillin (100 ug/mL)  
 ORF Nucleotide Sequence: >RG230921 representing NM\_001195637  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGTGCCTGTATTGCTGGGACATCGAGCCTTCCCAAGTCAACCCTGAAGGACCAAGACAACATCATCCTT  
 CAGAGGTCACTGAGCGGCAGCTTGCAAATAAACGTATTCAGAATATGCAACACCTAAAGAAAGAGAAGAG  
 GAGACTGAATAAAAGGTTTTCAAGGCCTTCTCTATTCCAGAACCAGGACTCCTATGGTCATCT

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG230921 representing NM\_001195637  
 Red=Cloning site Green=Tags(s)

MCLYCWDIEPSQVNPEGPRQHHPSEVTERQLANKRIQNMQHLKKEKRRLNKRFSRPSPIPEPGLLWSS


**TRTRPLE** - GFP Tag - V

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM\_001195637

ORF Size: 204 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM\_001195637.2  
 RefSeq Size: 800 bp  
 RefSeq ORF: 207 bp  
 Locus ID: 100500938  
 UniProt ID: H3BU77  
 Cytogenetics: 11p14.3

## Product images:



Circular map for RG230921