

Product datasheet for **RG231090**

OXR1 (NM_001198534) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: OXR1 (NM_001198534) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: OXR1
Synonyms: CHEGDD; Nbla00307; TLDC3
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG231090 representing NM_001198534
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTCGTCTCTGGTATGGGAAAAAGGGAGAAGACATCAACCAATTAATCATAAATACACTCTGGTAG
TGTCACTGGCTGAGTATCACCGCAGGATCGATGCTCTAAATACTGAAGAACTGCGCACACTCTGCAGACG
CCTCCAGATTACTACAAGGAAGACATAAATCAAAGCAGGTTGCTACAGTGAAGCAGACCTGGAGTCT
GAATCTTTTCGACCAAACCTAAGTGATCCCACTGAACTTTACTGCCAGATCAAATGAAAAGCTTACCA
AGCATCTCCACCAAGAACAATTGGCTATCCATGGACTTTGTTTATGGTACTGGAAAACATGGCACAAG
CTTGAAAACCTTTATCGAACAATGACAGGTTTAGACACCCAGTGTGATGGTGATTAAGACAGTGAT
GGACAGGTTTTGGTGCCTTAGCATCTGAGCCACTGAAAGTGAGTGATGGCTTTTATGGTACTGGAGAGA
CCTTTGTTTTTACATTCTGTCCGGAGTTTGGGCTTTAAGTGGACAGGAGATAATATGTTTTTATCAA
AGGAGACATGGATTCACTAGCTTTCCGGTGGTGGAGGAGGAGAAATTTGCGCTTTGGCTTGATGGAGATCTC
TACCATGGAAGAAGCCATTCTTGTAACGTTTGGGAATCGTACACTTTCTAAGAAGGAAGATTTCTTTA
TCCAAGATATTGAAATCTGGGCTTTTGAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG231090 representing NM_001198534
Red=Cloning site Green=Tags(s)

MSRLWYGKKGRRHQPINHKYTLVVSVAEYHRRIDALNTEELRTLCRRLQITTREDINSKQVATVKADLES
 ESFRPNLSDPSELLLPDQIEKLTKHLPPRTIGYPWTLVYGTGKHGTSKTLYRTMTGLDTPVLMVIKDS
 GQVFGALASEPLKVSDGFYGTGETFVFTECFEVEVFKWTGDNMFFIKGDMDSLAFGGGGGEFALWLDGDL
 YHGRSHSCKTFGNRTLKSKEDFFIQDIEIWAFF

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_001198534

ORF Size: 729 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001198534.1 , NP_001185463.1
RefSeq Size:	2635 bp
RefSeq ORF:	732 bp
Locus ID:	55074
UniProt ID:	Q8N573
Cytogenetics:	8q23.1
Gene Summary:	May be involved in protection from oxidative damage.[UniProtKB/Swiss-Prot Function]