

Product datasheet for **RG231171**

LIMS1 (NM_001193483) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LIMS1 (NM_001193483) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	LIMS1
Synonyms:	PINCH; PINCH-1; PINCH1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG231171 representing NM_001193483 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGGGCGTGGCGGCCGGGATGACCCACAGCAACATGGCCAACGCCCTGGCCAGCGCCACTTGCAGC
GCTGCAAGGGCGGCTTTGCGCCCGCTGAGAAGATCGTGAACAGTAATGGGGAGCTGTACCATGAGCAGTG
TTTCGTGTGCGCTCAGTGCTTCCAGCAGTCCAGAAAGGACTCTTCTATGAGTTGAAGGAAGAAAGTAC
TGTGAACATGACTTTCAGATGCTCTTTGCCCTTGCTGTCATCAGTGTGGTGAATTCATCATTGGCCGAG
TTATCAAAGCCATGAATAACAGCTGGCATCCGGAGTGCTCCGCTGTGACCTCTGCCAGGAAGTTCTGGC
AGATATCGGGTTTGTCAAGAATGCTGGGAGACACCTGTGTCGCCCTGTCATAATCGTGAGAAAGCCAGA
GGCCTTGGGAAATACATCTGCCAGAAATGCCATGCTATCATCGATGAGCAGCCTCTGATATTCAAGAACG
ACCCCTACCATCCAGACCATTTCAACTGCGCCAACCTGCGGGAAGGAGCTGACTGCCGATGCACGGGAGCT
GAAAGGGGAGCTATACTGCCTCCCATGCCATGATAAAATGGGGTCCCATCTGTGGTGTCTTGGCAGCG
CCCATCGAAGGGCGCGTGGTGAACGCTATGGGCAAGCAGTGCCATGTGGAGCATTTTGTGGTGGCAAGT
GTGAGAAACCTTTCTTGGACATCGCCATTATGAGAGGAAAGGCTGGCATATTGTGAAACTCACTATAA
CCAGCTATTTGGTGTGTTTGTCTTCCACTGCAATCGTGTATAGAAGGTGATGTGGTCTCTGCTCTTAAT
AAGGCTTGGTGCCTGAACTGCTTTGCCTGTTCTACCTGCAACATAAATAACACTCAAGAATAAGTTTG
TGGATTTGACATGAAGCCAGTCTGTAAGAAGTCTATGAGAAATTTCCATTGGAGCTGAAGAAAAGACT
TAAGAACTAGCTGAGACCTTAGGAAGGAAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG231171 representing NM_001193483
 Red=Cloning site Green=Tags(s)

MLGVAAGMTHSNMANALASATCERCKGGFAPAEEKIVNSNGELYHEQCFVCAQCFQFPEGLFYEFEGRKY
 CEHDFQMLFAPCCHQCGEFIIGRVIKAMNNSWHPECFRCDLCQEVLADIGFVKNAGRHLCRPCHNREKAR
 GLGKYICQKCHAIIDEQPLIFKNDPYHPDHFNCANCGKELTADARELKGELYCLPCHDKMGVPICGACRR
 PIEGRVNVAMGKQWHVEHFVCAKCEKPFLLGHRHYERKGLAYCETHYNQLFGDVCFHCNRVIEGDVVSALN
 KAWCVNCFACSTCNTKLLKNKFVEFDMKPVCKCYEKFPLELKKRLKLAETLGRK

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001193483

ORF Size: 1011 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001193483.3](#)

RefSeq Size: 4438 bp

RefSeq ORF: 1014 bp

Locus ID: 3987


UniProt ID: [P48059](#)

Cytogenetics: 2q12.3

Protein Families: Druggable Genome

Gene Summary: The protein encoded by this gene is an adaptor protein which contains five LIM domains, or double zinc fingers. The protein is likely involved in integrin signaling through its LIM domain-mediated interaction with integrin-linked kinase, found in focal adhesion plaques. It is also thought to act as a bridge linking integrin-linked kinase to NCK adaptor protein 2, which is involved in growth factor receptor kinase signaling pathways. Its localization to the periphery of spreading cells also suggests that this protein may play a role in integrin-mediated cell adhesion or spreading. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2010]

Product images:



Circular map for RG231171