

Product datasheet for **RG231303**

CPNE1 (NM_001198863) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CPNE1 (NM_001198863) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CPNE1
Synonyms:	COPN1; CPN1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG231303 representing NM_001198863
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCCCACTGCGTGACCTTGGTTTCAGCTGTCCATTTCTGTGACCATCTCATTGACAAGGACATCGGCT
 CCAAGTCTGACCACTCTGCGTCTTTTACAGGATGTGGGAGGGGCAGCTGGGCTGAGCTTGCCGGAC
 TGAACGGGTGCGAACTGCTCAAGCCCTGAGTTCTCCAAGACTCTACAGCTTGAGTACCGCTTTGAGACA
 GTCCAGAAGCTACGCTTTGGAATCTATGACATAGACAACAAGACGCCAGAGCTGAGGGATGATGACTTCC
 TAGGGGGTGTGAGTGTCCCTAGGACAGATTGTGTCCAGCCAGGTAAGTACTCTCCCTTGATGCTGAA
 GCCTGGAAAACCTGTGGGCGGGGACCATCACGGTCTCAGCTCAGGAATTAAGGACAATCGTGTAGTA
 ACCATGGAGGTAGAGGCCAGAACTAGATAAAGAGGACTTCTGGGAAAATCAGATCCATTTCTGGAGT
 TCTTCCGCCAGGGTATGGGAAATGGCACCTGGTGTACAGATCTGAGGTCATCAAGAACAACCTGAACCC
 TACATGGAAGCGTTTCTCAGTCCCGTTTACGATTTCTGTGGTGGAAACCCAGCACACCCATCCAGGTG
 CAATGCTCCGATTATGACAGTGACGGGTACATGATCTCATCGGTACCTCCACACCAGCTTGCCCGCAGC
 TGCAGGACGTCCCGCTGAGTTTGAATGCATCCACCTGAGAAGCAGCAGAAAAAGAAAGCTACAGAA
 CTCTGGAATACTCCGTGTCAAGATTTGTCCGGTAGAAACAGAGTACTCTTTCTGGACTATGTGATGGGA
 GGCTGTGAGTCAACTTCACTGTGGCGTGGACTTCACTGGTCCAATGGAGACCCTCCTCACCTGACT
 CCCTACACTACCTGAGTCCAACAGGGGTCAATGAGTACCTGATGGCACTGTGGAGTGTGGCAGCGTGGT
 TCAGGACTATGACTCAGACAAGCTGTCCCTGCATTTGGATTTGGGGCCAGGTTCCCTGACTGGCAG
 GTCTCGCATGAATTTGCCTTGAATTTCAACCCAGTAAACCCCTACTGTGCAGGCATCCAGGGCATTGTGG
 ATGCTACCGCCAAGCCCTGCCCAAGTTCGCCTCTATGGCCCTACCAACTTTGCACCCATCATCAACCA
 TGTGGCCAGGTTTGCAGCCAGGCTGCACATCAGGGGACTGCCTCGTACTTCATGCTGTTGCTGCTGACT
 GATGGTGTGTGACGGATGTGGAAGCCACAGTGGGCTGTGGTGCCTCGAACCTGCCATGTCAG
 TGATCATTGTGGGTGTGGGTGGTGTGACTTTGAGGCCATGGAGCAGCTGGACGCTGATGGTGGACCCCT
 GCATACACGTTCTGGCAGGCTGCTGCCCGGACATTGTGCAATTTGTACCCTACCGCGGTTCCAGAAT
 GCCCTCGGGAGGCAATTGGCACAGACCGTGTGCGAGAAGTGGCCACAACTGGTCTCATACTTCAGGG
 CCCAGGTTGGGCCCGCTCAAGCCACTCCACCCTCAGCCAAGGATCTGCACAGGCCCCCGAGGCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG231303 representing NM_001198863
 Red=Cloning site Green=Tags(s)

MAHCVTLVQLSISCDHLIDKDIGSKSDPLCVLLQDVGGGWAELGRTERVRNCSSPEFSKTLQLEYRFET
 VQKLRFGIYDIDNKTPELRDDDFLGAECSLGQIVSSQVLTPLMLKPGKPAGRGITVSAQELKDNRVV
 TMEVEARNLDDKDFLGSDFLEFFRQGDGKWHLYYRSEVIKNNLNPTWKRF SVPVQHFCCGNPSTPIQV
 QCSYDSDGSHDLIGTFHTSLAQLQAVPAEFECIHPEKQKQKSYKNSGTIRVKICRVETEYSFLDYVMG
 GCQINFYVGVDFTGSDPSSPDSLHYLSPTGVNEYLMALWSVGSVVQDYSDKLFPAFGFGAQVPPDWQ
 VSHEFALNFNPNPYCAGIQGIVDAYRQALPQVRLYGPTNFAPIINHVARFAAQAHQGTASYFMLLLLT
 DGAVTDVEATREAVVRASNLPM SVIIVGVGGADFEAMEQLDADGGPLHTRSGQAAARDIVQFVPPYRRFQN
 APREALAQTVLAEVPTQLVSYFRAQGWAPLKPLPPSAKDPAQAPQA

TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-MluI

ACCN:	NM_001198863
ORF Size:	1608 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001198863.2
RefSeq Size:	1977 bp
RefSeq ORF:	1611 bp
Locus ID:	8904
Cytogenetics:	20q11.22
Gene Summary:	Calcium-dependent membrane-binding proteins may regulate molecular events at the interface of the cell membrane and cytoplasm. This gene encodes a calcium-dependent protein that also contains two N-terminal type II C2 domains and an integrin A domain-like sequence in the C-terminus. However, the encoded protein does not contain a predicted signal sequence or transmembrane domains. This protein has a broad tissue distribution and it may function in membrane trafficking. This gene and the gene for RNA binding motif protein 12 overlap at map location 20q11.21. Alternate splicing results in multiple transcript variants encoding different proteins. [provided by RefSeq, Aug 2008]