

Product datasheet for **RG231332**

NT5C1B (NT5C1B-RDH14) (NM_001199104) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NT5C1B (NT5C1B-RDH14) (NM_001199104) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	NT5C1B-RDH14
Synonyms:	AIRP; cN-IB; cN1B; NT5C1B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG231332 representing NM_001199104
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAGTCAAACATCTCTCAAACAGAAAAAAGATGAGCCTGGAATGAGGTCCTCAAAGAGAGTCTAGAAG
 CAGAAAAAAGAAAGGAATCTGACAAAACAGGAGTTCGTCTGAGCAATCAGATGAGGCGTGCAGTCAATCC
 GAATCACTCGCTGAGATGTTGCCCTTCCAGGGTCACTCGTCGTGTAGACGCTGCCTTTGTGACGCTGAG
 GGAACAGCCCTTGGCCCTGCCACACAATACGTATTTATATTACATGTGCCTGTTGTGGGAGCAGGGCC
 AGCAGATCACCATGATGAGGGGATCACAAGAATCATCACTGCGGAAGACAGACTCTCGAGGGTACCTTGT
 GCGCAGTCAATGGTCTAGAATATCCCGGAGCCATCCACCAAGGCTCCATCCATAGATGAGCCTAGAAGC
 AGGAACACCAGTCTAAGCTCCCAGCAGCTCCACGAGCTCCCGACTCCATCCACCTCCCCAAGCCTGC
 ATGACTCCTCACCGCCGCGCTGTCCGGCAGCCCTCGTCCAGCCACCCGCGTCCGCCAGCTGCCCGG
 GTCGCTGGACTCGCGCCCTCCCACGCCCCAGAGCCCGATCCTGGCTCCCGGCGCAGCACAAAATGCAA
 GAGAATCCGGAGGCTGGGCCAAAGGCATCGTGCGGGAAATCCGCCAGACCCGGGACTCGCAGCCGCTGG
 AATATTCGCGCACGTCCCCACCGAGTGGAAAGTCTCCAGCCAGCGCAGGGGGATCTACCCCGCCTCCAC
 CCAGCTGGACCGCAACTCTGTCCGAGCAGCAGCAGCAGCAGCGGGAGGACGAGGACGACTACGAGGCT
 GCCTACTGGGCATCCATGAGGTCGTTCTACGAGAAGAACCCGAGCTGCTCGCGCCCTGGCCGCCAAAC
 CCAAGAACGCCATCACCATTGCTCTCTCATCTGCGCGCTCTTCAACATGGTGGACGGCAGGAAAATCTA
 CGAGCAAGAGGGTCTGGAAAAGTACATGGAGTATCAGCTACCAATGAGAACGTCATCTGACCCCGGGC
 CCGCGCTCCGCTTCGTCAAGGCACTACAGTATGTCAATGTAGACTCCGTGATCTATATCTGATGAAC
 AGGACTTATTTGATATTGTACTGATGACTAATAACCATGCCAAGTGGGAGTGGCGCTTATAAACACGCT
 CAATCACTACGGCTTACTGATTGACCGCTTCTGTCTGACCGGGGAAAAGACCCATTGGCTATTTGAAG
 GCATATCTTACCAACTTGTATATTGCTGCAGATTCTGAAAAAGTCAAGAGGCAATACAAGAAGGTATTG
 CCTCTGCGACAATGTTTGTGAGGCAAAAGACATGGCTTACTGTGACACTCAGCTCCGTGTAGCCTTTGA
 TGGGGATGCTGCTCTTCTCTGATGAGTCTGAACATTTTACCAAGGAGCATGGGCTCGACAAAATCTTC
 CAGTATGATACATTATGTGAAAGTAAGCCTCTTGCTCAGGGTCCCCTAAAAGGCTTTCTGGAAGATTAG
 GCAGACTGCAAAAGAAGTCTATGCCAAAATGAACGGTACTTTGCTCTATCAGGACCTACCTGGTTAC
 AGCTAGGAGTGCAGCCAGTTCAGGCGCCCGTGTGCTGAAGACCCTTCGACGCTGGGGTCTAGAGATAGAC
 GAAGCTCTTTTCTTGTGGAGCCCCAAAAGTCCCATCTTGGTGAAGATCCGGCCCCACATCTTCTTTG
 ATGACCACATGTTCCACATTGAAGGGGCACAGAGGAAGAGCCTAGGCTGGATGTCT

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG231332 representing NM_001199104
 Red=Cloning site Green=Tags(s)


MSQTSLKQKNEPGRSSKESLEAEKRKESDKTGVRLSNQMRRAVNPNSLRCCPFQGHSSCRRRCLCAAE
 GTALGPCHTIRIYIHMCLLWEQQQITMMRGSQESSLRKTDNRGVLVRSQWSRISRSPSTKAPSIDEP
 RNTSAKLPSSSTSSRTPSTSPSLHDSPPPLSGQPSLQPPASPQLPRSLDSRPPTPEPDPGSRRTKMQ
 ENPEAWAQGIVREIRQTRDSQPLEYSRTSPTWKSQRRGIYPASTQLDRNSLSEQQQQREDEDDYEA
 AYWASMSRFYEKNPSCSRPWPPKPKNAITIALSSCALFNMVDGRKIYEQEGLEKMEYQLTENVILTPG
 PAFRFVKALQYVNARLRDLYPDEQDLFDIVLMTNNHAQVGVRLINSVNHVGLLIDRFCLTGGKDP
 IGYLKYLAADSEKVEAIQEGIASATMFDGAKDMAYCDTQLRVAFDGDAVLFSDSEHFKEHGLDKFF
 QYDTLCEKPLAQGLKGFLEDLGRLQKFKYAKNERLLCPIRTYLVARSASSGARVLTLLRWGLEID
 EALFLAGAPKSPILVKIRPHIFFDDHMFHIEGAQRKSLGWMS

TRTRPLE – GFP Tag – V


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001199104

ORF Size: 1806 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001199104.1 , NP_001186033.1
RefSeq Size:	2982 bp
RefSeq ORF:	1809 bp
Locus ID:	100526794
UniProt ID:	Q96P26
Cytogenetics:	2p24.2
Gene Summary:	This locus represents naturally occurring read-through transcription between the neighboring NT5C1B (5'-nucleotidase, cytosolic IB) and RDH14 (retinol dehydrogenase 14) genes on chromosome 2. Alternative splicing results in multiple transcript variants, one of which encodes a fusion protein that shares sequence identity with the products of each individual gene. [provided by RefSeq, Nov 2010]