

Product datasheet for **RG231343**

MID1 (NM_033289) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MID1 (NM_033289) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MID1
Synonyms:	BBBG1; FXY; GBBB1; MIDIN; OGS1; OS; OSX; RNF59; TRIM18; XPRF; ZNFXY
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG231343 representing NM_033289
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAACACTGGAGTCAGAAGTACCTGCCCATTGTCTGGAGCTCTTTGAGGACCTCTTCTACTG
 CCTGCGCACACAGCCTCTGCTTCAACTGCGCCACCGCATCCTAGTATCACACTGTGCCACCAACGAGTC
 TGTGGAGTCCATCACCGCCTTCCAGTGCACCCACCTGCCGGCATGTCATCACCCCTAGCCAGCGAGGTCTA
 GACGGGCTCAAGCGCAACGTCACCCACAGAACATCATCGACAGGTTCCAGAAAGCATCAGTGAGCGGGC
 CCAACTCTCCAGCGAGACCCGTCGGGAGCGGGCCTTTGACGCCAACACCATGACCTCCGCCGAGAAGGT
 CCTCTGCCAGTTTTGTGACCAGGATCCTGCCAGGACGCTGTGAAGACCTGTGTCACCTGTGAAGTATCC
 TACTGTGACGAGTGCCTGAAAGCCACTACCCGAATAAGAAGCCCTTTACAGGCCATCGTCTGATTGAGC
 CAATCCGGACTCTCACATCCGGGGCTGATGTGCTTGGAGCATGAGGATGAGAAGCAAACTTAGAGAG
 TAACCTCACCAACCTTATTAAGAGGAACACAGAAGTGGAGACCTTTTGGCTAAACTCATCAAACCTGT
 CAACATGTTGAAGTCAATGCATCAGTCAAGAAGCCAAATTTGACAGAGGAGTGTGATCTTCTCATTGAGA
 TCATTCAGCAAAGACGACAGATTATTGGAACCAAGATCAAGAAGGGAAGGTGATGAGGCTTCGAAACT
 GGCTCAGCAGATTGCAAACGAGTGCATGAGCGGTGACGATCACTCATCTCCAAAGCGGAACAC
 TCTCTGAAGGAGAATGATCATGCGGTTTCTACAGACTGCTAAGAATATCACCGAGAGAGTCTCCATGG
 CAACTGCATCCTCCAGGTTCTAATTCCTGAAATCAACCTCAATGACACATTTGACACCTTTGCCTTAGA
 TTTTTCCGAGAGAAGAACTGCTAGAATGTCTGGATTACCTTACAGCTCCCAACCTCCACAATTAGA
 GAAGAGCTCTGCACAGCTTACATGACACCATCACTGTGCATTGGACCTCCGATGATGAGTTCAGCGTGG
 TCTCTACGAGCTCCAGTACACCATATTCACCGGACAAGCCAACGTCGTTAGTCTGTGTAATTCGGCTGA
 TAGCTGGATGATAGTACCAACATCAAGCAGAACCCTACACGGTGCACGGTCTGCAGAGCGGCACCAAG
 TACATCTTCATGGTCAAGGCCATCAACCAGGCGGCGAGCCGACGAGTGCCTGGGAAGTTGAAGACAA
 ACAGCCAACCATTTAAACTGGATCCCAAATCTGCTCATCGAAAAGTGAAGGTGTCCCATGATAAATTGAC
 AGTAGAACGTGATGAGTCATCATCAAGAAGAGTACACACCTGAACGCTTACCAGCCAGGGGAGCTAT
 GGAGTAGCTGAAATGTGTTTATTGATAGTGGCCGCATTATTGGGAAGTGGTATAAGTGAAGACAT
 GGTATGCCATTGGTCTTGCTTACAAATCAGCCCCGAAGCATGAATGGATTGGGAAGAACTGCTTCTCTG
 GCGCTCTGCCGCTGCAACAATAACTGGTGGTGGAGACACAATAGCAAGGAAATCCCATGAGCCTGCC
 CCCCACCTCCGGCGCTGGGCATCCTGCTGGACTATGATAACGGCTCTATCGCCTTTTATGATGCTTTGA
 ACTCCATCCACCTCTACACCTTCGACGTCGCATTTGCGCAGCCTGTTGCCCCACCTTACCGTGTGGAA
 CAAGTGTCTGACGATTATCACTGGGCTCCCTATCCAGACCATTGGACTGCACAGAGCAGCTGCCG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG231343 representing NM_033289
 Red=Cloning site Green=Tags(s)

METLESELTCPICLELFEPLLLPCAHSFCNCAHRILVSHCATNESVESITAFQCPTCRHVITLSQRGL
 DGLKRNVTLQNIIDRFQKASVSGPNPSETRRERAFDANTMTSAEKVLCQFCDQDPAQDAVKTCVTCEVS
 YCDECLKATHPNKPFTHRLIEPIPDHIRGLMCLHEHEDEKQNLNLTNLIKRNTLETLLAKLIQTC
 QHVEVNASRQEAKLTECDLLIEIIQRRQIIGTKIKEGKVMRLRKLAAQQIANCKQCIERSASLISQAEH
 SLKENDHARFLQAKNITERVSMATASSQVLIPEINLNDTFDFTALDFSREKKLLECLDYLTPNPPTIR
 EELCTASYDTITVHWTSDDEFVSVSYELQYTI FTGQANVVSLENSADSWMIVPNIKQNHVTVHGLQSGTK
 YIFMVKAINQAGSRSSPEGKLTNSQPFKLDPKSAHRKLVSHDNLTVRDESSSKSHTPERFTSQGSY
 GVAGNVFIDSGRHYWEVVISGSTWYAIGLAYKSAPKHEWIGKNSASWALCRCNNWVVRHNSKEIPIEPA
 PHLRRVGIILLDYDNGSIAFYDALNSIHL YTFDVAFAQVPCPTFTVWNKCLTIITGLPIPDLDCTEQLP

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_033289

ORF Size: 1887 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_033289.2](#)

RefSeq Size: 6349 bp

RefSeq ORF: 1890 bp

Locus ID: 4281


Cytogenetics: Xp22.2

Protein Families: Druggable Genome

Protein Pathways: Ubiquitin mediated proteolysis

Gene Summary: The protein encoded by this gene is a member of the tripartite motif (TRIM) family, also known as the 'RING-B box-coiled coil' (RBCC) subgroup of RING finger proteins. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. This protein forms homodimers which associate with microtubules in the cytoplasm. The protein is likely involved in the formation of multiprotein structures acting as anchor points to microtubules. Mutations in this gene have been associated with the X-linked form of Opitz syndrome, which is characterized by midline abnormalities such as cleft lip, laryngeal cleft, heart defects, hypospadias, and agenesis of the corpus callosum. This gene was also the first example of a gene subject to X inactivation in human while escaping it in mouse. Alternative promoter use, alternative splicing and alternative polyadenylation result in multiple transcript variants that have different tissue specificities. [provided by RefSeq, Dec 2016]

Product images:



Circular map for RG231343