

Product datasheet for **RG231360**

AGAP9 (NM_001190810) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AGAP9 (NM_001190810) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	AGAP9
Synonyms:	AGAP-9; bA301J7.2; CTGLF6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG231360 representing NM_001190810
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGGGAACATACTGACCTGTCGTGTGCACCCTAGCGTCAGCCTCGAGTTTGACCAGCAGCAGGGGTCGG
TGTGTCCCTCTGAATCTGAGATCTATGAGGCAGGAGCTGAGGACAGGATGGCAGGAGCGCCCATGGCTGC
TGCTGTACAGCCTGCTGAGGTGACTGTTGAAGTTGGTGAGGACCTCCACATGCACCAGGTTCTGTGACCGG
GAGATGCCTGAAGCTTTGGAGTTAACCTTTCTGCCAATCCAGAGGCAAGCACAATATCCAGAGGAACT
CTCAAACAGATGCTTTGGAGTTAACTCTTCTGCCAATCCAGAGGCAAGCACAATATCCAGAGGAACTC
TCAAACAGATGTTGTAGAAATAAGAAGAAGCAACTGTACAAACCATGTATCTACTGAGCGTTTCAGTCAA
CAATACAGCTCGTTCGACAATATTCCTTGATGACAGCACAGCCATCCAGCATTATCTTACAATGACAA
TAATATCTGTGACCTTGGAGATACCTCATCATATCACACAAAGAGATGCAGATAGATCTTTGAGCATAACC
TGATGAACAGTTACACTATTTGCGGTTTCCACCGTGCACATTATGAAGAAAAGAAATGGAGGTGGGAGT
TTAAATAACTATTCCTCCTCCATTCCACCGACTCCCAGCACCAGCCAGGAGGACCCTCAGTTCAGTGTTCC
CTCCCCTGCAACACACCCACCCCGTTTGCAAGCGGTCCATGCGCTGGTCCAACCTGTTTACATCTGA
GAAAGGGAGTGACCCAGACAAAGAGAGGAAAGCCCGGAGAATCATGCTGACACCATCGGGAGCGGCAGA
GCCATCCCATTAACAGGGCATGCTCTTAAAGCGAAGTGAAAAATGGCTGAAGACATGGAAAAAGAAAT
ACGTCACCCTGTGTTCCAATGGCGTGCTCACCTATTATTCAAGCTTAGGTGATTATGAAGAATATTCA
TAAAAAGAGATTGACCTTCGGACATCTACCATCAAAGTCCCAGGAAAGTGGCCATCCCTAGCCACATCA
GCCTGTGCACCCATCTCCAGCTCTAAAAGCAATGGCCTATCCAAGGACATGGACACCGGGCTGGTGACT
CCATATGCTTCAGCCCCGGTATCTCCAGCACACCAGCCCCAAGCTCAACCCGCCCTCTCCTCATGC
CAATAAAAAGAAACACCTAAAGAAGAAAAGCACCAACAACCTTTATGATTGTGTCTGCCACTGGCCAAACG
TGGCACTTTGAAGCCACGACGTATGAGGAGCGGGATGCCTGGGTCCAAGCCATCCAGAGCCAGATCCTGG
CCAGCCTGCAGTCATGCAAGAGCAGTAAAAGCAAGTCCCAGCTGACCAGCCAGAGTGAGGCCATGGCCCT
GCAGTCGATCCAAAACATGCGTGGGAACGCCACTGTGTGGACTGTGAGACCCAGAATCCTAAGTGGGCC
AGTTTGAAGTGGGAGTCTCATGTGATTGAATGCTCAGGTATCCACCGCAGTTTTGGCACCCGCTTT
CCCGTGTGCGATCTCTGGAGCTGGATGACTGGCCAGTTGAGCTCAGGAAGTTATGTCATCTATTGGCAA
TGAGCTAGCCAACAGCATCTGGGAAGGGAGCAGCCAGGGCAGACAAAACCTCAATAAAGTCCACGAGG
GAAGAGAAGGAATGGTGGATCCGTTCAAATATGAGGAGAAGCTTTTCTGGCCCACTACCTGCACTG
AGCTGTCCCTGGGCCAGCAGCTGCTGCGGGCCACCACTGATGAGGACCTGCAGACAGCCATCCTGTGCT
GGCACATGGCTCCCGTGAAGGAGTGAACGAGACCTGTGGGGAGGGAGACGGCTGCACGGCACTCCATCTG
GCCTGCCCAAGGGGAATGTGGTCTGGAGCAGCTCCTGACGGGGTGGACGTCATGGCCCGAGATGCCCA
CGGAACACAGCGC

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG231360 representing NM_001190810
 Red=Cloning site Green=Tags(s)


MGNILTCRVHPSVSLFDQQQGSVCPSESEIYEAGAEDRMAGAPMAAAVQPAEVTVEVGEDLHMHQVRDR
 EMPEALEFNLSANPEASTIFQRNSQTDALFNSSANPEASTIFQRNSQTDVVEIRRSNCTNHVSTERFSQ
 QYSSCSTIFLDDSTAIQHLYLMTIISVTLIPHHITQRDADRSLIPDEQLHSFAVSTVHIMKKRNGGGS
 LNNYSSSIPPTPSTSQEDPQFSVPPTANTPTPVCKRSMRWSNLFSEKGSDDPKERKAPENHADIIGSGR
 AIPKQGMMLKRSKGLKTKWKKYVTLCSNGVLTYYSSLGDYMKNIHKKEIDLRTSTIKVPGKWPSLATS
 ACAPISSSKSNGLSKDMDTGLGDSICFSPGISSTTSPKLNPPSPHANKKKHLKKSTNNFMIVSATGQT
 WHFEATTYEERDAWVQAIQSQILASLQCKSSKSKSQLTSQSEAMALQSIQNMGRNAHCVDCEQNPKWA
 SLNLGVLMCIECSGIHRSGTRLSRVRSLELDDWPVELRKMSSIGNELANSIWEGSSQGGQTKPSIKSTR
 EEKEWIRSKYEEKLFLAPLPCTELSLGQQLLRATTDDELQTAILLLAHGSREEVNETCGEGDGCTALHL
 ACRKGNVVLEQLLTGWTSWPEMPTGTQR

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_001190810

ORF Size: 1974 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001190810.1](#), [NP_001177739.1](#)

RefSeq Size: 2387 bp
RefSeq ORF: 1977 bp
Locus ID: 642517
UniProt ID: [Q5VTM2](#)
Cytogenetics: 10q11.22
Gene Summary: Putative GTPase-activating protein.[UniProtKB/Swiss-Prot Function]