

Product datasheet for **RG231373**

STRA6 (NM_001199041) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	STRA6 (NM_001199041) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	STRA6
Synonyms:	MCOPCB8; MCOPS9; PP14296
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG231373 representing NM_001199041
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGACTCTCAAGGGGACTGGGCTGCACAGGAGAAGGGCCAGAGAATGTCGTCCCAGCCAGCAGGGAACC
 AGACCTCCCCCGGGGCCACAGAGGACTACTCCTATGGCAGCTGGTACATCGATGAGCCCCAGGGGGCGA
 GGAGCTCCAGCCAGAGGGGAAGTGCCCTCCTGCCACACCAGCATAACCACCCGGCCTGTACCACGCCTGC
 CTGGCCTCGTGTCAATCCTTGTGCTGCTGCTCCTGGCCATGCTGGTGGGGCGCCGCCAGCTCTGGCCTG
 ACTGTGTGCGTGGCAGGCCCGCCCTGCCAGCCCTGTGGATTTCTTGGCTGGGGACAGGCCCCGGGCAGT
 GCCTGCTGCTGTTTTATGGTCTCCTGAGCTCCCTGTGTTTGTGCTCCCCGACGAGGACGCATTGCC
 TTCCTGACTCTCGCCTCAGCACCCAGCCAAGATGGGAAAAGTGGGCTCCAAGAGGGGCCTGGAAGATAC
 TGGGACTGTTCTATTATGCTGCCCTCTACTACCCTCTGGCTGCCTGTGCCACGGCTGGCCACACAGCTGC
 ACACCTGCTCGGCAGCAGCTGTCTGGGCCACCTTGGGGTCCAGGTCTGGCAGAGGGCAGAGTGTCCC
 CAGGTGCCAAGATCTACAAGTACTACTCCCTGCTGGCCTCCCTGCCTCTCTGCTGGGCCCTCGGATTCC
 TGAGCCTTTGGTACCCTGTGCAGCTGGTGAAGAAGCTTCAGCCGTAGGACAGGAGCAGGCTCCAAGGGGT
 GCAGAGCAGTACTCTGAGGAATATCTGAGGAACCTCCTTTCAGGAAGAAGCTGGGAAGCAGTACCAC
 ACCTCCAAGCATGGCTTCTGTCTGGGCCCGCTCTGCTTGAGACACTGCATCTACACTCCACAGCCAG
 GATTCATCTCCCGCTGAAGCTGGTGTCTTTCAGCTACACTGACAGGGACGGCCATTTACCAGGTGGCCCT
 GCTGCTGCTGGTGGCGTGGTACCCACTATCCAGAAGTGAGGGCAGGGGTACCACGGATGTCTCTAC
 CTGCTGGCCGGCTTTGGAATCGTCTCTCCGAGGACAAGCAGGAGGTGGTGGAGCTGGTGAAGCACCATC
 TGTGGGCTCTGGAAGTGTGCTACATCTCAGCCTTGGTCTTGTCTGCTTACTCACCTTCTGGTCTGAT
 GCGCTCACTGGTGACACACAGGACCAACCTTCGAGCTCTGCACCGAGGAGCTGCCCTGGACTTGAGTCCC
 TTGCATCGGAGTCCCATCCCTCCCGCCAAGCCATATTCTGTTGGATGAGCTTCAGTGCCTACCAGACAG
 CCTTTATCTGCCTTGGGCTCCTGGTGCAGCAGATCATCTTCTTCTGGGAACCACGGCCCTGGCCTTCT
 GGTGCTCATGCCTGTGCTCCATGGCAGGAACCTCCTGCTTTCGTTCCCTGGAGTCTCGTGGCCCTTC
 TGGCTGACTTTGGCCCTGGCTGTGATCCTGCAGAACATGGCAGCCATTGGGTCTTCTGGAGACTCATG
 ATGGACACCCACAGCTGACCAACCGCGAGTGCTCTATGCAGCCACCTTCTTCTTCTTCCCCTCAATGT
 GCTGGTGGGTGCCATGGTGGCCACCTGGCGAGTGCTCCTCTGCCCTCTACAACGCCATCCACCTTGGC
 CAGATGGACCTCAGCCTGCTGCCACCGAGAGCCGCACTCTGACCCCGGCTACTACACGTACCGAAACT
 TCTTGAAGATTGAAGTACGCCAGTCGCATCCAGCCATGACAGCCTTCTGCTCCCTGCTCCTGCAAGCGCA
 GAGCCTCTACCCAGGACCATGGCAGCCCCCAGGACAGCCTCAGACCAGGGGAGGAAGACGAAGGGATG
 CAGCTGTACAGACAAAGGACTCCATGGCCAAGGGAGCTAGGCCCGGGCCAGCCGCGGCAGGGCTCGCT
 GGGGTCTGGCCTACACGCTGCTGCACAACCAACCCTGCAGGTCTTCCGCAAGACGGCCCTGTTGGGTGC
 CAATGGTGGCCAGCCC

ACGCGTACGCGGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG231373 representing NM_001199041
 Red=Cloning site Green=Tags(s)

```
MDSQGDWAAQEKQRMSSQPAGNQTSPGATEDYSYGSWYIDEPQGGEEELQPEGEVPSCHTSIPPGLYHAC
LASLSILVLLLLAMLVRRRQLWPDCVRGRPLSPVDFLAGDRPRAVPAAVFMVLLSSCLLLPDEDALP
FLTLASAPSQDGKTEAPRGAWKILGLFYAALYPLAACATAGHTAAHLLGSTLSWAHLGVQVWQRAECP
QVPKIYKYSSLLASLPLLLGLGFLSLWYPVQLVRSFSRRTGAGSKGLQSSYSEEYLRNLLCRKKLGSSYH
TSKHGFLSWARVCLRHCIYTPQPGFHLPLKLVLSATLTGTAIYQVALLLLVGVVPTIQKVRAGVTTDVSY
LLAGFGIYVLSQDQEVVELVKHHLWALEVCYISALVLSCLLTFVLVLRSLVTHRTNLRALHRGAALDLS
LHRSPHPSRQAIFCWMFSAYQTAFCICLGLLVQQIIFFLGTALAFVLMVPLHGRNLLFRSLESSWPF
WLTALAVILQNMAAHWVLETHDGHPQLTNRRLVYAATFLLFPLNVLVGAMVATWRVLLSALYNAIHLG
QMDLSLLPPRAATLDPGYTYRNFLEKIEVSQSHPAMTAFCSLLQAQSLLPRTMAAPQDSLPRGEEDEGM
QLLQTKDSMAKGARPGASRGRARWGLAYTLLHNPTLQVFRKTALLGANGAQP
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001199041

ORF Size: 2046 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001199041.1](#), [NP_001185970.1](#)

RefSeq Size: 2993 bp

RefSeq ORF: 2049 bp

Locus ID: 64220


UniProt ID: [Q9BX79](#)

Cytogenetics: 15q24.1

Protein Families: Transmembrane

Gene Summary: The protein encoded by this gene is a membrane protein involved in the metabolism of retinol. The encoded protein acts as a receptor for retinol/retinol binding protein complexes. This protein removes the retinol from the complex and transports it across the cell membrane. Defects in this gene are a cause of syndromic microphthalmia type 9 (MCOPS9). Several transcript variants encoding a few different isoforms have been found for this gene. [provided by RefSeq, Dec 2008]

Product images:



Circular map for RG231373