

## Product datasheet for **RG231409**

### Enconsin (MAP7) (NM\_001198608) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Enconsin (MAP7) (NM_001198608) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Enconsin
Synonyms:	E-MAP-115; EMAP115
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RG231409 representing NM\_001198608  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**

ATGCCTGGATCAGCTACAGCTCTCCGACATGAGAGACTGAAGAAGACCAATGCAAGGCCAATTCCTCTTG  
 GTTTATTCACCATTAATGAGGAAGACGAAACAGCAAAAAGAATGAAAATTCAGAAGACAAAAGCACCCGA  
 CAGCTACAAAGTGCAAGATAAGAAAAATGCCTCCAGCCGCCCTGCCTCTGCAATTCAGGACAAAATAAC  
 AACCACTCAGGAAATAAACAGACCCTCCGCCTGTGTTACGTGTTGATGACCGGCAGCGGCTGGCCCGGG  
 AGCGACGTGAGGAACGGGAGAAACAGCTAGCTGCAAGAGAAATAGTGTGGTTAGAAAAGAGAAGAGCGAGC  
 CAGGCAGCACTACGAGAAGCACCTGGAAGAGCGGAAGAAGAGGTTGGAGGAGCAGAGGCAGAAGGAGGAG  
 CGGAGGAGGGCTGCTGTGGAGGAGAAGCGGAGGCAGAGACTTGAGGAGGACAAAGAACGCCACGAAGCTG  
 TTGTACGGCGCACAAATGAAAAGGAGCCAGAAGCCAAAACAGAAGCATAACCGTTGGTCGTGGGGAGGCTC  
 TCTCCATGGGAGCCCTAGCATCCACAGTGCAGATCCAGACAGGCGGTGAGTTCCACCATGAATCTTTTCG  
 AAATATGTTGATCCCGTCATTAGCAAGCGGCTCTCCTCTTCATCTGCAACTTTACTAAATTCCTCAGATA  
 GAGCTCGCCCGCTGCAGCTCAGCCCATGGGAGAGCAGCGTTGTTAACAGACTCCTGACGCCACACATTC  
 GTTCTGGCCAGAAGTAAAAGCACAGCTGCCTTGTCTGGAGAAGCAGCATCTTGACAGCCCATCATCATG  
 CCCTACAAAGCTGCACACTCTAGAAATTCGATGGATCGACAAAACCTCTTTGTAACACCCTGAGGGCT  
 CTTCTCGCAGGAGGATCATTTCATGGCACAGCGAGCTATAAAAAAGAAAGAGAGAGAGAAAATGACTCTT  
 CCTCACATCTGGCACCCGAAGGGCTGTATCTCCATCTAATCCCAAAGCAAGACAACCAGCTCGTCCCGA  
 CTTTGGCTTCCGTCCAAGTCTCTTCTCATTGCTGCGACACCCAGACCCGACATCCTCCTGCCACCCG  
 GCTCAGTCAAAGCTGCTCCTGCTCAGTCCGGCCCCATCCCCCGCAACATCCGCCCTGTCAAGAGGGA  
 AGTCAAAGTGGAGCCTGAGAAGAAAGATCCTGAGAAGGAACCTCAGAAAGTTGCCAATGAGCCCTACTA  
 AAGGGCAGAGCACCTTTAGTGAAGGTAGAAGAAGCCACAGTTGAAGAGCGGACACCTGCTGAACAGAAG  
 TTGGCCCTGCTGCTCCAGCCATGGCCCCAGCTCCAGCCTCGGCCCCAGCTCCAGCCTCGGCCCCAGCTCC  
 AGCCCCGGTCCCCACCCAGCCATGGTCTCAGCCCCGTCCACTGTGAATGCCAGTGCTTCTGTAAAG  
 ACTTCTGCAGGCACCACCCAGCCAGAGGAGGCCACAAGGCTTCTAGCTGAGAAGAGGCGGCTGGCCCGAG  
 AGCAGAGAGAAAAGGAAGAAAGGGAGAGGAGGGAGCAGGAAGAGCTTGAAGACAAAAGAGAGAGGAATT  
 GGCTCAACGTGTGGCTGAAGAGAGGACGACTCGCCGTGAGGAGGAGTCCGCGAGGCTGGAAGCCGAGCAG  
 GCCCGGGAGAAGGAGGAGCAGCTGCAGCGCAGCGGAGGAGCGGCGCTGCGCGAGCGGGAGGAGGCAG  
 AGCGCGCCAGAGGCAGAAAGAAGAAGAAGCTCGCGTTTCGTGAAGAAGCAGAGAGGGTCCGCGCAGGAACG  
 AGAGAAGCATTCCAGAGAGAAGAGCAAGAGCGCCTGGAGAGAAAAGAGCAGTTGAGGAGATTATGAAA  
 AGAACCAGGAGAACAGAGCTACAGATAAGAAAACCAAGTATCAGAGAAACGGTATATAGCCAAGGGAG  
 CTCTCACTGGAGGAACAGAGGTGTCTGCACTTCCATGTACAACAACCGCTCCGGGAAATGGAAAGCCAGT  
 TGGCAGCCACATGTGGTTACCTCACACCAGTCAAAGTGACAGTGGAGAGCACTCCCGATTTGGAAAAA  
 CAACCAATGAAAATGGTGTATCTGTTCCAGAAATGAAAATTTTGAAGAAATATAAACTTACCCATTGGAT  
 CTAACCATCCAGATTAGATGTCACCAACAGTGAAGAGCCAGAAATTCCTTTGAATCCAATTTGGCCCTT  
 TGATGATGAAGGGACACTTGGGCCCTGCCTCAGGTAGATGGTGTTCAGACACAGCAGACTGCAGAAGTT  
 ATA

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA


**Protein Sequence:** >RG231409 representing NM\_001198608  
 Red=Cloning site Green=Tags(s)

MPGSATALRHERLKKTNARPIPLGLFTINEEDEQQKNGNSRRPKAPDSYKVQDKKNASSRPASAI SGQNN  
 NHSGNKPDPPPVLRVDDRQLARERREEREKQLAAREIVWLEREERARQHYEKHLEERKKRLEEQRQKEE  
 RRRAAVEEKRRQRLEEDKERHEAVVRRTMERSQKPKQKHNRWSWGGSLHGSPSIHSADPDRRSVSTMNLS  
 KYVDPVISKRLSSSSATLLNSPDRARRLLQLSPWESSVVRLLTPTHSFLARSKSTAALSGEAASCSP IIM  
 PYKAAHSRNSMDRPFVTPPEGSSRRRIIHGTASYKKERENVLFLTSGTRRAVSPSNPKARQPARSR  
 LWLPSKSLPHLPGTPTSSLPPGSVKAAPAQVRPPSPGNIRPVKREVKVEPEKKDPEKEPQKVANEPSL  
 KGRAPLVKVEEATVEERTPAEPEVGAAPAMAPAPASAPAPASAPAPVPTPAMVSAPSSTVNASASVK  
 TSAGTTDPEEATRLLAEKRRLAREQREKEERERREQEELERQKREELAQRVAERTTRREEESRRLEAEQ  
 AREKEEQLRQAEEALREREEAERAQRQKEEEARVREEAERVQREREKHFQREEQERLERKKRLEEIMK  
 RTRRTEATDKKTSQQRNGDIAGALTGTEVSALPCTTNPAGNGKPVGSPHVVTSHQSKVTVESTPDLEK  
 QPNENGVSQNFEEIINLPIGSKPSRLDVTNSESPEIPLNPILAFDDEGLGPLPQVDGVQTQQTAEV  
 I

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001198608


**ORF Size:** 2313 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u>NM_001198608.1, NP_001185537.1</u>
<b>RefSeq Size:</b>	4512 bp
<b>RefSeq ORF:</b>	2316 bp
<b>Locus ID:</b>	9053
<b>UniProt ID:</b>	<u>Q14244</u>
<b>Cytogenetics:</b>	6q23.3
<b>Protein Families:</b>	Druggable Genome
<b>Gene Summary:</b>	The product of this gene is a microtubule-associated protein that is predominantly expressed in cells of epithelial origin. Microtubule-associated proteins are thought to be involved in microtubule dynamics, which is essential for cell polarization and differentiation. This protein has been shown to be able to stabilize microtubules, and may serve to modulate microtubule functions. Studies of the related mouse protein also suggested an essential role in microtubule function required for spermatogenesis. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2010]

### Product images:



Circular map for RG231409