

## Product datasheet for **RG231412**

### Enconsin (MAP7) (NM\_001198609) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Enconsin (MAP7) (NM_001198609) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Enconsin
Synonyms:	E-MAP-115; EMAP115
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RG231412 representing NM\_001198609  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGCCTGGATCAGCTACAGCTCTCCGACATGAGAGACTGAAGAAGACCAATGCAAGGCCAATTCCTCTTG  
 GTTTATTCACCATTAATGAGGAAGACGAAACAGCAAAAAGAATGGAAATTCAGAAGACAAAAGCACCCGA  
 CAGCTACAAAGTGCAAGATAAGAAAAATGCCTCCAGCCGCCCTGCCTCTGCAATTCAGGACAAAATAAC  
 AACCACTCAGGAAATAAACAGACCCTCCGCTGTGTTACGTGTTGATGACCGGCAGCGGCTGGCCCGGG  
 AGCGACGTGAGGAACGGGAGAAACAGCTAGCTGCAAGAGAAATAGTGTGGTTAGAAAAGAGAAGAGCGAGC  
 CAGGCAGCACTACGAGAAGCACCTGGAAGAGCGGAAGAAGAGGTTGGAGGAGCAGAGGCAGAAGGAGGAG  
 CGGAGGAGGGCTGCTGTGGAGGAGAAGCGGAGGCAGAGACTTGAGGAGGACAAAGAACGCCACGAAGCTG  
 TTGTACGGCGCACAAATGGAAGGAGCCAGAAGCCAAAACAGAAGCATAACCGTTGGTCTGGGGAGGCTC  
 TCTCCATGGGAGCCCTAGCATCCACAGTGCAGATCCAGACAGGCGGTGAGTTCCACCATGAATCTTTTCG  
 AAATATGTTGATCCCGTCATTAGCAAGCGGCTCTCCTCTTCATCTGCAACTTTACTAAATTTCTCCAGATA  
 GAGCTCGCCCGCTGCAGCTCAGCCCATGGGAGAGCAGCGTTGTTAACAGACTCCTGACGCCACACATTC  
 GTTCTGGCCAGAAGTAAAAGCACAGCTGCCTTGTCTGGAGAAGCAGTTATCCCCATTTGTCTCTGTTCA  
 GCATCTTGAGCCCATCATCATGCCCTACAAAGCTGCACACTCTAGAAATTCGATGGATCGACAAAAC  
 TCTTTGTAACACCACCTGAGGGCTCTTCTCGCAGGAGGATCATTTCATGGCACAGCGAGCTATAAAAAAGA  
 AAGAGAGAGAAAAATGTACTCTTCTCACATCTGGCACCCGAAGGGCTGTATCTCCATCTAATCCCAAA  
 GCAAGACAACAGCTCGCTCCCGACTTTGGCTTCCGTCCAAGTCTCTTCTCATTGCTGGCACACCCA  
 GACCGACATCCTCTTCCACCCGGCTCAGTCAAAGCTGCTCCTGCTCAGGTCCGGCCCCATCCCCCGG  
 CAACATCCGCCCTGTCAAGAGGGAAGTCAAAGTGGAGCCTGAGAAGAAAGATCCTGAGAAGGAACCTCAG  
 AAAGTTGCCAATGAGCCCTCACTAAAGGCAGAGCACCTTTAGTGAAGGTAGAAGGCCACAGTTGAAG  
 AGCGGACACCTGCTGAACCAGAAGTTGGCCCTGCTGCTCCAGCCATGGCCCCAGCTCCAGCCTCGGCCCC  
 AGCTCCAGCCTCGGCCCCAGCTCCAGCCCCGGTCCCCACCCAGCCATGGTCTCAGCCCCGTATCCACT  
 GTGAATGCCAGTGCTTCTGTTAAGACTTCTGCAGGCACCACCGACCCAGAGGAGGCCACAAGGCTTCTAG  
 CTGAGAAGAGGCGGCTGGCCCGAGAGCAGAGAGAAAAGGAAGAAAGGGAGAGGAGGGAGCAGGAAGAGCT  
 TGAAGACAAAAGAGAGAGGAATTGGCTCAACGTGTGGCTGAAGAGAGGACGACTCGCCGTGAGGAGGAG  
 TCGCGCAGGCTGGAAGCCGAGCAGGCCCGGAGAAAGGAGGAGCAGCTGCAGCGGCAGGCGGAGGAGCGGG  
 CGCTGCGCGAGCGGGAGGAGGCAGAGCGGCCAGAGGCAGAAAGAAGAAGCTCGCGTTCTGTAAGA  
 AGCAGAGAGGGTCCGGCAGGAACGAGAGAAGCATTTCAGAGAGAAGAGCAAGAGCGCCTGGAGAGAAAAG  
 AAGCGACTTGAGGAGATTATGAAAAGAACCAGGAGAACAGAAGCTACAGATAAGAAAACCAAGTATCAGA  
 GAAACGGTGATATAGCCAAGGGAGCTCTCACTGGAGGAACAGAGGTGTCTGCACTTCCATGTACAACAAA  
 CGCTCCGGGAAATGGAAGCCAGTTGGCAGCCACATGTGGTTACCTCACACCAGTCAAAAGTGACAGTG  
 GAGAGCACTCCCGATTTGGAAAAACAACCAATGAAAATGGTGTATCTGTTTCAAGTAAAAATTTGAAG  
 AAATTATAAACTTACCCATTGGATCTAAACCATCCAGATTAGATGTCACCAACAGTGAGAGCCCAGAAAT  
 TCCTTTGAATCCAATTTTGGCCTTTGATGATGAAGGGACACTTGGGCCCTGCCTCAGGTAGATGGTGT  
 CAGACACAGCAGACTGCAGAAGTTATA

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >RG231412 representing NM\_001198609  
 Red=Cloning site Green=Tags(s)


```

MPGSATALRHERLKKTNARPIPLGLFTINEEDEQQKNGNSRRPKAPDSYKVQDKKNASSRPASAI SGQNN
NHSGNKPPPPVLRVDDRQLARERREEREKQLAAREI VWLEREERARQHYEKHLEERKKRLEEQRQKEE
RRRAAVEEKRRQRLEEDKERHEAVVRRTMERSQPKQKHNRWSWGGSLHGSPSIHSADPDRRSVSTMNLS
KYVDPVISKRLSSSSATLLNSPDRARRLLQSPWESSVVNRLLP THTSFLARSKSTAALSGEAVIPI CPRS
ASCSP IIMPYKAAHSRNSMDRPKLFVTPPEGSSRRRI IHGTASYKKERENVLFLTSGTRRAVSPSNPK
ARQPARSRLWLP SKSLPHLPGTPRPTSSLPPGSVKAAPAQVRPPSPGNIRPVKREVKVEPEKDKPEKEPQ
KVANESLKGRAPLVKVEEATVEERTPAEPEVGPAA PAMAPAPASAPAPASAPAPVPTPAMVSAPSST
VNASASVKTSAGTTDPEEATRL LAEKRR LAREQREKEERERREQEELERQKREELAQRVAEERTTRREEE
SRRLAEQAREKEEQ LQRQAEERALREREEAERAQRQKEEEARVREEAERVQREREKHFQREEQERLERK
KRLEEIMKRTRRTEATDKKTS DQRNGDIAKGAL TGGTEVSALPCTT NAPGNGKPVGSPHVVTSHQSKVTV
ESTPDLEKQPNENGVS VQNFEEIINLP I GSKPSRLDVTNSE SPEIPLNPILAFDDEGTLGPLQVDGV
QTQQTAEVI
    
```

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001198609


**ORF Size:** 2337 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u>NM_001198609.1, NP_001185538.1</u>
<b>RefSeq Size:</b>	4536 bp
<b>RefSeq ORF:</b>	2340 bp
<b>Locus ID:</b>	9053
<b>UniProt ID:</b>	<u>Q14244</u>
<b>Cytogenetics:</b>	6q23.3
<b>Protein Families:</b>	Druggable Genome
<b>Gene Summary:</b>	The product of this gene is a microtubule-associated protein that is predominantly expressed in cells of epithelial origin. Microtubule-associated proteins are thought to be involved in microtubule dynamics, which is essential for cell polarization and differentiation. This protein has been shown to be able to stabilize microtubules, and may serve to modulate microtubule functions. Studies of the related mouse protein also suggested an essential role in microtubule function required for spermatogenesis. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2010]

### Product images:



Circular map for RG231412