

Product datasheet for **RG231584**

PKIB (NM_001270395) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: PKIB (NM_001270395) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: PKIB
Synonyms: PRKACN2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG231584 representing NM_001270395
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCACACCAGGATGTTGCTATGAGGACAGATTCATCAAAAATGACTGACGTGGAGTCTGGGGTCGCCA
ATTTTGCATCTTCAGCAAGGGCAGGCCCGGAATGCCTTACCAGACATCCAGAGTTCAGCTGCCACAGA
CGGAACCTCAGATTTGCCCTCAAACCTGGAGGCTCTCTCCGTGAAGGAAGATGCAAAGAGAAAGATGAA
AAAACAACACAAGACCAATTGAAAAGCCTCAAATGAAGAAAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG231584 representing NM_001270395
Red=Cloning site Green=Tags(s)
MSHQDVAMRTDSSKMTDVESGVANFASSARAGRNRNALPDIQSSAATDGTSDLPLKLEALSVKEDAKEKDE
KTTQDQLEKPQNEEK

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001270395

ORF Size: 255 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001270395.1</u> , <u>NP_001257324.1</u>
RefSeq Size:	1614 bp
RefSeq ORF:	258 bp
Locus ID:	5570
UniProt ID:	<u>Q9C010</u>
Cytogenetics:	6q22.31
Protein Families:	Druggable Genome
Gene Summary:	This gene encodes a member of the cAMP-dependent protein kinase inhibitor family. The encoded protein may play a role in the protein kinase A (PKA) pathway by interacting with the catalytic subunit of PKA, and overexpression of this gene may play a role in prostate cancer. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jul 2012]