

Product datasheet for **RG231595**

MAX (NM_001271068) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: MAX (NM_001271068) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: MAX
Synonyms: bHLHd4
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG231595 representing NM_001271068
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGCGATAACGATGACATCGAGGTGGAGAGCGACGCTGACAAACGGGCTCATCATAATGCACTGGAAC
GAAAACGTAGGGACCACATCAAAGACAGCTTTCACAGTTTGC GGACTCAGTCCCATCACTCCAAGGAGA
GAAGCTCTATTTCTCTTTGGAAATTGTGACTCCTGTCCTTCATCGTCAAAGTTTGATGCAGAAATGC
CACACCTTCATTTCAAGCTACCAAGTGCACAAGAAAAAGAATGCAAGATT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG231595 representing NM_001271068
Red=Cloning site Green=Tags(s)
MSDNDIEVESDADKRAHHNALERKRRDHKDSFHSRLRDSVPSLQGEKLYFLFWKLCPTVLRQSLMQKC
HTFISSYQVHKKECKI

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001271068

ORF Size: 261 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

| | |
|-------------------------------|--|
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_001271068.1 , NP_001257997.1 |
| RefSeq Size: | 910 bp |
| RefSeq ORF: | 264 bp |
| Locus ID: | 4149 |
| UniProt ID: | P61244 |
| Cytogenetics: | 14q23.3 |
| Protein Families: | Druggable Genome, Transcription Factors |
| Protein Pathways: | MAPK signaling pathway, Pathways in cancer, Small cell lung cancer |
| Gene Summary: | The protein encoded by this gene is a member of the basic helix-loop-helix leucine zipper (bHLHZ) family of transcription factors. It is able to form homodimers and heterodimers with other family members, which include Mad, Mxi1 and Myc. Myc is an oncoprotein implicated in cell proliferation, differentiation and apoptosis. The homodimers and heterodimers compete for a common DNA target site (the E box) and rearrangement among these dimer forms provides a complex system of transcriptional regulation. Mutations of this gene have been reported to be associated with hereditary pheochromocytoma. A pseudogene of this gene is located on the long arm of chromosome 7. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2012] |