

## Product datasheet for **RG231612**

### UBE2G2 (NM\_001202489) Human Tagged ORF Clone

#### Product data:

Product Type: Expression Plasmids  
Product Name: UBE2G2 (NM\_001202489) Human Tagged ORF Clone  
Tag: TurboGFP  
Symbol: UBE2G2  
Synonyms: UBC7  
Mammalian Cell Selection: Neomycin  
Vector: pCMV6-AC-GFP (PS100010)  
E. coli Selection: Ampicillin (100 ug/mL)  
ORF Nucleotide Sequence: >RG231612 representing NM\_001202489  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAGATTTACCTGTGAGATGTTTCATCCCAACATCTACCCTGATGGGAGAGTCTGCATTTCCATCCTCC  
ACGCGCCAGGCGATGACCCCATGGGCTACGAGAGCAGCGCGGAGCGGTGGAGTCTGTGCAGAGTGTGGA  
GAAGATCCTGCTGTCGGTGGTGTGAGCATGCTGGCAGAGCCCAATGACGAAAGTGGAGCTAACGTGGATGCG  
TCCAAAATGTGGCGCATGACCGGGAGCAGTTCTATAAGATTGCCAAGCAGATCGTCCAGAAGTCTCTGG  
GACTG

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG231612 representing NM\_001202489  
Red=Cloning site Green=Tags(s)  
MRFTCEMFHPNIYPDGRVCISILHAPGDDPMGYESSAERWSPVQSVEKILLSVVSMLAEPNDESGANVDA  
SKMWRDDREQFYKIAKQIVQKSLGL

**TRTRPLE** - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001202489

ORF Size: 285 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

|                               |  |
|-------------------------------|--|
| <b>OTI Annotation:</b>        | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).   |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>  |
| <b>RefSeq:</b>                | <a href="#">NM_001202489.2</a>   |
| <b>RefSeq Size:</b>           | 3318 bp  |
| <b>RefSeq ORF:</b>            | 288 bp   |
| <b>Locus ID:</b>              | 7327   |
| <b>UniProt ID:</b>            | <a href="#">P60604</a>   |
| <b>Cytogenetics:</b>          | 21q22.3  |
| <b>Protein Families:</b>      | Druggable Genome   |
| <b>Protein Pathways:</b>      | Parkinson's disease, Ubiquitin mediated proteolysis  |
| <b>Gene Summary:</b>          | The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. The encoded protein shares 100% sequence identity with the mouse counterpart. This gene is ubiquitously expressed, with high expression seen in adult muscle. Three alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jan 2011] |