

Product datasheet for **RG231668**

ZNF 559 (ZNF559) (NM_001202409) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ZNF 559 (ZNF559) (NM_001202409) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: ZNF559
Synonyms: NBLA00121
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG231668 representing NM_001202409
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGGCCATTGGTCTTCCCTTTGTGATCATCTTTCTCAAGATGTTTTCTGTCTTCATGAGTCAAAT
TTGAAGAGGAAAGGATGGTGGCTGGGTGGTTGACAAATTACTCTCAGGACTCAGTGACCTTTGAGGATGT
GGCTGTGGACTTCACCCAGGAGGAGTGGACTTTGCTGGATCAAACCTCAGAGAACTTATACAGAGATGTG
ATGCTGGAGAACTATAAGAATCTAGTTGCAGTAGATTGGGAGAGTCATATTAATACCAAATGGTCAGCAC
CTCAGCAGAATTTTTGCAGGGAAAACATCCAGTGTGGTGGAAATGAATTCAGAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG231668 representing NM_001202409
Red=Cloning site Green=Tags(s)
MLAIGLSLCDHLSQDVFCLHESKFEEERMVAGWLTNYSQDSVTFEDVAVDFTQEEWTLLDQTQRNLYRDV
MLENYKNLVAVDWESHINTKWSAPQQNFLQKTSVVMNSE

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001202409

ORF Size: 336 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001202409.1</u> , <u>NP_001189338.1</u>
RefSeq Size:	3309 bp
RefSeq ORF:	339 bp
Locus ID:	84527
Cytogenetics:	19p13.2
Protein Families:	Transcription Factors
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]